

# ABX MICROS 60

## TECHNICAL MANUAL


Part number :  
RAA 009 A Ind. B



# • *REVISIONS*

<i>INDEX</i>	<i>P/N REVISION</i>	<i>REVISION</i>	<i>SECTION</i>	<i>DATE</i>
A		Creation	all	30/03/98
B	RAH 581 AA	Miscellaneous corrections	1, 3, 4	07/10/98

# • *CONTENTS*

- 
1. HYDRAULIC & PNEUMATIC PRINCIPLES
  2. ELECTRIC & ELECTRONIC PRINCIPLES
  3. MAINTENANCE PROCEDURES
  4. OUTPUT FORMAT
  5. TRAINING SLIDES

# 1. HYDRAULIC & PNEUMATIC PRINCIPLES

## CONTENTS

1. GENERALITIES .....	2
2. MICROS 60 OT HYDRAULIC .....	4
2.1. With bottles .....	4
2.1.1. Transmission tubes list .....	4
2.1.2. Hydropneumatic connections .....	5
2.2. Pack .....	6
2.2.1. Transmission tubes list .....	6
2.2.2. Hydropneumatic connections .....	7
2.3. Hydraulic cycle description .....	8
2.3.1. Atmosphere circuit .....	8
2.3.2. Diluent circuit .....	8
2.3.3. Clean circuit .....	8
2.3.4. Lyse circuit .....	8
2.3.5. WBC/RBC counting circuit .....	8
2.3.6. Drain/bubbling circuit .....	8
3. MICROS 60 CT HYDRAULIC .....	15
3.1. With bottles .....	15
3.1.1. Transmission tubes list .....	15
3.1.2. Hydropneumatic connections .....	16
3.2. Pack .....	17
3.2.1. Transmission tubes list .....	17
3.2.2. Hydropneumatic connections .....	18
3.3. Hydraulic cycle description .....	19
3.3.1. Atmosphere circuit .....	19
3.3.2. Diluent circuit .....	19
3.3.3. Clean circuit .....	19
3.3.4. Lyse circuit .....	19
3.3.5. WBC/RBC counting circuit .....	19
3.3.6. Drain/bubbling circuit .....	19
4. PNEUMATIC DIAGRAMS .....	26
4.1. Micros 60 CT bottle version .....	26
4.2. Micros 60 CT pack version .....	26
4.3. Micros 60 OT bottle version .....	26
4.4. Micros 60 OT pack version .....	26

# 1. HYDRAULIC & PNEUMATIC PRINCIPLES

## 1. GENERALITIES

MICROS 60 instrument has been designed for simple mechanical operations.


4 stepper motors provide movements to mechanical assemblies.

Pressure and vacuum are provided by the vacuum/waste syringe up and down movements (diag 1).

Liquid movements are achieved either by means of mechanical assembly movements (liquid



Diag.1



Diag.2

syringes diag 2) or pressure/vacuum and simultaneous action of specific valves.

### • Dilution chambers (Diag .3)

WBC/HGB and RBC chambers are made of GRILAMID TR55 LY injected.

The diode and the cell of the spectrophotometer are glued on the WBC/HGB chamber.

Chamber positions can be modified in order to obtain the best sampling position possible.

#### *Dilutions :*

First dilution is carried out in the WBC/HGB chamber (with a bubbling phasis).

The RBC blood sample is aspirated from this dilution.

Lyse is sent from the drain nipple of the WBC/HGB chamber.

Counts have a duration of 2 x 6 seconds.

(see procedures RAS 188 A and RAS 187 A for cycle hydraulic details)



#### *Rinse :*

To obtain the best rinse in the counting heads, diluent is sent from the liquid syringes. This is carried out before, between and after the two counts.

# 1. HYDRAULIC & PNEUMATIC PRINCIPLES

## **IMPORTANT**

A window on the HGB/WBC chamber allows the needle to move down into the chamber and to inject reagents. As important light or variation of light can cause HGB result drifts, close instrument cover and door before running blood analyses.

### *Bubbling :*

Insulators avoid polluted liquid overflows during bubbling phasis. they also allows an accurate adjustment of the bubbling volume.

### *MICROS 60 CT specifics :*

- The piercing needle is equipped with two injectors to obtain a homogeneous diluent flow during needle rinsing phasis (see procedures RAS 188 A and RAS 187 A).
- Atmosphere is provided to sample tubes to allow a correct aspiration of blood.

# 1. HYDRAULIC & PNEUMATIC PRINCIPLES

## 2. MICROS 60 OT HYDRAULIC

Liquid circuits, hydropneumatic connections, as well as the transmission tubes used, are described in the following chart tables.

### 2.1. With bottles

#### 2.1.1. Transmission tubes list

DESIGNATION	PART NUMBER	DIAMETER	LENGTH	QUANTITY
SLEEVE HPS3	DBD005A	5-9		0.5
T CONNECTOR	EAB006B	2.3		3
T CONNECTOR	EAB032A	1.5		1
TUBE CAP	EAC017A	2.5		1
TYGON TUBE 0.051"	EAE006A	1.30	140	1
TYGON TUBE 0.051"	EAE006A	1.30	350	1
TYGON TUBE 0.060"	EAE007A	1.52	20	1
TYGON TUBE 0.060"	EAE007A	1.52	40	3
TYGON TUBE 0.060"	EAE007A	1.52	50	1
TYGON TUBE 0.060"	EAE007A	1.52	60	1
TYGON TUBE 0.060"	EAE007A	1.52	70	2
TYGON TUBE 0.060"	EAE007A	1.52	80	1
TYGON TUBE 0.060"	EAE007A	1.52	150	1
TYGON TUBE 0.060"	EAE007A	1.52	170	1
TYGON TUBE 0.060"	EAE007A	1.52	220	1
TYGON TUBE 0.060"	EAE007A	1.52	240	1
TYGON TUBE 0.060"	EAE007A	1.52	300	1
TYGON TUBE 0.060"	EAE007A	1.52	370	1
TYGON TUBE 0.060"	EAE007A	1.52	410	1
TYGON TUBE 0.060"	EAE007A	1.52	450	2
TYGON TUBE 0.060"	EAE007A	1.52	480	1
TYGON TUBE 0.081"	EAE008A	2.05	20	1
TYGON TUBE 0.081"	EAE008A	2.05	35	1
TYGON TUBE 0.081"	EAE008A	2.05	200	1
TYGON TUBE 0.081"	EAE008A	2.05	330	1
TYGON TUBE 0.081"	EAE008A	2.05	1080	1
TYGON TUBE 0.090"	EAE009A	2.28	20	2
TYGON TUBE 0.090"	EAE009A	2.28	50	1
TYGON TUBE 0.090"	EAE009A	2.28	60	2
TYGON TUBE 0.090"	EAE009A	2.28	120	1
TYGON TUBE 0.090"	EAE009A	2.28	140	1
TYGON TUBE 0.090"	EAE009A	2.28	150	1
TYGON TUBE 0.090"	EAE009A	2.28	190	1
BLUE TYGON TUBE 0.090"	EAE036A	2.28	1100	1
SLEEVE	GAL098A			30
TUBE SHIELD	GBC088A	4.4	30	1
TUBE SHIELD	GBC088A	4.4	60	1

# 1. HYDRAULIC & PNEUMATIC PRINCIPLES

## 2.1.2. Hydropneumatic connections

CIRCUIT	FROM	SLEEVE	DIAM	LENGTH	TO	SLEEVE
AIR	(atmosphere)		2.28	190	Liquid valve 2_2	
	Liquid valve 2_1	Y	2.28	50	Waste-chamb_1	Y
	Waste-chamb_2		2.28	20	Cap EAC017A	
DILUENT	Diluent input		2.05	330	Liquid valve11_1	
	Liquid valve11_3	Y	1.52	220	Temp. sensor xba281a	Y
	Temp. sensor. xba281a	Y	1.52	40	Liquid syringes_3	Y
	Liquid valve11_2	Y	1.52	40	Liquid valve10_3	Y
	Liquid valve10_1	Y	1.52	70	Liquid valve7_3	Y
	Liquid valve7_1	Y	1.52	370	needle rinsing block_2	
	Needle rinsing block_1		1.52	410	Liquid valve8_1	
	Liquid valve8_2	Y	1.52	300	Waste-chamb_3	
	Liquid valve 7_2	Y	1.52	240	Té 2.3_1	
	Té 2.3_2		1.52	50	Bac WBC/HGB_3	
	Liquid valve10_2	Y	1.30	140	Liquid syringes_1	Y
	Liquid syringes_2	Y	1.30	350	needle_1	
	Liquid syringes_2	Y			Needle_1	Y
CLEAN	Clean bottle		2.28	1100 blue	Liquid valve4_2	Y
	Liquid valve4_1	Y	1.52	450	T connector 2.3_3	
LYSE	Lyse bottle		2.05	1080	liquid valve 1_1	
	Liquid valve1_3		1.52	150	Liquid syringes_4	
	Liquid valve1_2		1.52	480	WBC grounding connector	
	WBC grounding connector		1.52	20	T connector 1.5_1	
WBC/RBC COUNTING	WBC/HGB_2 chamber		1.52	170	RBC_3 chamber	
	RBC_2 chamber		1.52	450	liquid valve 6_2	Y
	Liquid valve 6_1	Y	1.52	60	Waste-chamb_4	
DRAIN/BUBBLING	WBC/HGB_1 chamber		1.52	40	T connector 1.5_2	
	WBC/HGB_1 chamber		gbc088a	30	T connector 1.5_2	
	T connector 1.5_3		1.52	80	Insulator WBC_1	
	Insulator WBC_2		2.28	120	liquid valve12_2	Y
	Liquid valve12_1	Y	2.05	35	T connector 2.3_1	
	RBC_1 chamber		1.52	70	Insulator RBC_1	
	RBC_1 chamber		gbc088a	60	Insulator RBC_1	
	Insulator RBC_2		2.28	150	liquid valve13_2	Y
	Liquid valve13_1	Y	2.05	20	T connector 2.3_2	
	T connector 2.3_3		2.05	200	T connector 2.3_2	
	T connector 2.3_3		2.28	20	Cell xba199a	Y
	Cell xba199a	Y	2.28	60	Waste-chamb_5	Y
	T connector 2.3_1		2.28	60	Liquid valve5_2	Y
Liquid valve5_1		2.28	140	Waste ouput		

### NOTE

Read this table as follows in this example :  
The Liquid valve 7\_2 corresponds to the ouput 2 of the valve number 7 (see attached pneumatic diagram).

# 1. HYDRAULIC & PNEUMATIC PRINCIPLES

## 2.2. Pack

### 2.2.1. Transmission tubes list

DESIGNATION	PART NUMBER	DIAMETER	LENGTH	QUANTITY
SLEEVE HPS3	DBD005A	5-9		0.5
T CONNECTOR	EAB006B	2.3		3
T CONNECTOR	EAB032A	1.5		1
TUBE CAP	EAC017A	2.5		1
TYGON TUBE 0.051"	EAE006A	1.30	140	1
TYGON TUBE 0.051"	EAE006A	1.30	350	1
TYGON TUBE 0.060"	EAE007A	1.52	20	1
TYGON TUBE 0.060"	EAE007A	1.52	40	3
TYGON TUBE 0.060"	EAE007A	1.52	50	1
TYGON TUBE 0.060"	EAE007A	1.52	60	1
TYGON TUBE 0.060"	EAE007A	1.52	70	2
TYGON TUBE 0.060"	EAE007A	1.52	80	1
TYGON TUBE 0.060"	EAE007A	1.52	150	1
TYGON TUBE 0.060"	EAE007A	1.52	170	1
TYGON TUBE 0.060"	EAE007A	1.52	220	1
TYGON TUBE 0.060"	EAE007A	1.52	240	1
TYGON TUBE 0.060"	EAE007A	1.52	300	1
TYGON TUBE 0.060"	EAE007A	1.52	370	1
TYGON TUBE 0.060"	EAE007A	1.52	410	1
TYGON TUBE 0.060"	EAE007A	1.52	450	2
TYGON TUBE 0.060"	EAE007A	1.52	480	1
TYGON TUBE 0.081"	EAE008A	2.05	20	1
TYGON TUBE 0.081"	EAE008A	2.05	35	1
TYGON TUBE 0.081"	EAE008A	2.05	200	1
TYGON TUBE 0.081"	EAE008A	2.05	590	1
TYGON TUBE 0.090"	EAE009A	2.28	20	2
TYGON TUBE 0.090"	EAE009A	2.28	50	1
TYGON TUBE 0.090"	EAE009A	2.28	60	2
TYGON TUBE 0.090"	EAE009A	2.28	120	1
TYGON TUBE 0.090"	EAE009A	2.28	150	1
TYGON TUBE 0.090"	EAE009A	2.28	190	1
TYGON TUBE 0.090"	EAE009A	2.28	510	1
TYGON TUBE 0.090"	EAE009A	2.28	550	1
TYGON TUBE 0.090"	EAE009A	2.28	1100	1
SLEEVE	GAL098A			31
TUBE SHIELD	GBC088A	4.4	30	1
TUBE SHIELD	GBC088A	4.4	60	1
METALLIC SHEATH	GBC170A	5.2	35	3



# 1. HYDRAULIC & PNEUMATIC PRINCIPLES

## 2.2.2. Hydropneumatic connections

CIRCUIT	FROM	SLEEVE	DIAMETER	LENGTH	TO	SLEEVE
AIR	(atmosphere)		2.28	190	Liquid valve 2_2	
	Liquid valve 2_1	Y	2.28	50	Waste-chamb._1	Y
	Waste-chamb._2		2.28	20	cap EAC017A	
DILUENT	Pack_3 (Diluent)		2.28	550	Liquid valve 11_1	Y
	Pack_3 (Diluent)		gbc170a			
	Liquid valve 11_3	Y	1.52	220	temp. sensor xba281a	Y
	temp. sensor xba281a	Y	1.52	40	Liquid syringes _3	Y
	Liquid valve 11_2	Y	1.52	40	Liquid valve 10_3	Y
	Liquid valve 10_1	Y	1.52	70	Liquid valve 7_3	Y
	Liquid valve 7_1	Y	1.52	370	Needle rinsing block_2	
	needle rinsing block._1		1.52	410	Liquid valve 8_1	
	Liquid valve 8_2	Y	1.52	300	Waste-chamb._3	
	Liquid valve 7_2	Y	1.52	240	T connector 2.3_1	
	T connector 2.3_2		1.52	50	Bac WBC/HGB_3	
	Liquid valve 10_2	Y	1.30	140	Liquid syringes _1	Y
	Liquid syringes _2	Y	1.30	350	NEEDLE_1	
	Liquid syringes _2	Y			NEEDLE_1	Y
CLEAN	Pack-1 (Clean)		2.28	510	Liquid valve 4_2	Y
	Pack-1 (Clean)		gbc170a			
	Liquid valve 4_1	Y	1.52	450	T connector 2.3_3	
LYSE	Pack_2 (Lyse)		2.05	590	Liquid valve 1_1	
	Pack_2 (Lyse)		gbc170a			
	Liquid valve 1_3		1.52	150	Liquid syringes _4	
	Liquid valve 1_2		1.52	480	WBC grounding connector	
WBC/RBC counting	WBC/HGB chamber _2		1.52	170	RBC chamber _3	
	RBC chamber _2		1.52	450	Liquid valve 6_2	Y
	Liquid valve 6_1	Y	1.52	60	Waste-chamb._4	
DRAIN / BUBBLING	WBC/HGB chamber _1		1.52	40	T connector 1.5_2	
	WBC/HGB chamber _1		gbc088a	30	T connector 1.5_2	
	T connector 1.5_3		1.52	80	Insulator WBC_1	
	Insulator WBC_2		2.28	120	Liquid valve 12_2	Y
	Liquid valve 12_1	Y	2.05	35	T connector 2.3_1	
	Bac RBC_1		1.52	70	Insulator RBC_1	
	Bac RBC_1		gbc088a	60	Insulator RBC_1	
	Insulator RBC_2		2.28	150	Liquid valve 13_2	Y
	Liquid valve 13_1	Y	2.05	20	T connector 2.3_2	
	T connector 2.3_3		2.05	200	T connector 2.3_2	
	T connector 2.3_3		2.28	20	Cellule xba199a	Y
	Cell xba199a	Y	2.28	60	Waste-chamb._5	Y
	T connector 2.3_1		2.28	60	Liquid valve 5_2	Y
	Liquid valve 5_1		2.28	1100	Pack_4 (Waste)	

### NOTE

Read this table as follows in this example :

The Liquid valve 7\_2 corresponds to output 2 of the valve number 7 (see attached pneumatic diagram).

# 1. HYDRAULIC & PNEUMATIC PRINCIPLES

## 2.3. Hydraulic cycle description

*2.3.1. Atmosphere circuit*

*2.3.2. Diluent circuit*



*2.3.3. Clean circuit*

*2.3.4. Lyse circuit*

*2.3.5. WBC/RBC counting circuit*

*2.3.6. Drain/bubbling circuit*

ATMOSPHERE CIRCUIT





Translation inch to mm (Ø int.)			
0.010"	0.25mm	0.045"	1.14mm
0.020"	0.50mm	0.051"	1.30mm
0.025"	0.64mm	0.060"	1.52mm
0.030"	0.76mm	0.085"	2.05mm
0.035"	0.89mm	0.095"	2.28mm
0.040"	1.02mm	0.100"	2.54mm

- C: tube Cristal
- P: tube Polyuréthane
- : tube Tygon
- S: tube Silicone
- V: tube Viton

DILUENT CIRCUIT



- - - - - COUNTING HEAD RINSE CIRCUIT
- - - - - OUTER NEEDLE RINSE CIRCUIT
- - - - - INNER NEEDLE CIRCUIT



Translation inch to mm [Ø int.]			
0.010"	0.25mm	0.045"	1.14mm
0.020"	0.50mm	0.051"	1.30mm
0.025"	0.64mm	0.060"	1.52mm
0.030"	0.76mm	0.081"	2.05mm
0.035"	0.89mm	0.090"	2.28mm
0.040"	1.02mm	0.100"	2.54mm

- C: tube Cristal
- P: tube Polyuréthane
- : tube Tygon
- S: tube Silicone
- V: tube Viton



CLEAN CIRCUIT



Translation inch to mm (Ø int.)			
0.010"	0.25mm	0.045"	1.14mm
0.020"	0.50mm	0.051"	1.30mm
0.026"	0.64mm	0.060"	1.52mm
0.030"	0.76mm	0.085"	2.05mm
0.035"	0.89mm	0.095"	2.26mm
0.040"	1.02mm	0.100"	2.54mm

- C: tube Cristal
- P: tube Polyuréthane
- : tube Tygon
- S: tube Silicone
- V: tube Viton



LYSE CIRCUIT



Translation inch to mm (Ø int.)			
0.010"	0.25mm	0.045"	1.14mm
0.020"	0.50mm	0.051"	1.30mm
0.025"	0.64mm	0.060"	1.52mm
0.030"	0.76mm	0.081"	2.05mm
0.035"	0.89mm	0.090"	2.28mm
0.040"	1.02mm	0.100"	2.54mm

- C: tube Cristal
- P: tube Polyuréthane
- : tube Tygon
- S: tube Silicone
- V: tube Viton


WBC / RBC COUNTING PHASIS



Translation inch to mm (Ø int.)			
0.010"	0.25mm	0.045"	1.14mm
0.020"	0.50mm	0.051"	1.30mm
0.026"	0.64mm	0.060"	1.52mm
0.030"	0.76mm	0.081"	2.05mm
0.035"	0.89mm	0.090"	2.28mm
0.040"	1.02mm	0.100"	2.54mm

- C: tube Cristal
- P: tube Polyuréthane
- : tube Tygon
- S: tube Silicone
- V: tube Viton

DRAIN / BUBBLING CIRCUIT





# 1. HYDRAULIC & PNEUMATIC PRINCIPLES

## 3. MICROS 60 CT HYDRAULIC

### 3.1. With bottles

#### 3.1.1. Transmission tubes list

DESIGNATION	PART NUMBER	DIAMETER	LENGTH	QUANTITY
T CONNECTOR	EAB006B	2.3		4
STRAIGHT CONNECTOR	EAB015A	1.5/2.5		3
T CONNECTOR	EAB032A	1.5		2
TUBE CAP	EAC017A	2.5		1
TYGON TUBE 0.040"	EAE005A	1.02	110	1
TYGON TUBE 0.040"	EAE005A	1.02	335	1
TYGON TUBE 0.060"	EAE007A	1.52	15	3
TYGON TUBE 0.060"	EAE007A	1.52	20	2
TYGON TUBE 0.060"	EAE007A	1.52	40	4
TYGON TUBE 0.060"	EAE007A	1.52	50	2
TYGON TUBE 0.060"	EAE007A	1.52	60	1
TYGON TUBE 0.060"	EAE007A	1.52	70	2
TYGON TUBE 0.060"	EAE007A	1.52	80	1
TYGON TUBE 0.060"	EAE007A	1.52	100	1
TYGON TUBE 0.060"	EAE007A	1.52	150	1
TYGON TUBE 0.060"	EAE007A	1.52	170	1
TYGON TUBE 0.060"	EAE007A	1.52	220	1
TYGON TUBE 0.060"	EAE007A	1.52	240	1
TYGON TUBE 0.060"	EAE007A	1.52	420	2
TYGON TUBE 0.060"	EAE007A	1.52	450	2
TYGON TUBE 0.060"	EAE007A	1.52	480	1
TYGON TUBE 0.081"	EAE008A	2.05	20	1
TYGON TUBE 0.081"	EAE008A	2.05	35	1
TYGON TUBE 0.081"	EAE008A	2.05	200	1
TYGON TUBE 0.081"	EAE008A	2.05	330	1
TYGON TUBE 0.081"	EAE008A	2.05	650	1
TYGON TUBE 0.081"	EAE008A	2.05	1080	1
TYGON TUBE 0.090"	EAE009A	2.28	20	2
TYGON TUBE 0.090"	EAE009A	2.28	50	1
TYGON TUBE 0.090"	EAE009A	2.28	60	3
TYGON TUBE 0.090"	EAE009A	2.28	120	1
TYGON TUBE 0.090"	EAE009A	2.28	140	1
TYGON TUBE 0.090"	EAE009A	2.28	150	1
TYGON TUBE 0.090"	EAE009A	2.28	190	1
SILICONE TUBE	EAE025A	1.5/3.5	50	2
BLUE TYGON TUBE 0.090"	EAE036A	2.28	1100	1
SLEEVE	GAL098A			32
TUBE SHIELD	GBC088A	4.4	30	1
TUBE SHIELD	GBC088A	4.4	60	1

# 1. HYDRAULIC & PNEUMATIC PRINCIPLES

## 3.1.2. Hydropneumatic connections

CIRCUIT	FROM	SLEEVE	DIAMETER	LENGTH	TO	SLEEVE
AIR	(atmosphere)		2.28	190	Liquid valve 2_2	
	Liquid valve 2_1	Y	2.28	50	Waste-chamb._1	Y
	Waste-chamb._2		2.28	20	cap EAC017A	
	(atmosphere)		2.28	60	Liquid valve 3_2	
	Liquid valve 3_1		1.52	40	T connector 2.3_1	
DILUENT	Diluent input		2.05	330	Liquid valve 11_1	
	Liquid valve 11_3	Y	1.52	220	temp sensor. xba281a	Y
	temp sensor. xba281a	Y	1.52	40	liquid syringes_3	Y
	Liquid valve 11_2	Y	1.52	40	Liquid valve 10_3	Y
	Liquid valve 10_1	Y	1.52	70	Liquid valve 7_3	Y
	Liquid valve 7_1	Y	1.52	50	Liquid valve 9_3	Y
	Liquid valve 9_1	Y	1.52	420	Needle rinsing block._2	
	Liquid valve 9_2	Y	1.52	420	T connector 1.5_2	
	T connector 1.5_1		S1.5/3.5	50	Needle rinsing block._3	
	T connector 1.5_3		S1.5/3.5	50	Needle rinsing block._4	
	Liquid valve 7_2	Y	1.52	240	T connector 2.3_1	
	T connector 2.3_2		1.52	50	WBC/HGB chamber_3	
	Liquid valve 10_2	Y	1.52	15	Connector 1.5/2.5	
	Connector 1.5/2.5		1.02	110	Connector 1.5/2.5	
	Connector 1.5/2.5		1.52	15	liquid syringes_1	Y
	liquid syringes_2	Y	1.52	15	Connector 1.5/2.5	
	Connector 1.5/2.5		1.02	335	needle_1	
	Connector 1.5/2.5	Y			needle_1	Y
	Needle rinsing block._1		1.52	100	T connector 2.3_3	
	T connector 2.3_2		1.52	20	Liquid valve 8_1	
Liquid valve 8_2		2.05	650	Waste-chamb._3		
CLEAN	Clean bottle		2.28	1100 blue	Liquid valve 4_2	Y
	Liquid valve 4_1	Y	1.52	450	T connector 2.3_3	
LYSE	Lyse bottle		2.05	1080	Liquid valve 1_1	
	Liquid valve 1_3		1.52	150	liquid syringes_4	
	Liquid valve 1_2		1.52	480	WBC grounding connector	
	WBC grounding connector		1.52	20	T connector 1.5_1	
WBC/RBC counting	WBC/HGB chamber_2		1.52	170	RBC chamber_3	
	RBC chamber_2		1.52	450	Liquid valve 6_2	Y
	Liquid valve 6_1	Y	1.52	60	Waste-chamb._4	
DRAIN / BUBBLING	WBC/HGB Chamber_1		1.52	40	T connector 1.5_2	
	WBC/HGB Chamber_1		gbc088a	30	T connector 1.5_2	
	T connector 1.5_3		1.52	80	insulator WBC_1	
	insulator WBC_2		2.28	120	Liquid valve 12_2	Y
	Liquid valve 12_1	Y	2.05	35	T connector 2.3_1	
	RBC_1 chamber		1.52	70	insulator RBC_1	
	RBC_1 chamber		gbc088a	60	insulator RBC_1	
	insulator RBC_2		2.28	150	Liquid valve 13_2	Y
	Liquid valve 13_1	Y	2.05	20	T connector 2.3_2	
	T connector 2.3_3		2.05	200	T connector 2.3_2	
	T connector 2.3_3		2.28	20	Cell xba199a	Y
	Cell xba199a	Y	2.28	60	Waste-chamb._5	Y
	T connector 2.3_1		2.28	60	Liquid valve 5_2	Y
Liquid valve 5_1		2.28	140	Waste ouput		

### NOTE

Read this table as follows in this example :  
The Liquid valve 7\_2 corresponds to output 2 of the valve number 7 (see attached pneumatic diagram).

# 1. HYDRAULIC & PNEUMATIC PRINCIPLES

## 3.2. Pack

### 3.2.1. Transmission tubes list

DESIGNATION	PART NUMBER	DIAMETER	LENGTH	QUANTITY
T CONNECTOR	EAB006B	2.3		4
STRAIGHT CONNECTOR	EAB015A	1.5/2.5		3
T CONNECTOR	EAB032A	1.5		2
TUBE CAP	EAC017A	2.5		1
TYGON TUBE 0.040"	EAE005A	1.02	110	1
TYGON TUBE 0.040"	EAE005A	1.02	335	1
TYGON TUBE 0.060"	EAE007A	1.52	15	3
TYGON TUBE 0.060"	EAE007A	1.52	20	2
TYGON TUBE 0.060"	EAE007A	1.52	40	4
TYGON TUBE 0.060"	EAE007A	1.52	50	2
TYGON TUBE 0.060"	EAE007A	1.52	60	1
TYGON TUBE 0.060"	EAE007A	1.52	70	2
TYGON TUBE 0.060"	EAE007A	1.52	80	1
TYGON TUBE 0.060"	EAE007A	1.52	100	1
TYGON TUBE 0.060"	EAE007A	1.52	150	1
TYGON TUBE 0.060"	EAE007A	1.52	170	1
TYGON TUBE 0.060"	EAE007A	1.52	220	1
TYGON TUBE 0.060"	EAE007A	1.52	240	1
TYGON TUBE 0.060"	EAE007A	1.52	420	2
TYGON TUBE 0.060"	EAE007A	1.52	450	2
TYGON TUBE 0.060"	EAE007A	1.52	480	1
TYGON TUBE 0.081"	EAE008A	2.05	20	1
TYGON TUBE 0.081"	EAE008A	2.05	35	1
TYGON TUBE 0.081"	EAE008A	2.05	200	1
TYGON TUBE 0.081"	EAE008A	2.05	590	1
TYGON TUBE 0.081"	EAE008A	2.05	650	1
TYGON TUBE 0.090"	EAE009A	2.28	20	2
TYGON TUBE 0.090"	EAE009A	2.28	50	1
TYGON TUBE 0.090"	EAE009A	2.28	60	3
TYGON TUBE 0.090"	EAE009A	2.28	120	1
TYGON TUBE 0.090"	EAE009A	2.28	150	1
TYGON TUBE 0.090"	EAE009A	2.28	190	1
TYGON TUBE 0.090"	EAE009A	2.28	510	1
TYGON TUBE 0.090"	EAE009A	2.28	550	1
TYGON TUBE 0.090"	EAE009A	2.28	1100	1
SILICONE TUBE	EAE025A	1.5/3.5	50	2
SLEEVE	GAL098A			33
TUBE SHIELD	GBC088A	4.4	30	1
TUBE SHIELD	GBC088A	4.4	60	1
METALLIC SHEATH	GBC170A	5.2	35	3

# 1. HYDRAULIC & PNEUMATIC PRINCIPLES

## 3.2.2. Hydropneumatic connections

CIRCUIT	FROM	SLEEVE	DIAMETER	LENGTH	TO	SLEEVE
AIR	(atmosphere)		2.28	190	Liquid valve 2_2	
	Liquid valve 2_1	Y	2.28	50	Waste-chamb._1	Y
	Waste-chamb._2		2.28	20	cap EAC017A	
	(atmosphere)		2.28	60	Liquid valve 3_2	
	Liquid valve 3_1		1.52	40	T Connector 2.3_1	
DILUENT	Pack_3 (Diluent)		2.28	550	Liquid valve 11_1	Y
	Pack_3 (Diluent)		gbc170a			
	Liquid valve 11_3	Y	1.52	220	Temp. sensor xba281a	Y
	Temp. sensor xba281a	Y	1.52	40	Liquid syringes_3	Y
	Liquid valve 11_2	Y	1.52	40	Liquid valve 10_3	Y
	Liquid valve 10_1	Y	1.52	70	Liquid valve 7_3	Y
	Liquid valve 7_1	Y	1.52	50	Liquid valve 9_3	Y
	Liquid valve 9_1	Y	1.52	420	Needle rinsing block_2	
	Liquid valve 9_2	Y	1.52	420	T Connector 1.5_2	
	T Connector 1.5_1		S1.5/3.5	50	Needle rinsing block_3	
	T Connector 1.5_3		S1.5/3.5	50	Needle rinsing block_4	
	Liquid valve 7_2	Y	1.52	240	T Connector 2.3_1	
	T Connector 2.3_2		1.52	50	Bac WBC/HGB_3	
	Liquid valve 10_2	Y	1.52	15	connector1.5/2.5	
	connector1.5/2.5		1.02	110	connector1.5/2.5	
	connector 1.5/2.5		1.52	15	Liquid syringes_1	Y
	Liquid syringes_2	Y	1.52	15	connector 1.5/2.5	
	connector 1.5/2.5		1.02	335	Needle_1	
	connector 1.5/2.5	Y			Needle_1	Y
	Needle rinsing block_1		1.52	100	T Connector 2.3_3	
T Connector 2.3_2		1.52	20	Liquid valve 8_1		
Liquid valve 8_2		2.05	650	Waste-chamb._3		
CLEAN	Pack_1 (Clean)		2.28	510	Liquid valve 4_2	Y
	Pack_1 (Clean)		gbc170a			
	Liquid valve 4_1	Y	1.52	450	T Connector 2.3_3	
LYSE	Pack_2 (Lyse)		2.05	590	Liquid valve 1_1	
	Pack_2 (Lyse)		gbc170a			
	Liquid valve 1_3		1.52	150	Liquid syringes_4	
	Liquid valve 1_2		1.52	480	WBCgrounding connec.	
Comptage WBC/RBC	WBC grounding connect.		1.52	20	T Connector 1.5_1	
	Bac WBC/HGB_2		1.52	170	RBC chamber_3	
	RBC chamber_2		1.52	450	Liquid valve 6_2	Y
EVACUAT <sup>9</sup> / Bullage	Liquid valve 6_1	Y	1.52	60	Waste-chamb._4	
	Bac WBC/HGB_1		1.52	40	T Connector 1.5_2	
	Bac WBC/HGB_1		gbc088a	30	T Connector 1.5_2	
	T Connector 1.5_3		1.52	80	insulator WBC_1	
	insulator WBC_2		2.28	120	Liquid valve 12_2	Y
	Liquid valve 12_1	Y	2.05	35	T Connector 2.3_1	
	RBC chamber_1		1.52	70	insulator RBC_1	
	RBC chamber_1		gbc088a	60	insulator RBC_1	
	insulator RBC_2		2.28	150	Liquid valve 13_2	Y
	Liquid valve 13_1	Y	2.05	20	T Connector 2.3_2	
	T Connector 2.3_3		2.05	200	T Connector 2.3_2	
	T Connector 2.3_3		2.28	20	CELL xba199a	Y
	CELL xba199a	Y	2.28	60	Waste-chamb._5	Y
T Connector 2.3_1		2.28	60	Liquid valve 5_2	Y	
Liquid valve 5_1		2.28	1100	Pack_4 (Waste)		

### NOTE

Read this table as follows in this example :

The Liquid valve 7\_2 corresponds to output 2 of the valve number 7 (see attached pneumatic diagram)

# 1. HYDRAULIC & PNEUMATIC PRINCIPLES

## 3.3. Hydraulic cycle description

*3.3.1. Atmosphere circuit*

*3.3.2. Diluent circuit*



*3.3.3. Clean circuit*

*3.3.4. Lyse circuit*

*3.3.5. WBC/RBC counting circuit*

*3.3.6. Drain/bubbling circuit*

ATMOSPHERE CIRCUIT





Translation inch to mm (Ø int.)			
0.010"	0.25mm	0.045"	1.14mm
0.020"	0.50mm	0.051"	1.30mm
0.025"	0.64mm	0.060"	1.52mm
0.030"	0.75mm	0.068"	1.73mm
0.035"	0.89mm	0.090"	2.28mm
0.040"	1.02mm	0.100"	2.54mm

- C: tube Cristal
- P: tube Polyuréthane
- T: tube Tygon
- S: tube Silicone
- V: tube Viton

DILUENT CIRCUIT



- - - - - COUNTING HEAD RINSE CIRCUIT
- - - - - OUTER SAMPLING NEEDLE RINSE CIRCUIT
- - - - - SAMPLING NEEDLE INNER CIRCUIT
- - - - - PIERCING NEEDLE RINSE CIRCUIT



Translation inch to mm (Ø int.)			
0.010"	0.25 mm	0.046"	1.14 mm
0.020"	0.50 mm	0.051"	1.30 mm
0.025"	0.64 mm	0.060"	1.52 mm
0.030"	0.76 mm	0.084"	2.05 mm
0.035"	0.89 mm	0.090"	2.28 mm
0.040"	1.02 mm	0.100"	2.54 mm

- C: tube Cristal
- P: tube Polyuréthane
- : tube Tygon
- S: tube Silicone
- V: tube Viton

CLEAN CIRCUIT




Translation inch to mm (Ø int.)			
0.010"	0.25mm	0.046"	1.14mm
0.020"	0.50mm	0.051"	1.30mm
0.025"	0.64mm	0.060"	1.52mm
0.030"	0.76mm	0.081"	2.05mm
0.035"	0.89mm	0.090"	2.28mm
0.040"	1.02mm	0.100"	2.54mm


- C: tube Cristal
- P: tube Polyuréthane
- : tube Tygon
- S: tube Silicone
- V: tube Viton



LYSE CIRCUIT



XBA199A




3/2-3T








2/2 NC

Translation inch to mm (Ø int.)			
0.010"	0.25 mm	0.046"	1.14 mm
0.020"	0.50 mm	0.051"	1.30 mm
0.025"	0.64 mm	0.060"	1.52 mm
0.030"	0.76 mm	0.081"	2.05 mm
0.035"	0.89 mm	0.090"	2.28 mm
0.040"	1.02 mm	0.100"	2.54 mm

- C: tube Cristal
- P: tube Polyuréthane
- : tube Tygon
- S: tube Silicone
- V: tube Viton





WBC / RBC COUNTING CIRCUIT



Translation inch to mm (Ø int.)			
0.010"	0.25mm	0.046"	1.14mm
0.020"	0.50mm	0.051"	1.30mm
0.025"	0.64mm	0.060"	1.52mm
0.030"	0.76mm	0.081"	2.05mm
0.035"	0.89mm	0.090"	2.28mm
0.040"	1.02mm	0.100"	2.54mm

- C: tube Cristal
- P: tube Polyuréthane
- : tube Tygon
- S: tube Silicone
- V: tube Viton

DRAIN / BUBBLING CIRCUIT



Translation inch to mm (Ø int.)			
0.010"	0.25 mm	0.046"	1.14 mm
0.020"	0.50 mm	0.051"	1.30 mm
0.025"	0.64 mm	0.060"	1.52 mm
0.030"	0.76 mm	0.061"	2.05 mm
0.035"	0.89 mm	0.090"	2.28 mm
0.040"	1.02 mm	0.100"	2.54 mm

- C: tube Cristal
- P: tube Polyuréthane
- : tube Tygon
- S: tube Silicone
- V: tube Viton

# 1. HYDRAULIC & PNEUMATIC PRINCIPLES

## 4. PNEUMATIC DIAGRAMS

- 4.1. Micros 60 CT bottle version
- 4.2. Micros 60 CT pack version
- 4.3. Micros 60 OT bottle version
- 4.4. Micros 60 OT pack version

# 2. ELECTRIC & ELECTRONIC PRINCIPLES


## CONTENTS

1. POWER SUPPLY ASSEMBLY .....	2
1.1. Principle .....	2
1.2. Power supply fan .....	4
1.3. Secondary supply cable .....	5
1.4. Main supply filter .....	6
2. MOTHER BOARD .....	7
2.1. Configuration .....	7
3. LCD BOARD .....	8
4. COAXIALS .....	9
4.1. RBC coaxial .....	9
4.2. WBC coaxial .....	11
5. SENSORS .....	12
5.1. Drain detection sensor .....	12
5.2. Carriage & needle sensors .....	13
5.3. Syringe sensors .....	14
6. MISCELLANEOUS .....	15
6.1. Needle carriage motor .....	15
6.2. Barcode cable .....	16
6.3. Temperature sensor .....	17
6.4. HGB Chamber assembly .....	18
6.5. Chip card reader cable .....	19
6.6. CT Piercing assembly cable .....	20
6.7. CT twin valve cable .....	21
7. FLAT CABLES .....	22
8. ELECTRICAL SYNOPTICS .....	26

# 2. ELECTRIC & ELECTRONIC PRINCIPLES


## 1. POWER SUPPLY ASSEMBLY

### 1.1. Principle



Désignation: Schéma de principe		XBA274D
ENSEMBLE ALIM* MICROS_60 JAPON Ind Modif:A		
REF. du document: XBA274DS1	Nbr.:	
Plan dessiné le: 10/02/98 Par: A. DA COSTA		
HEMATOLOGIE		
Ce plan est la propriété de la société ABX et ne peut être reproduit ou communiqué sans autorisation		
In. DM	Date	
A 3781	10/02/98	CREATION POUR JAPON: NORME IEC 601-1/IEC 1000-4-5 (IMMUNITÉ AUX ONDES DE CHOC: Foudre)
Modification		

# 2. ELECTRIC & ELECTRONIC PRINCIPLES




Designation: Schéma de définition XBA274D  
 ENSEMBLE ALIM. MICROS 80 JAPON Inc Modif: A  
 REF DU DOCUMENT: XBA274DSE2  
 Nbr:   
 Plan Gess.1n6 16: 10/02/98 Per: A. DA COSTA  
 Ce plan est le propriété de la société AXI et ne peut être reproduit ni communiqué sans autorisation

A 3781 10/02/98 CREATION POUR JAPAN NOME IEC 801-2/IEC 1000-4-5 (IMMEDIATE AUX ONDES DE CHOC. FOURME)  
 In Date MODIFICATION

# 2. ELECTRIC & ELECTRONIC PRINCIPLES

## 1.2. Power supply fan





# 2. ELECTRIC & ELECTRONIC PRINCIPLES

## 1.3. Secondary supply cable

**VERSION CARTE PU65**

**J2-PU65 SWITCHING SUPPLY BOARD**

**J26-MASTER BOARD**

P26

1 2 3 4 5 6 7 8

+5V +5V +5V GND +24V +24V +12V -12V

18 RED 18 RED 18 RED 18 GREEN 18 BLUE 18 BLUE 18 VIOLET 18 BROWN

48 MEGAMIDE SLEEVE

5195-8 (REF: 09-50-1081) CAK024A

CONTACTS: 5194 (REF: 09-70-1030) CAK025A

MANCHON HP53 COLLIER 2-4#32

450 mm

400 mm

P26

1

**VERSION CARTE US50**

**J2-US50 SWITCHING SUPPLY BOARD**

**J26-MASTER BOARD**

P26

1 2 3 4 5 6 7 8

+5V +5V +5V GND +24V +24V +12V -12V

18 RED 18 RED 18 RED 18 GREEN 18 BLUE 18 BLUE 18 VIOLET 18 BROWN

48 MEGAMIDE SLEEVE

5195-8 (REF: 09-50-1081) CAK024A

CONTACTS: 5194 (REF: 09-70-1030) CAK025A

MANCHON HP53 COLLIER 2-4#32

450 mm

400 mm

P26

1

NOTA : LES FILS DE MEME FONCTION (TENSION, JAUGE, COULEUR) PEUVENT ETRE INTERVERTIS.  
PLACER UN DETROMPEUR EN 5 DE P2.

Désignation: TORON ALIMENTATION SECONDAIRE		XBA276A
MICROS		Ind Modif: B
REF du document: XBA276AS		Nbr:
Plan dessiné le: 01/01/94		Par: A. DA COSTA
Ce plan est la propriété de la société Roche ABX et ne peut être reproduit ou communiqué sans autorisation		HEMATOLOGIE

**VERSION CARTE PU65**

**J2-PU65 SWITCHING SUPPLY BOARD**

**J26-MASTER BOARD**

P26

1 2 3 4 5 6 7 8

+5V +5V +5V GND +24V +24V +12V -12V

18 RED 18 RED 18 RED 18 GREEN 18 BLUE 18 BLUE 18 VIOLET 18 BROWN

48 MEGAMIDE SLEEVE

5195-8 (REF: 09-50-1081) CAK024A

CONTACTS: 5194 (REF: 09-70-1030) CAK025A

MANCHON HP53 COLLIER 2-4#32

450 mm

400 mm

P26

1

**VERSION CARTE US50**

**J2-US50 SWITCHING SUPPLY BOARD**

**J26-MASTER BOARD**

P26

1 2 3 4 5 6 7 8

+5V +5V +5V GND +24V +24V +12V -12V

18 RED 18 RED 18 RED 18 GREEN 18 BLUE 18 BLUE 18 VIOLET 18 BROWN

48 MEGAMIDE SLEEVE

5195-8 (REF: 09-50-1081) CAK024A

CONTACTS: 5194 (REF: 09-70-1030) CAK025A

MANCHON HP53 COLLIER 2-4#32

450 mm

400 mm

P26

1

NOTA : LES FILS DE MEME FONCTION (TENSION, JAUGE, COULEUR) PEUVENT ETRE INTERVERTIS.  
PLACER UN DETROMPEUR EN 5 DE P2.

Désignation: TORON ALIMENTATION SECONDAIRE		XBA276A
MICROS		Ind Modif: B
REF du document: XBA276AS		Nbr:
Plan dessiné le: 01/01/94		Par: A. DA COSTA
Ce plan est la propriété de la société Roche ABX et ne peut être reproduit ou communiqué sans autorisation		HEMATOLOGIE

**VERSION CARTE PU65**

**J2-PU65 SWITCHING SUPPLY BOARD**

**J26-MASTER BOARD**

P26

1 2 3 4 5 6 7 8

+5V +5V +5V GND +24V +24V +12V -12V

18 RED 18 RED 18 RED 18 GREEN 18 BLUE 18 BLUE 18 VIOLET 18 BROWN

48 MEGAMIDE SLEEVE

5195-8 (REF: 09-50-1081) CAK024A

CONTACTS: 5194 (REF: 09-70-1030) CAK025A

MANCHON HP53 COLLIER 2-4#32

450 mm

400 mm

P26

1

**VERSION CARTE US50**

**J2-US50 SWITCHING SUPPLY BOARD**

**J26-MASTER BOARD**

P26

1 2 3 4 5 6 7 8

+5V +5V +5V GND +24V +24V +12V -12V

18 RED 18 RED 18 RED 18 GREEN 18 BLUE 18 BLUE 18 VIOLET 18 BROWN

48 MEGAMIDE SLEEVE

5195-8 (REF: 09-50-1081) CAK024A

CONTACTS: 5194 (REF: 09-70-1030) CAK025A

MANCHON HP53 COLLIER 2-4#32

450 mm

400 mm

P26

1

NOTA : LES FILS DE MEME FONCTION (TENSION, JAUGE, COULEUR) PEUVENT ETRE INTERVERTIS.  
PLACER UN DETROMPEUR EN 5 DE P2.

Désignation: TORON ALIMENTATION SECONDAIRE		XBA276A
MICROS		Ind Modif: B
REF du document: XBA276AS		Nbr:
Plan dessiné le: 01/01/94		Par: A. DA COSTA
Ce plan est la propriété de la société Roche ABX et ne peut être reproduit ou communiqué sans autorisation		HEMATOLOGIE

**VERSION CARTE PU65**

**J2-PU65 SWITCHING SUPPLY BOARD**

**J26-MASTER BOARD**

P26

1 2 3 4 5 6 7 8

+5V +5V +5V GND +24V +24V +12V -12V

18 RED 18 RED 18 RED 18 GREEN 18 BLUE 18 BLUE 18 VIOLET 18 BROWN

48 MEGAMIDE SLEEVE

5195-8 (REF: 09-50-1081) CAK024A

CONTACTS: 5194 (REF: 09-70-1030) CAK025A

MANCHON HP53 COLLIER 2-4#32

450 mm

400 mm

P26

1

**VERSION CARTE US50**

**J2-US50 SWITCHING SUPPLY BOARD**

**J26-MASTER BOARD**

P26

1 2 3 4 5 6 7 8

+5V +5V +5V GND +24V +24V +12V -12V

18 RED 18 RED 18 RED 18 GREEN 18 BLUE 18 BLUE 18 VIOLET 18 BROWN

48 MEGAMIDE SLEEVE

5195-8 (REF: 09-50-1081) CAK024A

CONTACTS: 5194 (REF: 09-70-1030) CAK025A

MANCHON HP53 COLLIER 2-4#32

450 mm

400 mm

P26


1

NOTA : LES FILS DE MEME FONCTION (TENSION, JAUGE, COULEUR) PEUVENT ETRE INTERVERTIS.  
PLACER UN DETROMPEUR EN 5 DE P2.

Désignation: TORON ALIMENTATION SECONDAIRE		XBA276A
MICROS		Ind Modif: B
REF du document: XBA276AS		Nbr:
Plan dessiné le: 01/01/94		Par: A. DA COSTA
Ce plan est la propriété de la société Roche ABX et ne peut être reproduit ou communiqué sans autorisation		HEMATOLOGIE

# 2. ELECTRIC & ELECTRONIC PRINCIPLES

## 1.4. Main supply filter



# 2. ELECTRIC & ELECTRONIC PRINCIPLES

## 2. MOTHER BOARD

### 2.1. Configuration

## JUMPERS FUNCTION

E3: 6BHCL1.MODB\_sional\_opt100 on: MODB = 0 (test mode, future use)  
off: MODB = 1 (default value)

E4: on: clock connected (default value)  
off: clock no connected

E6: BS232\_CONFIGURATION

E7: on: RTS/CTS are connected to DSR  
off: RTS/CTS are not connected to DSR (default value)

## CONFIGURATION

ABX Recherche & Développement  
128 Rue du Caducee, 34184 MONTPELLIER


Title: MICRO5 MASTER Board  
Size Document Number: XAA355BR  
REV: A

12/10/95  
DATE: 1995  
BY: October 18, 1995  
Sheet: 1 of 1

E11 EPROM address area opt100		AREA	
EPROM	JUMPERS	AREA	AREA
128Kx8	FSM...NC	1-3	don't
	SAIPE		francais
256Kx8	FSM...AAT	1-3	4-6
		1-3	2-4
	ALB...AAT		francais
		3-5	4-6
		3-5	2-4
		1-3	4-6
		1-3	2-4
			espagnoi

# 2. ELECTRIC & ELECTRONIC PRINCIPLES

## 3. LCD BOARD



ABX R&D  
 Rue du Caduceus, BP7290  
 34184 MONTPELLIER Cedex 4

Title Clavier, Afficheur, LEDs MICROS

Size Document Number XAA354B5

REV B


DATE	IND	DESCRIPTION	BY	Date	Sheet	of
22/11/96	B	Application DM3401	OD			
14/12/95	A	Creation suite DM2533	OD			
				November 22, 1996		

(H1) (H2) (H3) (H4)

# 2. ELECTRIC & ELECTRONIC PRINCIPLES


## 4. COAXIALS

### 4.1. RBC coaxial



Désignation: PREPARATION		XBA371A
COAX GR MICROS 60		Ind Modif: A
REF du document: XBA371AZ		Nbr:
Plan dessiné le: 15/07/97		Ech: 1
Par: D. BOISSIER		
Ce plan est la propriété de la société ABX et ne peut être reproduit ou communiqué sans autorisation		HEMATOLOGIE
A 3881	15/07/97	Création
In. DM	Date	Modification

# 2. ELECTRIC & ELECTRONIC PRINCIPLES



Désignation: CORDON COAXIAL GR  
MICROS-60  
REF du document: XBA365AS

XBA365A  
Ind Modif.: A

Nbr.:


Plan dessiné le: 06/08/97 Par: A. DA COSTA

Ce plan est la propriété de la société ABX et ne peut être reproduit ou communiqué sans autorisation

A	3632	06/08/97	Creation suivant prescrip.
In.	DM	Date	Modification

# 2. ELECTRIC & ELECTRONIC PRINCIPLES

## 4.2. WBC coaxial




Designation: CORDON COAXIAL GB		XBA272B
MICROS-60		Ind Modif: A
REF du document: XBA272BS		Nbr:
Plan dessiné le: 06/08/97		Par: A. DA COSTA
Ce plan est la propriété de la société ABX et ne peut être reproduit ou communiqué sans autorisation		HEMATOLOGIE

A 3631	07/08/97	Creation suivant preserie.
In. DW	Date	Modification

# 2. ELECTRIC & ELECTRONIC PRINCIPLES

## 5. SENSORS

### 5.1. Drain detection sensor




Ce plan est la propriété de la société ABX et ne peut être reproduit ou communiqué sans autorisation



# 2. ELECTRIC & ELECTRONIC PRINCIPLES


## 5.2. Carriage & needle sensors



Ce plan est la propriété de la société ABX et ne peut être reproduit ou communiqué sans autorisation

# 2. ELECTRIC & ELECTRONIC PRINCIPLES

## 5.3. Syringe sensors



C	3745	03/03/98	SECONDE SOURCE CELLULE IR: HONEYWELL HOA0890-T55. APPLICATION SUR LC270.
B	2636	05/05/95	Montage de la cellule sous-traité. Remplace la cellule XBA228AS sur ARGOS, HELIOS
A	2358	30/11/94	Création.

In.	DM	Date	Modification	
Désignation:		CELLULE IR 2 PATTES CABLEE	XBA319A	
		L=260mm AG/HE/VE/MS/LC220/270	Ind Modif: C	
REF du document:		XBA319AS		
			Nbr:	
Plan dessiné le:		30/11/94	Par:	A. DA COSTA




HÉMATOLOGIE

Ce plan est la propriété de la société ABX et ne peut être reproduit ou communiqué sans autorisation

# 2. ELECTRIC & ELECTRONIC PRINCIPLES

## 6. MISCELLANEOUS

### 6.1. Needle carriage motor




A	30/08/95	Création	
In.	DM	Date	Modification
Désignation: MOTEUR CHARIOT AIGUILLE CABLE		DAL008A	
MICROS		Ind Modif: A	
REF du document: DAL008AS		Nbr:	
Plan dessiné le: 30/08/95		Par: A. DA COSTA	

Ce plan est la propriété de la société Roche ABX et ne peut être reproduit ou communiqué sans autorisation

# 2. ELECTRIC & ELECTRONIC PRINCIPLES


## 6.2. Barcode cable



Ce plan est la propriété de la société ABX et ne peut être reproduit ou communiqué sans autorisation

# 2. ELECTRIC & ELECTRONIC PRINCIPLES


## 6.3. Temperature sensor



Ce plan est la propriété de la société ABX et ne peut être reproduit ou communiqué sans autorisation


# 2. ELECTRIC & ELECTRONIC PRINCIPLES

## 6.4. HGB Chamber assembly




# 2. ELECTRIC & ELECTRONIC PRINCIPLES

## 6.5. Chip card reader cable



# 2. ELECTRIC & ELECTRONIC PRINCIPLES

## 6.6. CT Piercing assembly cable




Désignation: TORON PERCUTEUR TUBE FERME XBA275A  
 Ind Modif: B  
 REF du document: XBA275AS Nbr: 1  
 Plan dessiné le: 01/01/94 Par: A. DA COSTA  
 Ce plan est la propriété de la société Roche ABX et ne peut être reproduit ou communiqué sans autorisation  
 HENATOLOGIE





# 2. ELECTRIC & ELECTRONIC PRINCIPLES

## 6.7. CT twin valve cable



# 2. ELECTRIC & ELECTRONIC PRINCIPLES

## 7. FLAT CABLES



**NOTA:** Se référer à la procédure de fabrication des limandes : RCF013A


C	11/05/94	Compatibilité entre version standard et version surmoulée.
B	08/03/94	Modification longueur de coupe.
A	03/02/94	Création.
In. DM	Date	Modification
Désignation: LIMANDE SORTIE RS232 MICROS		DAD078A Ind Modif: C
REF du document: DAD078AS		
		Nbr: 1
Plan dessiné le: 01/01/94		Par: A. DA COSTA




HÉMATOLOGIE

Ce plan est la propriété de la société Roche ABX et ne peut être reproduit ou communiqué sans autorisation

# 2. ELECTRIC & ELECTRONIC PRINCIPLES





**NOTA:** Se référer à la procédure de fabrication des limandes : RCF013A

D	2533	25/10/95	Normes CE: Allonger la limande (pour montage d'une ferrite)
C		11/05/94	Compatibilité entre version standard et version surmoulée.
B		08/03/94	Modification longueur de coupe.
A		03/02/94	Création.
In.	DM	Date	Modification
Désignation: LIMANDE SORTIE IMPRIMANTE		DAD079A	
MICROS		Ind Modif: D	
REF du document: DAD079AS		Nbr: 1	
Plan dessiné le: 01/01/94		Par: A. DA COSTA	
			 <b>HEMATOLOGIE</b>


Ce plan est la propriété de la société Roche ABX et ne peut être reproduit ou communiqué sans autorisation

# 2. ELECTRIC & ELECTRONIC PRINCIPLES




**NOTA:** Se référer à la procédure de fabrication des limandes : RCF013A

C	2584	14/04/95	Rallonger la limande de 30 mm, pour une utilisation commune sur MICROS et LC220
B		11/05/94	Compatibilité entre version standard et version surmoulée.
A		03/02/94	Création.
In.	DM	Date	Modification
Désignation: LIMANDE 5/6 ELECTROVANN MICROS 4/6 ELECTROVANNES LC220			DAD075A Ind Modif: C
REF du document: DAD075AS			Nbr:
Plan dessiné le: 01/01/94			Par: A. DA COSTA

  
**HEMATOLOGIE**

Ce plan est la propriété de la société Roche ABX et ne peut être reproduit ou communiqué sans autorisation

# 2. ELECTRIC & ELECTRONIC PRINCIPLES



# 2. ELECTRIC & ELECTRONIC PRINCIPLES

## 8. ELECTRICAL SYNOPTICS

# 3. M A I N T E N A N C E

## ***CONTENTS***

1. MAINTENANCE .....	2
1.1. Introduction .....	2
1.2. Daily customer maintenance .....	2
1.3. Weekly customer maintenance .....	2
2. Maintenance kits .....	3
3. PROCEDURES .....	4
3.1. Procedure chart tables .....	5
3.2. Required tools and products .....	7

# 3. MAINTENANCE

## 1. MAINTENANCE

### 1.1. Introduction

## **WARNING !**



Customer maintenance has to be carried out according to the recommended frequency chart table and after having performed an ABX approved customer training course.

The system warranty may be affected if damage occurs after a non trained technician intervenes or if replaced spare parts and consumables do not come from an ABX approved origin.

### 1.2. Daily customer maintenance

No special adjustments or maintenance has to be done on your equipment if the recommended startup and shutdown procedures are explicitly respected.

See the **ABX MICROS 60** User Manual for the daily rinsing and cleaning of the system.

### 1.3. Weekly customer maintenance

An overall check for cleanliness of the system is recommended every week.

All traces of blood or reagent have to be wiped off as soon as possible using a piece of cloth and distilled water.

## **CAUTION**

Never use solvent or abrasive cleaning material to clean the system.



# 3. MAINTENANCE

## 2. MAINTENANCE KITS

O R I N G   K I T   :   X E A   3 2 8   A S	
FAA 036 A	1 Lyse dispenser O ring 1 Piercing needle O ring (MICROS CT only)
FAA 029 A	1 Diluent dispenser O ring
FAA 055 A	2 Sampling syringe O rings
FAA 049 A	2 Aperture O rings
FAA 046 A	2 Coaxial O rings
FAA 054 A	2 Sampling needle O rings (CT)
FAA 053 A	1 Sampling needle O ring (OT)
FAA 017 A	1 Vacuum/waste syringe O ring
KAM 022 A	6 Board holder clips (MICROS 45 only)
FAA 023 A	1 WBC/HGB chamber cap O ring
GBC 030 A	1 Air Syringe piston

Part number	SPARE PARTS KIT X E A   4 5 8   A S	Quantity
CAE 006 A	tube holder detection switch	1
EAE 005 A	Tygon tube 1,02mm (0,040")	2
EAE 007 A	Tygon tube 1,52mm (0,060")	2
EAE 008 A	Tygon tube 2,05mm (0,081")	2
EAE 009 A	Tygon tube 2,29mm (0,090")	2
FAK 001 A	Aperture 50µ	1
FAK 003 A	Aperture 80µ	1
FBR 002 A	Notched belt 290	1
FBR 003 A	Notched belt 380	1
GBC 052 A	C.T sampling needle	1
GBC 069 A	O.T sampling needle	1
GBC 189 A	CT Piercing needle	1
GBC 193 A	RBC chamber	1
XBA 199 A	Drain detection sensor	1
XBA 250 A	Carriage/needle sensor	1
XBA 272 B	WBC coaxial	1
XDA 472 B	WBC/HGB chamber	1
XDA 481 B	Liquid valve 2W NC without solenoid	1
XDA 483 B	Liquid valve 3W without solenoid	1
XBA 365 A	RBC coaxial	1
XBA 319 A	Liquid/air syringe sensor	1
XEA 328 AS	Maintenance kit	1

# 3. MAINTENANCE

## 3. PROCEDURES

### **NOTE**

Maintenance and adjustments that need to be done on **ABX MICROS 60** are divided into «procedures» according to concerned assemblies. This should make any updating easier as all interventions can be done with the corresponding «procedure» on its own.

#### CONCERNING PARTS :

- Hydraulic maintenance and adjustments.
- Pneumatic maintenance and adjustments.
- Electrical maintenance and adjustments.
- Power supply maintenance and adjustments.
- Electronic maintenance and adjustments.
- Printer maintenance and setup.

### **NOTE**

Each procedure has to be read entirely before beginning the intervention.

## **WARNING !**



When cleaning instruments, disposable gloves should be worn.

# 3. MAINTENANCE

## 3.1. Procedure chart tables

<b>M I C R O S 6 0 O T</b>		
<b>P / N</b>	<b>PROCEDURE</b>	<b>CONCERNS</b>
RAS 165 A	<b>MICROS 60 OT</b> INSTALLATION AND STARTUP	Unpacking - Working conditions - Visual checks - Reagent connection - Printer & instrument connections - Priming & Startup
RAS 168 A	SAMPLING NEEDLE MAINTENANCE <b>MICROS OT</b>	Needle replacement - O ring replacement
RAS 169 A	CHAMBER MAINTENANCE	RBC, WBC/HGB chamber cleaning - Aperture O ring replacement - Coaxial O ring replacement
RAS 170 A	LIQUID VALVE MAINTENANCE	Liquid valve assy replacement - Valve body replacement
RAS 171 A	POWER SUPPLY CHECK/REPLACEMENT	Voltage supply check - Power supply module replacement Fan operation check.
RAS 172 A	TECHNICIAN FUNCTIONS <b>MICROS 60 OT</b>	Version display. - Adjustments : HGB photometer calibration, Aperture voltage, pressure check, WBC gain , RBC & PLT gain, Sensor, Needle height and motion, bubbling. - Temperature sensor adjustment - Run mode - Reagent pack - Serial number - Cycle number - Burning
RAS 173 A	MECHANIC FUNCTIONS	Sensor replacement - Needle motion check - Carriage motion check - Liquid syringe motion check - Vac/Waste syringe motion check - Valve operation check - LCD contrast - Piercing mechanism check
RAS 174 A	DRAIN DETECTION ADJUSTMENT	Drain detection sensor adjustment
RAS 175 A	PCB VOLTAGE CHECKS	Voltage supply check/adj. - Aperture voltage check/adj. - RBC threshold check/adj. - WBC threshold check/adj. - PLT threshold check/adj. - HGB blank voltage check - Stepper motor voltage adjustment
RAS 177 A	LX300 PRINTER	Configuration - Control panel - Control LEDs and keys - Description
RAS 178 A	LIQUID SYRINGE MAINTENANCE	Lyse dispenser O ring replacement - Diluent dispenser O ring replacement - Sampling needle dispenser O ring replact. - Lubrification
RAS 179 A	VACUUM/WASTE SYRINGE MAINTENANCE	O ring replacement
RAS 180 A	CHANGING THE INSTRUMENT LANGUAGE	Changing the instrument language
RAS 181 A	REAGENT PACK	Connector O ring replacement
RAS 182 A	BARCODE READER SETUP	Reader configuration
RAS 187 A	HYDRAULIC CYCLE CHECKUP <b>MICROS 60 OT</b>	Step by step control of the hydraulic cycle
RAS 191 A	OVERALL MAINTENANCE	Instrument maintenance step by step
RAS 192 A	DECONTAMINATION	Instrument decontamination
RAS 197 A	DRAIN & RINSE	Instrument rinse and drain for an extended shutdown

# 3. MAINTENANCE

<b>M I C R O S 6 0 C T</b>		
<b>P / N</b>	<b>PROCEDURE</b>	<b>CONCERNS</b>
RAS 166 A	<b>MICROS 60 CT</b> INSTALLATION AND STARTUP	Unpacking - Working conditions - Visual checks - Reagent connection - Printer & instrument connections - Priming & Startup
RAS 167 A	SAMPLING NEEDLE MAINTENANCE <b>MICROS CT</b>	Needle O ring replacement - Sampling needle replacement Piercing needle replacement
RAS 169 A	CHAMBER MAINTENANCE	RBC, WBC/HGB chamber cleaning - Aperture O ring replacement - Coaxial O ring replacement
RAS 170 A	LIQUID VALVE MAINTENANCE	Liquid valve assy replacement - Valve body replacement
RAS 171 A	POWER SUPPLY CHECK/REPLACEMENT	Voltage supply check - Power supply module replacement Fan operation check.
RAS 176 A	TECHNICIAN FUNCTIONS <b>MICROS 60 CT</b>	Version display. - Adjustments : HGB photometer calibration, Aperture voltage, pressure check, WBC gain , RBC & PLT gain, Sensor, Needle height and motion, bubbling. - Temperature sensor adjustment - Run mode - Reagent pack - Serial number - Cycle number - Burning
RAS 173 A	MECHANIC FUNCTIONS	Sensor replacement - Needle motion check - Carriage motion check - Liquid syringe motion check - Vac/Waste syringe motion check - Valve operation check - LCD contrast - Piercing mechanism check
RAS 174 A	DRAIN DETECTION ADJUSTMENT	Drain detection sensor adjustment
RAS 175 A	PCB VOLTAGE CHECKS	Voltage supply check/adj. - Aperture voltage check/adj. - RBC threshold check/adj. - WBC threshold check/adj. - PLT threshold check/adj. - HGB blank voltage check - Stepper motor voltage adjustment
RAS 177 A	LX300 PRINTER	Configuration - Control panel - Control LEDs and keys - Description
RAS 178 A	LIQUID SYRINGE MAINTENANCE	Lyse dispenser O ring replacement - Diluent dispenser O ring replacement - Sampling needle dispenser O ring replact. - Lubrification
RAS 179 A	VACUUM/WASTE SYRINGE MAINTENANCE	O ring replacement
RAS 180 A	CHANGING THE INSTRUMENT LANGUAGE	Changing the instrument language
RAS 181 A	REAGENT PACK	Connector O ring replacement
RAS 182 A	BARCODE READER SETUP	Reader configuration
RAS 188 A	HYDRAULIC CYCLE CHECKUP <b>MICROS 60 CT</b>	Step by step control of the hydraulic cycle
RAS 191 A	OVERALL MAINTENANCE	Instrument maintenance step by step
RAS 192 A	DECONTAMINATION	Instrument decontamination
RAS 197 A	DRAIN & RINSE	Instrument rinse and drain for an extended shutdown
RAS 198 A	PIERCING BLOCK	Description - Maintenance - Sampling position

# 3. MAINTENANCE

## 3.2. Required tools and products

T O O L S		P R O D U C T S	
DESIGNATION	PART NUMBER	DESIGNATION	PART NUMBER
HEXAGONAL KEYS		EMPTY SAMPLE TUBES	
DYNAMOMETRIC SCREW DRIVER A302	MAG 019 A	SILICONE GREASE	LAM 004 A
DYNAMOMETRIC SCREW DRIVER A301	MAG 020 A	GREASE FOR MECHANICAL ASSEMBLIES	XEA 381 AS
DYNAMOMETRIC SCREW DRIVER A300	MAG 013 A	SOFT TISSUE	
CLAMPS		LIQUID SOAP	
SCALPEL		DISTILLED WATER	
CUTTING PLIERS		MICROPIPETTE TIP	
PAIR OF SCISORS		FLAT PIECE OF STIFF PLASTIC	
VOLTMETER		LATEX WBC	LAD 001 AS
FLAT SCREWDRIVER		LATEX RBC	LAD 002 AS
BARFLEX		FELT PEN	
THERMOMETER		SYRINGE 5ML	
TORX KEYS			

# INSTALLATION



Hématologie

01/10/98

## • CONCERNS

- 1 - Unpacking
- 2 - Working conditions
- 3 - Visual checks
- 4 - Reagent connections
- 5 - Printer and instrument connections
- 6 - Priming and startup

## • REQUIRED TOOLS

Hexagonal keys

## • REQUIRED PRODUCTS

- MICROS 60 Reagents : Bottles or Pack.
- Waste container (for bottle model).

## • INTERVENTION TIME

- 30 minutes

## • FREQUENCY

## • SPECIFIC KIT OR CONSUMABLES

- PACK installation kit : XEA 314 A
- or
- Bottle installation kit : XEA 332 A

RAS 165 A Ind.B

MICROS 60 OT

- PROCEDURE

### 1 - Unpacking



*Diag.1*

The instrument is enveloped in a special, protective foam before being placed in a cardboard box. Cut the four angles of the box to unpack the system.

Remove the cardboard box containing the **ABX MICROS 60** installation kit from its location (see Diag.1).

### 2 - Working conditions

- Environment

**ABX MICROS 60-OT** should be operated in an indoor location only. Operation at an altitude over 2000 meters is not recommended. Instrument is designed to be safe for transient voltages according to INSTALLATION CATEGORY II and POLLUTION DEGREE 2.

- Location

**ABX MICROS 60-OT** should be placed on a clean and level table or work station. Please note that **ABX MICROS 60-OT**, printer and reagents weigh approximately 30 kilograms (66 lbs). Avoid exposure to sunlight. Proper ventilation requires that a space of at least 20 cm (8 inches) must be left behind the apparatus.

- Grounding

Proper grounding is required. Check that the wall ground (earth) plug is correctly connected to the laboratory grounding electricity installation. If there is no ground then use a ground stake. Current electricity norms must be applied.

- Humidity and temperature conditions

**ABX MICROS 60-OT** can function between 18 to 32°C (65 to 90°F), with relative humidity, meaning less than 80% with no condensation. If it is kept at a temperature less than 10°C (50°F), the instrument should be allowed to sit for an hour at the correct room temperature before use.

# INSTALLATION

## 3 - Visual checks




Diag.2

Using the key contained into the installation kit, turn the locker as shown on the Diag. 2 to open the pneumatic protection door.

Unscrew the 3 cover fixation screws (Diag. 3) and loosen the 2 tightening screws under the reagent flap (Diag 4).


Remove the cover : pull it backward and lift it up to the rear of the instrument .



Diag.3



Diag.4



Diag.5

Push the black plastic carriage locking clip (A) as far as possible to the left and place the sample needle carriage (B) as far forward as possible to the right-hand side, as shown in Diag. 5.


Check that the needle is not bent and make sure it is in its upper position.

MICROS 600T



## INSTALLATION

Unscrew slightly the 2 screws of the WBC/HGB chamber protection cover (diag 6). Remove the cover and check that both chambers (RBC/PLT, WBC/HGB) are fixed properly in their clips and the electrode blocks are attached firmly to the chambers (Diag. 7).




Diag.6



Diag.7


Re-install the HGB/WBC chamber cover.



Diag.8

Check that the connectors on the printed circuit board are securely in place (Diag. 8).

Re-install the instrument cover.



Diag.9

Remove the fuse holder from its location on the rear panel pressing on the holder lock (Diag.9) and check the fuse characteristics : they should be 1 Ampere, 250 Volts Slow-Blow.

# INSTALLATION

## 4 - Reagent connections

- Bottle connections



Diag.10

Lyse and cleaning reagents are placed inside the reagent compartment as shown in the Diag. 10. Install the reagent straws and the bottle stoppers. Connect the blue tube to the MINICLEAN bottle and the white tube to the MINILYSE bottle. Close the compartment cover.


## IMPORTANT

When the ABX MICROS 60-07 is set up with the 16 or 18 parameters mode, it is mandatory to use specific MINILYSE LMG and MINIDIL LMG reagents.

## CAUTION

The Diluent container will be located on the bench at the same level than the instrument.

- Waste connection



Diag.11

Install the male connectors included in the installation kit at the liquid input and output located at the bottom of the instrument rear panel (Diag.11).

Connect the diluent container (see CAUTION above) using the diluent straw and a 3x6 cristal tube (1 meter maximum) on the diluent input located at the bottom of the instrument rear panel (Diag.11). Connect the waste container using the cristal tube 3 x 6 on the waste output, and place the waste container below the instrument level (under the bench).

## WARNING !




Always follow the recommended procedures for waste disposal. Never connect the instrument wastes directly to the laboratory drain pipes. For each waste container, follow the neutralization procedure as described in the user manual.

MICROS 6007

## INSTALLATION

- Reagent pack connection

Remove the reagent output protections, as well as the waste input protection (Diag.12 & 13)



Diag.12



Diag.13

Install the pack directly into the compartment of the instrument as shown on the Diag. 14, 15 & 16. Push the pack down in order to plug correctly the pack on the male connectors.



Diag.14



Diag.15



Diag.16

The free male connector (see Diags. 14, 15 & 16) must be plugged on the pack upper valve in order to receive the waste liquids.

### CAUTION

In order to avoid leak problems it is recommended not to unplug several times the same reagent pack.


# INSTALLATION

## 5 - Printer and instrument connections

- Instrument connection.


If the instrument has to be connected to a laboratory computer, use the plug RS232.

Connect the power cable to the plug located on the rear left-hand side of the device (Diag. 17).



Diag.17

- Setting up the printer
  - Remove all the package protections
  - Install the paper feed-knob




Diag.18




Diag.19

Open the ink ribbon access door at the top of the printer and install the ribbon as shown in diag 20 & 21 :

- Slide the printer head to the middle of the printer.
- Insert the ribbon cartridge into the printer
- Guide the ribbon between the print head and ribbon guide.
- Slide the printer head from side to side to make sure it moves smoothly.




Diag.20



Diag.21


## INSTALLATION

- Install the paper supports (Diag 22) for single sheets paper use only.



*Diag.22*

- The printer is connected to **ABX MICROS 60-OT** with the cable delivered with the instrument. Lock the connector in place by tightening the 2 screws on each end of the connector to the **ABX MICROS 60-OT**. Attach the other end of the cable to the printer and lock the printer connector in place by the means of the 2 clips located on the connector itself (Diag. 23).




*Diag.23*

- For continuous paper, introduce it in the slot at the back of the printer and use the sprocket covers to load paper, or else feed the paper forward when using single sheets (see printer user manual).




*Diag.24*



*Diag.25*

# INSTALLATION

- Printer command keys



Diag.26

**LED PAUSE** : The orange LED PAUSE lights when the printer stops printing. During each power ON, this LED blinks for few seconds and 4 audible beeps occur. When the printer runs out of paper, the LED blinks and 3 audible beeps occur. This LED lights also when the paper is in its tear off position. When a problem occurs, this LED lights ON and 5 audible beeps occur.

**LEDS FONT 1 and FONT 2** : These 2 green LEDs indicate the selected font. Refer to the printer user's manual to select the font.

**Key FONT** : During normal operation, the FONT key allows the font selection. For each pressure on this key, the selection is modified. Refer to the printer user's manual to select the font. When this key is pressed during the printer power ON, the printer setup menu is entered.


**Key LF/FF** : During normal operation, a quick pressure on this key allows a ligne feed of the paper. Keep the pressure on this key to feed a whole page. This key can be used to load or eject the paper.

When this key is pressed during the printer power ON, the printing test starts. **Key PAUSE** : When this key is pressed during the printing, the printout stops. Press again on this key to restart the printout.

**PAPER PARK** : If Z folded paper is used, the paper can be driven to its parking position when pressing simultaneously on the keys LF/FF and FONT.

**MICRO ADJUST** : This function allows to adjust the loading paper position. See the user's manual for details.

- Printer Configuration  
Switch on the printer when pressing the <FONT>. The configuration should be the following :



Diag.29

# INSTALLATION

Each parameter can be modified by the corresponding parameter chart.  
 Each chart is accessible using the keys <PAUSE>, <FONT> and <LF/FF> according to the control LED combinations (Diag 30).

Lights				Table B	
1	2	PAUSE	Setting	Go to submenu	
BLINKS	OFF	OFF	Character spacing	Table C	
BLINKS	ON	OFF	Shape of zero	Table D	
OFF	BLINKS	OFF	Skip-over-perforation	Table E	
ON	BLINKS	OFF	Character table	Table F	
BLINKS	BLINKS	OFF	Auto line feed	Table G	
BLINKS	OFF	ON	Page length	Table H	
BLINKS	ON	ON	Auto tear off	Table I	
OFF	BLINKS	ON	Tractor	Table J	
ON	BLINKS	ON	Interface	Table K	
BLINKS	BLINKS	ON	Bit rate	Table L	
OFF	OFF	BLINKS	Parity	Table M	
BLINKS	OFF	BLINKS	Data length	Table N	
ON	OFF	BLINKS	ETX/ACK	Table O	
ON	ON	BLINKS	Software	Table P	
BLINKS	ON	BLINKS	Auto CR	Table Q	

Diag.30

## 6 - Instrument startup

- Reagent priming

When the **ABX MICROS 60-0T** is first installed, it contains no reagents. All the reagents have to be primed now. Turn ON instrument by pressing the ON/OFF switch located on the rear panel. When the instrument turns on, the display shows :

PLEASE WAIT FOR 3 MIN  
ESCAPE : ESC

This time is required at the startup for the instrument initialization and stabilization, specifically for the HGB diode to reach its operational temperature. Press ESC several times in order to abort the cycle : the LED of the front panel turns from red to green and the display shows the following :

STARTUP NOT INITIATED  
PRESS A KEY TO CONTINUE...

This message appears when the instrument is setup with the manual startup cycle to prevent any analysis cycle before running a startup cycle. Press any key, the main menu is displayed :

MAIN MENU	1 RESULTS
HH : MM	2 QC <span style="float: right;">▼</span>

# INSTALLATION

• *Bottles and containers set :*

Install the reagent bottle and carry out a PRIME cycle to clear the reagent line of air bubbles. This procedure should be done whenever a new bottle of reagent is installed.

From the MAIN MENU, move the cursor to the function (4) SERVICE and press ENTER. The service menu is displayed :

SERVICE	1 BACK FLUSH	▼
HH : MM	2 DRAIN CHAMBERS	

Move the cursor to function (3) PRIME REAGENTS and press ENTER. The PRIME menu is displayed :

PRIME	1 ALL REAGENTS	▼
HH : MM	2 DILUENT	

Select either the function (1) prime ALL REAGENTS or move the cursor next to the required reagent and press ENTER.

The priming cycle starts while the following menu is displayed :

ALL REAGENTS	WAIT FOR 2 MN 3 S
*****	

## IMPORTANT

**Before analyzing samples, visually inspect reagent lines and pumps for air bubbles. Repeat priming if air bubbles are still present. Call the *ABX representative service department* if priming does not eliminate air bubbles.**

**Never initiate two Lyse prime cycles back-to-back. This causes excessive foaming in the waste chamber. Run a blank cycle between each Lyse prime cycle.**

- Run a STARTUP cycle.



# INSTALLATION

## • Reagent pack

From the MAIN MENU, move the cursor to the function **4** SERVICE and press ENTER. The service menu is displayed :

SERVICE HH : MM	> 2 DRAIN CHAMBER 3 PRIME REAGENTS	▲ ▼
--------------------	---------------------------------------	--------

Move the cursor to **3** PRIME REAGENTS and press the ENTER key. Select the function **1** CHANGE PACK and follow the instructions given by the LCD in order to install the pack.

REAGENT PACK HH : MM	> 1 CHANGE PACK 2 CBC LEFT < 150>	▼
-------------------------	--------------------------------------	---

Once the new PACK is installed a priming cycle will be automatically carried out and the following menu is displayed.

PRIME *****	WAIT FOR 2 MIN 3 S
----------------	--------------------

## IMPORTANT

**Before analyzing samples, visually inspect reagent lines and pumps for air bubbles. Repeat priming if air bubbles are still present. Call the *ABX representative service department* if priming does not eliminate air bubbles.**

From the REAGENT PACK menu, the function **2** "CBC LEFT" displays the number of analysis cycles left to run with the same pack.

It is also possible to run a priming cycle at any time using the selection **3** "PRIME REAGENTS" of the SERVICE menu.

## IMPORTANT

**It is recommended not to remove the pack several times before the reagents are totally used in order to avoid leak problems.**

- Run a STARTUP cycle.

# INSTALLATION

Once the instrument ready for the analyses, remove the adhesive protection from the front panel (diag 31)



*Diag.31*

MICROS 60 OT

# INSTALLATION



Hématologie

01/10/98

## • CONCERNS

- 1 - Unpacking
- 2 - Working conditions
- 3 - Visual checks
- 4 - Reagent connections
- 5 - Printer and instrument connections
- 6 - Priming and startup

## • REQUIRED TOOLS

Hexagonal keys

## • REQUIRED PRODUCTS

Reagents : Bottles or Pack.  
Waste container (for bottle model).

## • INTERVENTION TIME

30 minutes

## • FREQUENCY

## • SPECIFIC KIT OR CONSUMABLES

PACK installation kit : XEA 317 A or  
Bottle installation kit : XEA 335 A or

MICROS 60 CT

RAS 166 A Ind.B

- PROCEDURE

**1 - Unpacking****Diag.1**

The instrument is enveloped in a special, protective foam before being placed in a cardboard box. Cut the four angles of the box to unpack the system.

Remove the cardboard box containing the instrument installation kit from its location (see Diag.1).

**2 - Working conditions**

- Environment

**ABX MICROS 60-CT** should be operated in an indoor location only. Operation at an altitude over 2000 meters is not recommended. Instrument is designed to be safe for transient voltages according to INSTALLATION CATEGORY II and POLLUTION DEGREE 2.

- Location

**ABX MICROS 60-CT** should be placed on a clean and level table or work station. Please note that **ABX MICROS 60-CT**, printer and reagents weigh approximately 30 kilograms (66 lbs). Avoid exposure to sunlight. Proper ventilation requires that a space of at least 20 cm (8 inches) must be left behind the apparatus.

- Grounding

Proper grounding is required. Check that the wall ground (earth) plug is correctly connected to the laboratory grounding electricity installation. If there is no ground then use a ground stake. Current electricity norms must be applied.

- Humidity and temperature conditions

**ABX MICROS 60-CT** can function between 18 to 32°C (65 to 90°F), with relative humidity, meaning less than 80% with no condensation. If it is kept at a temperature less than 10°C (50°F), the instrument should be allowed to sit for an hour at the correct room temperature before use.

# INSTALLATION



## 3 - Visual checks



Using the key contained into the installation kit, turn the locker as shown on the Diag. 2 to open the pneumatic protection door.


Unscrew the 3 cover fixation screws (Diag. 3) and loosen the 2 tightening screws under the reagent flap (Diag 4).  
Remove the cover : pull it backward and lift it up to the rear of the instrument.

Diag.2



Diag.4

Diag.3



Push the black plastic carriage locking clip (A) as far as possible to the left and place the sample needle carriage (B) as far forward as possible to the right-hand side, as shown in Diag. 5.  
Check that the needle is not bent and make sure it is in its upper position.

Diag.5

MICROS 60 CT

## INSTALLATION

Unscrew slightly the 2 screws of the WBC/HGB chamber protection cover (diag 6). Remove the cover and check that both chambers (RBC/PLT, WBC/HGB) are fixed properly in their clips and the electrode blocks are attached firmly to the chambers (Diag. 7).



Diag.6



Diag.7

Re-install the HGB/WBC chamber cover.



Diag.8

Check that the connectors on the printed circuit board are securely in place (Diag. 8).

Re-install the instrument cover.



Diag.9

Remove the fuse holder from its location on the rear panel pressing on the holder lock (Diag.9) and check the fuse characteristics : they should be 1 Ampere, 250 Volts Slow-Blow.

# INSTALLATION

## 4 - Reagent connections

- Bottle connections



Diag.10

Lyse and cleaning reagents are placed inside the reagent compartment as shown in the Diag. 10. Install the reagent straws and the bottle stoppers. Connect the blue tube to the MINICLEAN bottle and the white tube to the MINILYSE bottle. Close the compartment cover.


## IMPORTANT

When the ABX MICROS 60-CT is set up with the 16 or 18 parameters mode, it is mandatory to use specific MINILYSE LMG and MINIDIL LMG reagents.

## CAUTION

The Diluent container will be located on the bench at the same level than the instrument.

- Waste connection



Diag.11

Install the male connectors included in the installation kit at the liquid input and output located at the bottom of the instrument rear panel (Diag.11).

Connect the diluent container (see CAUTION above) using the diluent straw and a 3x6 cristal tube (1 meter maximum) on the diluent input located at the bottom of the instrument rear panel (Diag.11). Connect the waste container using the cristal tube 3 x 6 on the waste output, and place the waste container below the instrument level (under the bench).

## WARNING !




Always follow the recommended procedures for waste disposal. Never connect the instrument wastes directly to the laboratory drain pipes. For each waste container, follow the neutralization procedure as described in the user manual.

## INSTALLATION

- Reagent pack connection

Remove the reagent output protections, as well as the waste input protection (Diag.12 & 13)



Diag.12



Diag.13

Install the pack directly into the compartment of the instrument as shown on the Diag. 14, 15 & 16. Push the pack down in order to plug correctly the pack on the male connectors.



Diag.14



Diag.15



Diag.16

The free male connector (see Diags. 14, 15 & 16) must be plugged on the pack upper valve in order to receive the waste liquids.

## CAUTION

In order to avoid leak problems it is recommended not to unplug several times the same reagent pack.




# INSTALLATION

## 5 - Printer and instrument connections

- Instrument connection.


If the instrument has to be connected to a laboratory computer, use the plug RS232.

Connect the power cable to the plug located on the rear left-hand side of the device (Diag. 17).



Diag.17

- Setting up the printer
  - Remove all the package protections
  - Install the paper feed-knob




Diag.18




Diag.19

Open the ink ribbon access door at the top of the printer and install the ribbon as shown in diag 20 & 21 :

- Slide the printer head to the middle of the printer.
- Insert the ribbon cartridge into the printer
- Guide the ribbon between the print head and ribbon guide.
- Slide the printer head from side to side to make sure it moves smoothly.




Diag.20



Diag.21


## INSTALLATION

- Install the paper supports (Diag 22) for single sheets paper use only.



*Diag.22*

- The printer is connected to **ABX MICROS 60-CT** with the cable delivered with the instrument. Lock the connector in place by tightening the 2 screws on each end of the connector to the **ABX MICROS 60-CT**. Attach the other end of the cable to the printer and lock the printer connector in place by the means of the 2 clips located on the connector itself (Diag. 23).




*Diag.23*

- For continuous paper, introduce it in the slot at the back of the printer and use the sprocket covers to load paper, or else feed the paper forward when using single sheets (see printer user manual).




*Diag.24*



*Diag.25*

# INSTALLATION

- Printer command keys



Diag.26

**LED PAUSE** : The orange LED PAUSE lights when the printer stops printing. During each power ON, this LED blinks for few seconds and 4 audible beeps occur. When the printer runs out of paper, the LED blinks and 3 audible beeps occur. This LED lights also when the paper is in its tear off position. When a problem occurs, this LED lights ON and 5 audible beeps occur.

**LEDS FONT 1 and FONT 2** : These 2 green LEDs indicate the selected font. Refer to the printer user's manual to select the font.

**Key FONT** : During normal operation, the FONT key allows the font selection. For each pressure on this key, the selection is modified. Refer to the printer user's manual to select the font. When this key is pressed during the printer power ON, the printer setup menu is entered.


**Key LF/FF** : During normal operation, a quick pressure on this key allows a ligne feed of the paper. Keep the pressure on this key to feed a whole page. This key can be used to load or eject the paper.

When this key is pressed during the printer power ON, the printing test starts. **Key PAUSE** : When this key is pressed during the printing, the printout stops. Press again on this key to restart the printout.

**PAPER PARK** : If Z folded paper is used, the paper can be driven to its parking position when pressing simultaneously on the keys LF/FF and FONT.

**MICRO ADJUST** : This function allows to adjust the loading paper position. See the user's manual for details.

- Printer Configuration  
Switch on the printer when pressing the <FONT>. The configuration should be the following :



Diag.29

# INSTALLATION

Each parameter can be modified by the corresponding parameter chart.  
 Each chart is accessible using the keys <PAUSE>, <FONT> and <LF/FF> according to the control LED combinations (Diag 30).

Lights				Table B	
1	2	PAUSE	Setting	Go to submenu	
BLINKS	OFF	OFF	Character spacing	Table C	
BLINKS	ON	OFF	Shape of zero	Table D	
OFF	BLINKS	OFF	Skip-over-perforation	Table E	
ON	BLINKS	OFF	Character table	Table F	
BLINKS	BLINKS	OFF	Auto line feed	Table G	
BLINKS	OFF	ON	Page length	Table H	
BLINKS	ON	ON	Auto tear off	Table I	
OFF	BLINKS	ON	Tractor	Table J	
ON	BLINKS	ON	Interface	Table K	
BLINKS	BLINKS	ON	Bit rate	Table L	
OFF	OFF	BLINKS	Parity	Table M	
BLINKS	OFF	BLINKS	Data length	Table N	
ON	OFF	BLINKS	ETX/ACK	Table O	
ON	ON	BLINKS	Software	Table P	
BLINKS	ON	BLINKS	Auto CR	Table Q	

Diag.30

## 6 - Instrument startup

- Reagent priming

When the **ABX MICROS 60-CT** is first installed, it contains no reagents. All the reagents have to be primed now. Turn ON instrument by pressing the ON/OFF switch located on the rear panel. When the instrument turns on, the display shows :

PLEASE WAIT FOR 3 MIN  
 ESCAPE : ESC

This time is required at the startup for the instrument initialization and stabilization, specifically for the HGB diode to reach its operational temperature. Press the ESC key in order to abort the cycle : the LED of the front panel turns from red to green and the display shows :

STARTUP NOT INITIATED  
 PRESS A KEY TO CONTINUE...

This message appears when the instrument is setup with the manual startup cycle to prevent any analysis cycle before running a startup cycle. Press any key, the main menu is displayed :

MAIN MENU	1 RESULTS
HH : MM	2 QC <span style="float: right;">▼</span>

# INSTALLATION

## Bottles and containers set :

Install the reagent bottle and carry out a PRIME cycle to clear the reagent line of air bubbles. This procedure should be done whenever a new bottle of reagent is installed.

From the MAIN MENU, move the cursor to the function (4) SERVICE and press ENTER. The service menu is displayed :

SERVICE	1 BACK FLUSH	▼
HH : MM	2 DRAIN CHAMBERS	

Move the cursor to function (3) PRIME REAGENTS and press ENTER. The PRIME menu is displayed :

PRIME	1 ALL REAGENTS	▼
HH : MM	2 DILUENT	

Select either the function (1) prime ALL REAGENTS or move the cursor next to the required reagent and press ENTER.

The priming cycle starts while the following menu is displayed :

ALL REAGENTS	WAIT FOR 2 MN 3 S
*****	

## IMPORTANT

Before analyzing samples, visually inspect reagent lines and pumps for air bubbles. Repeat priming if air bubbles are still present. Call the *ABX representative service department* if priming does not eliminate air bubbles.

Never initiate two Lyse prime cycles back-to-back. This causes excessive foaming in the waste chamber. Run a blank cycle between each Lyse prime cycle.

- Run a STARTUP cycle.

# INSTALLATION

## Reagent pack

From the MAIN MENU, move the cursor to the function **4** SERVICE and press ENTER. The service menu is displayed :

SERVICE HH : MM	> 2 DRAIN CHAMBER 3 PRIME REAGENTS	▲ ▼
--------------------	---------------------------------------	--------

Move the cursor to **3** PRIME REAGENTS and press the ENTER key. Select the function **1** CHANGE PACK and follow the instructions given by the LCD in order to install the pack.

REAGENT PACK HH : MM	> 1 CHANGE PACK 2 CBC LEFT < 150>	▼
-------------------------	--------------------------------------	---

Once the new PACK is installed a priming cycle will be automatically carried out and the following menu is displayed.

PRIME *****	WAIT FOR 2 MIN 3 S
----------------	--------------------

## IMPORTANT

**Before analyzing samples, visually inspect reagent lines and pumps for air bubbles. Repeat priming if air bubbles are still present. Call the *ABX representative service department* if priming does not eliminate air bubbles.**

From the REAGENT PACK menu, the function **2** "CBC LEFT" displays the number of analysis cycles left to run with the same pack.

It is also possible to run a priming cycle at any time using the selection **3** "PRIME REAGENTS" of the SERVICE menu.

## IMPORTANT

**It is recommended not to remove the pack several times before the reagents are totally used in order to avoid leak problems.**

- Run a STARTUP cycle.

# INSTALLATION

Once the instrument ready for the analyses, remove the adhesive protection from the front panel (diag 31)



*Diag.31*

MICROS 60 CT

## SAMPLING NEEDLE MAINTENANCE



Hématologie

18/06/98

### • CONCERNS

- Needle O ring replacement
- Sampling needle replacement
- Piercing needle replacement

### • REQUIRED TOOLS

- Hexagonal keys
- Dynamometric screw driver A302 : MAG 019 A

### • REQUIRED PRODUCTS

- Silicone grease : LAM 004 A

### • INTERVENTION TIME

- 15min

### • FREQUENCY

- See frequency chart table for cleaning.
- O ring replacement : 1/year
- Needle replacement : On request only


### • SPECIFIC KIT OR CONSUMABLES

- O ring kit : XEA 328 AS
- Sampling needle : GBC 052 A
- Piercing needle : GBC 189 A



# SAMPLING NEEDLE MAINTENANCE

## • PROCEDURE



Diag.1


### A - Needle or O ring replacement

- Disconnect the tube (1) from the top of the needle (2) (diagram 1).

Remove the clip (3) (diagram 2) and manually lift up the sampling needle (Diagram 1).

- Unscrew the 2 screws (3) (Diagram 1) in order to freed the needle rining block (4) from the carriage frame

(5).



Diag.2

- Remove the rining block/needle assy from the carriage taking care not to bend the needle.

- Lift up the O ring holder (1) and replace the O rings (2) by new ones previously greased. Wipe all excess of grease away.

- If necessary clean the inner surface of the rining block with a little piece of paper.


Proceed the same way to replace the needle if necessary

- Reassemble in the reverse order. Use a dynamometric screw driver to tighten the screws (3) (Diag.1) to **700 mN.m (99.4 Ozf.in)**.

# SAMPLING NEEDLE MAINTENANCE

## B - Piercing needle replacement

- Disconnect the tube (1) from the piercing needle (4) (Diag.3).



Diag.3

- Lift the needle (2) in the upper position.

- Loosen the 2 screws (3) just enough to enable the rinsing block to be lifted up of about 5 mm.


- Pull the piercing needle (4) (forward) and replace it by a new one if necessary.

- Replace the piercing needle O ring (5) (FAA 036 A) by a new one.

- Reassemble in the reverse order. Apply the same torque : **700mN.m (99.4 Ozf.in)** to tighten the two screws (3).

## NOTE

It is recommended to check the correct motion of the needle. Proceed as follows :  
Enter the «SERVICE» menu and then the sub menu «MECHANIC» and perform a «NEEDLE U/D» cycle.



Diag.4

Blockage problems may occur on some instruments during the needle or carriage motions giving some motor error messages.

Before replacing the concerned motor, it is necessary to check the correct rotation of the free pulley located at the end of the notched belt.

Remove the axle screw of the pulley and clean its 2 parts and the washer. Reinstall the pulley assy, the rounded edge facing the pulley. Tighten the screw with a torque of **400mN.m (56.8 Ozf.in)**.

Check that the pulley turns freely after the tightening. Add a drop of oil (LAM 007 A) if necessary.

- 1 - FAG 011 A : Autolub. axle
- 2 - GBC 146 A : Free pulley
- 3 - GBC 147 A : Pulley holder

## SAMPLING NEEDLE MAINTENANCE



Hématologie

10/04/98

### • CONCERNS

- Needle replacement
- O ring replacement

### • REQUIRED TOOLS

- Hexagonal keys : 2,5
- Dynamometric screw driver A300 : MAG 013 A

### • REQUIRED PRODUCTS

- Silicone grease : LAM 004 A

### • INTERVENTION TIME

- 15 min

### • FREQUENCY

- Needle replacement : On request only
- O ring replacements : 1/year


### • SPECIFIC KIT OR CONSUMABLES

- O ring kit : XEA 328 AS
- Spare parts kit : XEA 458 AS

**SAMPLING NEEDLE MAINTENANCE**

## • PROCEDURE

- Disconnect the tube (7) from the top of the needle. Manually lift up the sampling needle (Diagram 1).
- Unscrew the 2 screws (5) in order to freed the needle rinsing block (2) from the carriage frame. Remove the clip (6).



Diag.1

- Remove the rinsing block/needle assy from the carriage taking care not to bend the needle.
- If necessary clean the inner surface of the rinsing block by means of a little piece of soft paper. Spread a little amount of grease in between the rinsing block (2) and its support.
- Lift up the O ring holder (4) and replace the O rings (3) by new ones previously greased. Wipe all excess of grease away.

**NOTE**

Proceed the same way to replace the needle if necessary

# SAMPLING NEEDLE MAINTENANCE

- Reassemble in the reverse order. Use a dynamometric screw driver to tighten the screws <sup>⑤</sup> (Diag.1) to 100 mN.m.

## NOTE

It is recommended to check the correct motion of the needle. Proceed as follows :  
Enter the «SERVICE» menu and then the sub menu «MECHANIC» and perform a «NEEDLE UP/DOWN» cycle.


Blockage problems may occur on some instruments during the needle or carriage motions giving some motor error messages.

Before replacing the concerned motor, it is necessary to check the correct rotation of the free pulley located at the end of the notched belt.

Remove the axle screw of the pulley and clean its 2 parts and the washer. Reinstall the pulley assy, the rounded edge facing the pulley. Tighten the screw with a torque of **400mN.m (99.4 Ozf.in)**.

Check that the pulley turns freely after the tightening. Add a drop of oil (LAM 007 A) if necessary.

- ① - FAG 011 A : Autolub. axle
- ② - GBC 146 A : Free pulley
- ③ - GBC 147 A : Pulley holder



Diag.2

# CHAMBER MAINTENANCE



Hématologie  
10/04/98

## • CONCERNS

- RBC & WBC/HGB cleaning
- Aperture O ring replacement
- Coaxial O ring replacement

## • REQUIRED TOOLS

- Hexagonal keys : 2,5
- Soft paper
- Dynamometric screw driver : MAG 013 A
- Cutting pliers

## • REQUIRED PRODUCTS

- Liquid soap
- Distilled water
- A scalpel
- A Micropipette tip

## • INTERVENTION TIME

- 30 min

## • FREQUENCY

- RBC & WBC cleaning : 2 (type 1 & 2) or 3/year (type 3).
- Aperture O ring replacement : 1/year
- Electrode O ring replacement : 1/year

## • SPECIFIC KIT OR CONSUMABLES


- O ring kit : XEA 328 AS

RAS 169 A Ind.A

M I C R O S 6 0


### • PROCEDURE

#### 1 - RBC chamber cleaning :



Diag.1

- Run a drain chamber cycle (SERVICE menu, DRAIN CHAMBERS Sub menu).
- Record the tube positions before dismantling the chambers.
- Disconnect the chamber tubes.
- Unclip the RBC chamber.
- Dismantle the electrode
- ⑥ loosening the 2 fixation screws ⑦ and the terminal holding screw ③ (Diagram 1).



Diag.2

- *Electrode O ring replacement :*
- Use a previously cut micropipette tip to replace the electrode O ring as shown on Diag.2.
- *Aperture O ring replacement :*
- Install the chamber over a piece of white paper or cloth.
- Carefully remove the counting head ④ and plunge the aperture ⑤ in distilled water.
- Replace the O ring ② by a new one.

- Clean the chamber and the counting head with liquid soap, do not introduce any sharp instruments inside so as to avoid damaging the inside of the chamber and the aperture.

## CAUTION

Do not manipulate the aperture using hard instruments. Clean the aperture with a piece of soft paper or preferably, in between 2 fingers.


- Rinse thoroughly with distilled water
- Dry the exterior of the chamber with a soft paper.

## CAUTION

- Do not apply too much pressure on the electrode fixation screws, as it can break the aperture (tightening torque = 100mN.m / 14.2 Ozf.in).
- It is recommended to reconnect the tubes on the counting head before reassembling the "electrode/chamber/counting head" assy in order to avoid applying constraint on the chamber.

- Position the chamber in its fixation clips.
- Reconnect the tubes

### 2 - WBC/HGB chamber cleaning :



Diag.3

- Run a drain chamber cycle (SERVICE menu, DRAIN CHAMBERS Sub menu).
- Loosen the cover screws of the WBC/HGB chamber and remove the cover.
- Record the tube positions before dismantling the chambers.
- Disconnect the chamber tubes.
- Unclip the RBC chamber.
- Dismantle the electrode (6) loosening the 2 fixation screws (7) and the terminal holding screw (3).
- Proceed as described in 1 - **RBC chamber** to clean the chamber and to replace the electrode and aperture O rings.

## NOTE

The spectrophotometer can not be dismantled from the chamber. If this one has been damaged it is necessary to replace the whole chamber assy. When cleaning the spectrophotometer, make sure to thoroughly rinse it in order to obtain a correct HGB blank measure.



# LIQUID VALVE MAINTENANCE



Hématologie

14/04/98

## • CONCERNS

Liquid valve assy replacement  
Valve body replacement

## • REQUIRED TOOLS

Hexagonal keys  
Pair of pliers

## • REQUIRED PRODUCTS

Soft paper

## • INTERVENTION TIME

15 min

## • FREQUENCY

On request only

## • SPECIFIC KIT OR CONSUMABLES

6 Valve assembly (**MICROS 60 CT**) : XDA 579 CS  
6 Valve assembly (**MICROS 60 OT**) : XDA 578 CS  
5 Valve assembly : XDA 580 CS  
2 ways NC liquid valve without solenoïd : XDA 481 B  
3 ways liquid valve without solenoïd : XDA 483 B  
Solenoid 24V 4W : EAZ 004 A


RAS 170 A Ind.A

M I C R O S 6 0

## LIQUID VALVE MAINTENANCE

### • PROCEDURE

- Several reasons can enable the correct operations of the valves :
  - Leaks on the valve bodies
  - Defective Solenoids
  - Liquid discharge on valves
  - Corrosion traces on the axis, ect...
- Run a «DRAIN CHAMBERS» cycle.
- Switch off the instrument.
- Note the tube positions on the valve assembly that requires to be dismantled.
- Unscrew the valve assembly fixation screws (Diag.1).



*Diag.1*



*Diag.2*

- Disconnect the supplying flat cable from the valve assembly (diag.2).

### **NOTE**

When leaks occur on valves it is recommended to replace the entire valve assembly by a new one.


- If only one valve has been damaged it is possible to dismantle the valve body on its own as shown on the diagram 3 : use a pair of pliers to disconnect the valve holder clip and remove the body.

### **CAUTION**

When replacing one valve only, it is important to check if the «neighboured» valves have not been damaged too.

# LIQUID VALVE MAINTENANCE

When replacing one valve only, it is important to check if the «neighboured» valves have not been damaged too.



Diag.3

## CAUTION

The solenoid can not be dismantled unless unsoldering it.  
If this one is suspected solder a new one or replace the entire valve assembly.

- Reconnect the tubes on the nipples (with the sleeves).
- Re-install in the reverse order. Switch on the instrument.
- Control the watertightness of the valves and check for the correct operations :  
Go to «SERVICE» menu, then to «MECHANIC» sub menu (5) and require a «VALVES» (6) test (see Mechanic functions : RAS 173 A).
- Check the calibration too.

### • REMARKS

On **MICROS 60 - CT**, when piercing several times the same tube cap, some pieces of cork may be dragged along towards the WBC/HGB chamber and then the liquid valve 12.  
This may damaged the operation of the valve.

M  
I  
C  
R  
O  
S  
  
6  
0

# POWER SUPPLY CHECK OR REPLACEMENT



Hématologie

10/04/98

## • CONCERNS

- Voltage supply check
- Power supply module replacement
- Fan operation check

## • REQUIRED TOOLS

- Hexagonal keys
- Flat screw driver
- Voltmeter

## • REQUIRED PRODUCTS

- None

## • INTERVENTION TIME

- 15 min

## • FREQUENCY

- 1/year

## • SPECIFIC KIT OR CONSUMABLES

- None

RAS 171 A Ind.A

M I C R O S 6 0

# POWER SUPPLY CHECK OR REPLACEMENT

## • PROCEDURE

### 1 - Supply voltage check

## CAUTION

The supply voltage check has to be done with the power supply module connected only.


- The supply voltage check is done on the mother board (XAA 355 A) according to the following chart table :

\* Ground on TP31, TP 30 or TP 29

TEST POINTS	VOLTAGES
TP 20	-12V $\pm$ 0,5V
TP 22	24V + 1.5V - 0V
TP 23	5V + 0,3V - 0V
TP21	12V $\pm$ 0,5V

## CAUTION


The two screws on the top of the coaxial cover holding the board on the MICROS frame are the grounding connection of the mother board (diag 1). It is mandatory to check the correct tightening of these screws to obtain correct voltage values.



Diag.1

# POWER SUPPLY CHECK OR REPLACEMENT

The test points are located as shown on the diagram below :



Diag.2

## NOTE


If the voltages values are not correct or among the ranges no adjustment can be carried out either on the board or on the power supply module. Replace the power supply module as described below :

MICROS 60


# POWER SUPPLY CHECK OR REPLACEMENT

## 2 - Power supply module replacement

- Switch off the instrument.
- Disconnect the main supply voltage cable and the printer cable.
- Disconnect the power supply cable from the mother board, connector J26 (see Diag.2).
- Unscrew the 4 screws ① (Diag.3) and start to move out the power supply module. Disconnect the grounding wire (Yellow/green wire see diag 4) from the rear panel of the module. Route the cable (from J26) down to make the removing of the module easier.



Diag.3



Diag.4

## WARNING !



The power supply module internal fuse is not to be replaced even when this one has blown down.

- Replace the power supply module by a new one and reinstall in the reverse order.
- Switch on the instrument and check the voltages on the mother board as described in the previous paragraph.

### • REMARKS

The four leds in front of the test points are lit to indicate a voltage presence but whatever its value!!!

Check the operation of the fan as following : When the fan has stopped, move the cursor to the 4 SERVICE menu and press ENTER : the fan should start.

# TECHNICIAN FUNCTION



Hématologie

16/06/98

## • CONCERNS

- |  |                                   |
|--|-----------------------------------|
| 1 - Version display.                         | 3 - Temperature sensor adjustment |
| 2 - Adjustments                              | 4 - Run mode                      |
| HGB photometer calibration, Aperture voltage | 5 - Reagent pack                  |
| Vacuum check, WBC gain, RBC & PLT gain       | 6 - Serial number                 |
| Sensor, Needle height, Needle motion         | 7 - Cycle number                  |
| Bubbling                                     | 8 - Burn-in                       |

## • REQUIRED TOOLS

- Flat screw driver
- thermometer
- Barflex
- Hexagonal keys
- Voltmeter

## • REQUIRED PRODUCTS

- WBC latex : LAD 001 AS
- RBC and PLT latex : LAD 002 AS
- Soft paper
- Flat piece of stiff plastic

## • INTERVENTION TIME

- 60 min

## • FREQUENCY

- See maintenance chart table.

## • SPECIFIC KIT OR CONSUMABLES


- Needle position tool : GBC 218 A

MICROS 600T

RAS 172 A Ind.A







• PROCEDURE



**MICROS 60 OT**

# TECHNICIAN FUNCTION

Once entering the «SERVICE» sub menu, move to «TECHNICIAN FUNCTIONS» and press . A specific password (421) is required to enter the sub menus.

Move the cursor by means of  and  and choose the required menus pressing the  key.

I - The version number is displayed.

## II - Adjustments :

### 1 - HGB photometer calibration

- Dismantle the WBC/HGB chamber cover.
- Check the general cleanliness of the WBC chamber/spectrophotometer assy.
- Re-install the chamber cover.

**NOTE** If the WBC chamber has been dismantled previously make sure no liquid has flown in between the spectrophotometer and the chamber. Clean the inner surfaces of the spectrophotometer as well as the chamber. Reassemble the assy and tighten the two screws to the following torque : 400mN.m (see RAS 169 A : Chamber maintenance)

- Run the CAL PHOTOMETER function (function 1 of the «ADJUSTMENTS» menu) : diluent is delivered to the WBC/HGB chamber twice.

An HGB channel is displayed on the LCD screen :

VALUE XXX
--------------


## IMPORTANT

- 1 - The HGB photometer calibration must be done 20min at least after the instrument has been switched on.
- 2 - This adjustment must be done with the WBC chamber cover installed!!!

- By means of **R97** (See diagram 1) adjust the HGB channel according to the room temperature using the chart table given on next page.

**NOTE** After 40 seconds approximately, the function is automatically exited.

- Run the CAL PHOTOMETER function again to verify the adjustment.



Diag.1

ROOM TPT (°C)	CHANNEL		
	Mini.	Nominal	Maxi.
15	240	245	250
16	240	245	250
17	239	244	249
18	238	243	248
19	237	242	247
20	236	241	246
21	235	240	245
22	234	239	244
23	234	239	244
24	233	238	243
25	232	237	242
26	231	236	241
27	230	235	240
28	229	234	239
29	228	233	238
30	228	233	238
31	227	232	237
32	226	231	236
33	225	230	235
34	224	229	234
35	223	228	233


# TECHNICIAN FUNCTION

## 2 - Aperture voltage


- Once entering the «APERTURE VOLTAGE» menu, connect a voltmeter between the ground (TP30) and TP19.
- Check that the value is **60V +2.8V -1.5V** (Diagram 1). The aperture voltage is not adjustable.
- Press any key to escape.

## 3 - Vacuum check

Enter the 2 - ADJUSTMENTS / 3 - CHECK PRESSURE menu.



Diag.2

- Disconnect the tube from the vacuum/waste syringe coming from the valve  (see diagram 2).

- Follow the instructions given on the LCD screen :

«PLEASE PLUG BARFLEX ON AIR SYRINGE» (On the free nipple).

«CHECK PRESSURE : -200mB ± 10mB»  
(The piston has raised in order to create a vacuum in the syringe body).

- Check the stability of the vacuum during 30 secondes : The vacuum drop down must be  $\leq 2$  mbar.

- If the results are not correct check the O ring and the tubing watertightness.

«PLEASE RECONNECT PREVIOUS TUBE» (disconnect the Barflex and replugin the tube instead).

## 4 - WBC adjust

- Put the WBC latex to mix on a Vortex during 1min or shake thoroughly

## IMPORTANT

As the WBC gain is a factory adjustment it is mandatory not to readjust it without taking the following precautions :

Carry out previously, an autoconcentrated cleaning to make sure of the cleanliness of the WBC counting circuit.

If necessary clean the WBC chamber aperture as described in RAS 169 A.

Make sure the Latex has been thoroughly mixed before.

- Run a blank cycle to check the cleanliness of the instrument.
- Enter the 4 - WBC ADJUST sub menu and then 1 - WBC ADJUST

## TECHNICIAN FUNCTION



- Present the vial of Latex to the open probe as shown on diagram 3 and press the sampling bar located behind the sampling needle : an analysis cycle begins.

During the cycle measuring phasis (around 1 minute) the **Lymphocyte** and the **Granulocyte** volumes are displayed on the screen every 3 seconds as shown below :

Diag.3

LYM <57 +/- 1> 57	GRA <180 +/- 2> 180
----------------------	------------------------

- Wait for several results to be displayed and check the stability of both values.


## IMPORTANT

It is mandatory not to operate the gain adjustment as long as the lympho and granulo values are not stable.

After one minute the last volume values displayed on the screen are printed out as well as the WBC, PLT and RBC histograms and the CBC's results. Check that these printed values correspond to the following target values :

- LYM =  $57 \pm 1$
- GRA =  $180 \pm 2$

- If not rerun a «WBC ADJUST» cycle after having previously mixed the latex vial again. By means of R74 (see diagram 4) adjust the volumes to the target values during the measuring phasis.



Diag.4

# TECHNICIAN FUNCTION

## NOTE

Both sub menus 2 - # LYM < > and 3 - # GRA < > allow the technician to change the Latex target values if the latex run on the instrument different from the Latex recommended above.

### 5 - RBC PLT adjust

- Put the RBC and PLT latex to mix on a Vortex during 1min or shake thoroughly

## IMPORTANT

As the RBC/PLT gain is a factory adjustment it is mandatory not to readjust it without taking the following precautions :

Carry out previously an autoconcentrated cleaning to make sure of the cleanliness of the RBC/PLT counting circuit.

If necessary clean the RBC/PLT chamber aperture as described in the procedure RAS 169 A.

Make sure the Latex has been thoroughly mixed before.

- Run a blank cycle to check the cleanliness of the instrument.
- Enter the «5 - RBC PLT ADJUST» sub menu.
- Present the vial of Latex to the open probe as shown on diagram 3 and press the sampling bar located behind the sampling needle : the needle directly delivers the latex sample in the RBC chamber dilution and a measuring phasis begins.
- During the cycle measuring phasis (around 1 minute) the **Platelet** and the **Red Blood cell** volumes are displayed on the screen every 3 seconds as shown below :

RBC <74 +/- 1> 74	PLT <59 +/-1> 59
----------------------	---------------------

- Wait for several results to be displayed and check the stability of both values.

## IMPORTANT

It is mandatory not to operate the gain adjustment as long as the platelet and RBC values are not stable.

After one minute the last volume values displayed on the screen are printed out as well as the PLT and RBC histograms and the CBC's results.

- Check that these printed values correspond to the following target values :

- RBC = 74 ± 1
- PLT = 59 ± 1

## NOTE

From the latex lot # 980311 included, balls having a different size, a drift of MPV peak has been noticed, i.e. a modification of the PLT gain target value : it becomes 64 instead of 59.

The program default value will be modified in the next MICROS version.

Both sub menus 2 - MCV < > and 3 - MPV < > allow the technician to change the Latex target values. If the lot # ≥ 980311, modify the target values and proceed the same way to adjust the PLT gain.

MICROS 600T

- If not, rerun a «RBC PLT ADJUST» cycle after having previously mixed the latex vial again.
- Adjust the volumes to the target values during the measuring phasis by means of the following potentiometers (see Diagram 4) :

- RBC gain with R81
- PLT gain with R87

### 6 - Needle sensor adjust

If the needle detector has been replaced by a new one or dismantled for any reason, it is mandatory to re-position it at the right heighth. Proceed as following :



Diag.5

- Install the piece of plastic (diagram 5) underneath the needle rinsing block.

- Once entering the menu «6 - SENSOR ADJUST», enter the sub menu «1 - NEEDLE SENSOR».

- Push the sampling needle downward until it stops against the piece of plastic and press any key in order to raise the needle back in its upper position.

The current number of steps, the mini and maxi values are displayed as well as the way to move the sensor (shown by an arrow) if the current value is out of ranges (see below).

CURRENT : 62

MIN : 65

MAX : 75



Diag.6

- For a current number of steps out of ranges, unloosen the 2 cell fixation screws (diagram 6) and gently move the sensor

- upward if the current value is too low
- downward if the current value is too high.

- Tighten the screws and rerun a «NEEDLE SENSOR» cycle. Check that the current value is correct.

## NOTE

10 steps correspond to around 1 mm. The target number of steps is  $70 \pm 5$

# TECHNICIAN FUNCTION

## 7 - Position adjustment

### • Needle height adjustment

When replacing a needle, it is mandatory to adjust the height of the needle in the chambers. Proceed as follows :

- Remove the WBC/HGB chamber cover.
- Enter the menu 7 - POSITION ADJUST / 1 - NEEDLE / 2 - ADJUSTMENT.
- Manually pull down the needle until it comes into contact with the edge of the RBC chamber (Diag. 7).
- Press a key to continue : the needle moves up to the initial position ; the value is stored.



Diag.7



Diag.8

### • Needle motion adjustment

The needle position in the WBC chamber can be adjusted as follows :

- Enter the menu 7 - POSITION ADJUST / 2 - CARRIAGE / 2 - ADJUSTMENT.
- Position the tool (P/N : GBC 218 A) over the RBC and WBC chambers (Diag 8).
- Manually lower the needle into the WBC chamber.
- Press a key to continue : the needle moves up, and the carriage comes back to the initial position ; the value is stored.
- Carry out a 3 - CHECK cycle : the needle comes down to the WBC chamber.
- Check that the needle is centered in the hole.

If not, enter the 2 - CARRIAGE / 1 - CARRIAGE < > menu ; the display shows the current number of steps carriage motion.

If the needle goes too far on the right, add 1 step to the current value for 0.1mm.

If the needle is too much on the left, decrease the current value of 1 step for 0.1mm.

CARRIAGE ? :  
CURRENT : 893

EXIT : ESC  
SAVE : ENTER

Carry out a 2 - CARRIAGE / 3 - CHECK cycle again to control the needle position.

MICROS 600T



**8 - Bubbling adjustment**

An overflow protection tank is installed on the drain circuit of each chamber. This one prevents from polluted liquid overflow during bubbling phasis (Diag .9).



*Diag.9*

Two bubbling phasis are adjustable :

- "BUBBLING 1" is the first dilution (WBC/HGB chamber) bubbling value.
- "BUBBLING 2" is the second dilution (WBC/HGB chamber + lyse) value and RBC chamber bubbling value.

Both values correspond to a number of steps carried out by the waste/vacuum syringe.

Default values are        BUBBLING 1 : 175  
                                       BUBBLING 2 : 120

## CAUTION

These values are factory adjusted (and may be different from the default values shown above) and should be modified only when hematologic erroneous results are given by the instrument : If values are too important, liquid overflows can occur or if bubbling is too low homogeneity of the dilution can be decreased.

Ranges :    150 < BUBBLING 1 < 200  
                   80 < BUBBLING 2 < 140

To modify the bubbling values, enter the menu :

- 2 - ADJUSTMENTS / 8 - BUBBLING ADJUST. / BUBBLING 1 < >
- 2 - ADJUSTMENTS / 8 - BUBBLING ADJUST. / BUBBLING 2 < >

and type in new step value.


Carry out a 3 - CHECK to control the adjustment.

# TECHNICIAN FUNCTION

## III - Temperature

When entering the «temperature adjustment» menu the following sub menus are displayed :

### 1 - Uncal. temperature

When pressing the  key the sensor temperature value **uncalibrated** is displayed. This value should be close to the diluent temperature.


### 2 - Calibrate :

The temperature must be calibrated according to the diluent temperature :



- Plunge a thermometer directly into the diluent container and leave it for a while until stabilization.

## NOTE

For a pack equipped instrument, the thermometer must be plunged in the WBC/HGB chamber and the temperature must be note as soon as possible.

- Run 2 diluent primes («SERVICE» menu, «PRIME» sub menu, «DILUENT» selection).
- Enter the «CALIBRATE» menu. Note the temperature of the diluent and type in the value (if it is different from the previous on the instrument).
- Press  to save the new value.

### 3 - Temperature :

- 1 - **Visible** : Press  to validate the temperature display.  
The value displayed is the calibrated temperature.
- 2 - **Invisible** : Press  to cancel the temperature display.  
(Temperature invisible by default)

## IV - Run mode

The instrument must be configured in the «USER» mode (configured by default in the «USER « mode), the «CONTROL» mode intends for a factory use.

## V - Reagent pack

This function is used to update the instrument from a bottle mode to a pack mode.

## VI - Serial #

Displays the instrument serial number.

## VII - Cycle #

Displays the :  
Startup number since the first use of the instrument.  
Stand by number since the first use of the instrument.  
CBC number since the first use of the instrument.

## IMPORTANT

The startup, stand by and CBC numbers are adjustable in this menu but it is mandatory to keep the initial values (useful for maintenance schedules).

## VIII - Burn-in

This function which allows the burn-in of the instrument is intended for a factory use.

# MECHANIC FUNCTION



Hématologie

08/04/98

## • CONCERNS

- Sensor replacements
- Needle motion check
- Carriage motion check
- Liquid syringe motion check
- Vacuum/waste syringe motion check
- Valve operation check
- LCD contrast adjustment
- Piercing mechanism check (**MICROS 60 CT**)

## • REQUIRED TOOLS

- Hexagonal keys
- Felt-pen

## • REQUIRED PRODUCTS

- None

## • INTERVENTION TIME

- 30 minutes

## • FREQUENCY

- See maintenance chart table.

## • SPECIFIC KIT OR CONSUMABLES


- Vacuum/waste syringe sensor : XBA 319 AS
- Liquid syringe sensor : XBA 319 AS
- Carriage and needle sensors : XBA 250 A

RAS 173 A Ind.A

M I C R O S 6 0

• PROCEDURE

The MECHANIC function arborescence is as follows :



From the main menu, enter the 4 - SERVICE menu and move to 5 - MECHANIC sub menu.

**1 - Sensors**

Press  ;

<b>NEEDLE SENSOR</b>	<b>0</b>
<b>CARRIAGE SENSOR</b>	<b>0</b>

should be displayed on the LCD screen.

Manually raise the needle support in the upper position (see diag 1)

The «0» should switch to «1». This indicates the correct operation of the needle sensor (diagram 1).

If nothing happened try to move the needle up and down again to get the commutation from «0» to «1».

If the test is still wrong, check the correct connection of the sensor on the connector J7 of the mother board.

• **Replacing the needle sensor**

Switch off the instrument.

Use a felt-pen to mark the sensor position.

Unscrew the 2 sensor fixation screws and remove the cell holder .

Disconnect the wire from the J7 connector (See diagram 2).

Replace the previous sensor by a new one and reassemble in the reverse order.


**CAUTION**

**The new sensor must be installed exactly on the same position.**


Proceed as described in **RAS 172 A (MICROS 60 OT)** or **RAS 176 A (MICROS 60 CT)** «Sensor adjustment» to check the correct position of the sensor.

MICROS 60 OT

MCROS 60 CT



Diag.1



Diag.2

• **Replacing the carriage sensor**

Use the 1 - SENSOR test to check the carriage motion detection :

Move rightward the needle carriage in order to perform the sensor detection (diagram 3).

To replace it, note its position with a felt-pen, unscrew the cell fixation screws and remove the cell holder.

Disconnect the wire from the connector J8 (See diagram 2)

Replace the previous sensor by a new one and reassemble in the reverse order.

## CAUTION

The new sensor must be installed exactly on the same position.

Proceed as described in **RAS 172 A (MICROS 60 OT)** or **RAS 176 A (MICROS 60 CT)** « Sensor adjustment » to check the correct position of the sensor.

MICROS 60 OT

MICROS 60 CT



*Diag.3*

### 2 - Liquid syringes motion check

Run a 4 - LIQUID SYRINGE cycle : The syringes (diag 4) are moved upward and downward once.



If the sensor (diag 5) detection is not correct a message «ERROR : SYRINGE MOTOR» is displayed.

Three parts are involved :

- The motor ———> Check the connection on the board (connector J3)
- The syringes ———> Check the motion up/down
- The sensor ———> Check the connection on the board or replace it.

SYRINGE FRONT SIDE VIEW

SYRINGE BACK SIDE VIEW



*Diag.4*

*Diag.5*

### • Replacing the sensor :

- Switch off the instrument.
- Unscrew the 2 screws (Diag.5) and disconnect the wire from the connector J9 (Diag.2)
- Replace the previous sensor by a new one and re-install in the reverse order.
- Switch on the instrument.
- Run a 4 - LIQUID SYRINGE cycle again to control the correct detection.

# MECHANIC FUNCTION

## 3 - Vacuum/waste syringe motion check

Run a 5 - PRESSURE SYRINGE cycle : The syringe piston moves upward and downward once.

If the sensor (Diag.6) detection is not correct a message «ERROR : PRESSURE SYRINGE MOTOR» is displayed.

Three parts are involved :

- |               |       |  |
|---------------|-------|--|
| - The motor   | ————> | Check the connection on the board (connector J4) |
| - The syringe | ————> | Check the motion up/down                         |
| - The sensor  | ————> | Check the connection on the board or replace it. |



**Diag.6**

### • Replacing the sensor :

- Switch off the instrument.
- Unscrew the 2 cell fixation screws (Diag.6) and disconnect the wire from the connector J10 (Diag.2)
- Replace the previous sensor by a new one and re-install in the reverse order.
- Switch on the instrument.
- Run a 5 - PRESSURE SYRINGE cycle again to control the correct detection.

## 4 - Valve operation check

Run a 6 - VALVES cycle. The number of the valves from 1 to 13 is displayed and each valve is activated once.

Check their correct operations. If a valve is suspected proceed as described in the procedure : RAS 170 A.

## 5 - LCD contrast adjustment


Enter the sub menu «7 - CHG. CONTRAST». Use the arrows to modify the contrast and press



to validate the adjustment.

### NOTE

If the LCD screen happens to be unreadable because of a bad contrast adjustment it is possible to get back a correct contrast pressing at the same time both keys :



## 6 - Piercing mechanism check (MICROS 60 CT)

Enter the sub menu 9 - PIERCING. Follow the instructions : «PLEASE CLOSE TUBE HOLDER DOOR» : a piercing operation is simulated.

This function gives the sample tube holder position and the height of the needle in its lower position. See **RAS 172 A (MICROS 60 OT)** or **RAS 176 A (MICROS 60 CT)**

### NOTE

If the sample tube holder has been removed the following message is displayed : «ERROR : NO SAMPLE TUBE HOLDER».  
If the sample tube holder has been turned in between two piercing positions (it means that the tube is not in front of the piercing needle) the following message is displayed : «ERROR : TUBE HOLDER POSITION».



# DRAIN DETECTION



Hématologie

05/10/98

## • CONCERNS

- Drain detection sensor adjustment/replacement

## • REQUIRED TOOLS

- Voltmeter
- Flat screw driver
- 5ml syringe

## • REQUIRED PRODUCTS

- None

## • INTERVENTION TIME

- 15 min

## • FREQUENCY

- On request or once a year

## • SPECIFIC KIT OR CONSUMABLES

- Drain detection sensor : XBA 199 A.

RAS 174 A Ind.B

M I C R O S 6 0

• PROCEDURE

The Vacuum/Waste syringe and the chambers drains are controlled by an infrared sensor located below the syringe.

Controls carried out by the cell during a cycle are as follows :

- *Control of the correct operation of the cell*

At the first chamber drain, air must be detected in the cell within defined timeout, and followed by liquids.

If this switch "air-liquid" has been successful, the adjustment of the cell (see below) is validated. If not, the cycle is stopped and the following message is triggered : "sensor error or diluent empty" (check the connection or the adjustment of the cell).

- *Control of the chamber drains*

A first measure controls that liquids circulate in the cell during a defined timeout.

The second measure checks that air has replaced liquid.

If so, the drain phasis is validated.

If not, the following message is displayed : «.....» (the instrument carries on the current cycle)


It means that the sensor always detects liquids.

(Check the watertightness of the syringe.)

• **Drain detection sensor adjustment**

- Raise the piston up and press the valve <2>.

- Manually perform a syringe drain pressing the valve <5> and pulling down the syringe piston (Diag.1).



**Diag.1**


- Make sure the cell is perfectly drained (no bubble).

- Connect a voltmeter between the ground **TP30** and **TP6** (See Diag.2) and adjust the voltage to **4,5V ± 0,3V** by means of **R6**.


# DRAIN DETECTION

Fill the syringe up with diluent and connect it on the detection cell as shown on diagram 3.

## MICROS 60 MOTHER BOARD



Diag.2



Diag.3

- Push diluent through the sensor and check that the voltage falls down below **1Volt** on the voltmeter.

## NOTE

Check the commutation from 4,5Volts to  $\approx$  0Volt once again pushing and drawing alternately liquid and air through the sensor by means of the syringe.

- If nothing happens switch off the instrument and disconnect the connector from J13.
- Replace the drain detection sensor by a new one.
- Switch on the instrument
- Carry out the new sensor adjustment as described above.

## NOTE

However it is possible to control the correct operation of the sensor by means of the LED located next to **R6**. Indeed the LED should be lit when liquid circulates through the sensor and switched off with air.

# PCB VOLTAGE CHECKS/ ADJUSTMENTS



Hématologie

01/10/98

## • CONCERNS

- Voltage supply check
- Aperture voltage check
- RBC threshold check/adjustment
- WBC threshold check/adjustment
- PLT threshold check/adjustment
- HGB blank voltage check.
- Stepper motor voltages adjustment

## • REQUIRED TOOLS

- Voltmeter
- Flat screw driver

## • REQUIRED PRODUCTS

- None

## • INTERVENTION TIME

- 20min

## • FREQUENCY

- See maintenance chart table.

## • SPECIFIC KIT OR CONSUMABLES

- None

RAS 175 A Ind.B

M I C R O S 6 0

# PCB VOLTAGE CHECKS/ADJUSTMENTS

## • PROCEDURE

### 1 - RBC, PLT, WBC threshold checks/adjustments

- Ground on TP31.
- Adjust the thresholds according to the below chart table :

THRESHOLDS	TEST POINTS	VOLTAGE	POTENTIOMETERS
WBC	TP 10	350 mV $\pm$ 7	R68
RBC	TP13	400 mV $\pm$ 7	R75
PLT	TP16	180 mV $\pm$ 3	R82

(See Diag.1)

### 2 - Aperture voltage check

- See RAS 172 A : **MICROS 60 OT** Technician function
- RAS 176 A : **MICROS 60 CT** Technician function

### 3 - Voltage supply check

- See RAS 171 A

### 4 - HGB blank voltage check

- See RAS 172 A : **MICROS 60 OT** Technician function
- RAS 176 A : **MICROS 60 CT** Technician function

### 5 - Stepper motor voltage checks/adjustments


- Ground on TP31.
- Adjust the motor voltages according to the below chart table :

MOTORS	TEST POINTS	VOLTAGES	POTENTIOMETERS
Air syringe	TP36	2.50 V $\pm$ 0.05 V	R195
Liquid syringe	TP35	2.00 V $\pm$ 0.05 V	R186
Horizontal carriage	TP37	1.50 V $\pm$ 0.05 V	R177
Vertical carriage	TP38	1.00 V $\pm$ 0.05 V	R168

(See Diag.1)

# PCB VOLTAGE CHECKS/ADJUSTMENTS

MICROS 60



Diag.1

# TECHNICIAN FUNCTION



Hématologie

16/06/98

## • CONCERNS

- |   |  |
|---|--|
| 1 - Version display.  | 3 - Temperature sensor adjustment  |
| 2 - Adjustments :<br>HGB photometer calibration, Aperture voltage<br>Vacuum check, WBC gain, RBC & PLT gain<br>Sensor, Needle height, Needle motion<br>Bubbling | 4 - Run mode<br>5 - Reagent pack<br>6 - Serial number<br>7 - Cycle number<br>8 - Burn-in |

## • REQUIRED TOOLS

- Flat screw driver
- thermometer
- Barflex
- Hexagonal keys
- Voltmeter

## • REQUIRED PRODUCTS

- WBC latex : LAD 001 AS
- RBC and PLT latex : LAD 002 AS
- Soft paper
- Flat piece of stiff plastic

## • INTERVENTION TIME

- 60 min

## • FREQUENCY

- See maintenance chart table.


## • SPECIFIC KIT OR CONSUMABLES

- None

MICROS 60 CT

RAS 176 A Ind.A

• PROCEDURE



**MICROS 60 CT**



# TECHNICIAN FUNCTION

Once entering the «SERVICE» sub menu, move to «TECHNICIAN FUNCTIONS» and press



Enter

. A specific password (421) is required to enter the sub menus.

Move the cursor by means of  and  and choose the required menus pressing

the  key.

**I - The version number is displayed.**

**II - Adjustments :**

## 1 - HGB photometer calibration

- Dismantle the WBC/HGB chamber cover.
- Check the general cleanliness of the WBC chamber/spectrophotometer assy.
- Re-install the chamber cover.

### NOTE

If the WBC/HGB chamber has been dismantled previously make sure no liquid has flown in between the spectrophotometer and the chamber. Clean the inner surfaces of the spectrophotometer as well as the chamber. Reassemble the assy and tighten the two screws to the following torque : 400mN.m (see RAS 169 A : Chamber maintenance)

- Run the CAL PHOTOMETER function (selection 1 of the «ADJUSTMENTS» menu) : diluent is delivered to the WBC/HGB chamber twice.

An HGB channel is displayed on the LCD screen :

VALUE  
XXX

### IMPORTANT


- 1 - The HGB photometer calibration must be done 20min at least after the instrument has been switched on.
- 2 - This adjustment must be done with the WBC chamber cover installed!!!

- By means of **R97** (See diagram 1) adjust the HGB channel according to the room temperature using the chart table given on next page.

### NOTE

After 40 seconds approximately, the function is automatically exited.

- Run the CAL PHOTOMETER function again to verify the adjustment.



Diag.1


ROOM TPT (°C)	CHANNEL		
	Mini.	Nominal	Maxi.
15	240	245	250
16	240	245	250
17	239	244	249
18	238	243	248
19	237	242	247
20	236	241	246
21	235	240	245
22	234	239	244
23	234	239	244
24	233	238	243
25	232	237	242
26	231	236	241
27	230	235	240
28	229	234	239
29	228	233	238
30	228	233	238
31	227	232	237
32	226	231	236
33	225	230	235
34	224	229	234
35	223	228	233


# TECHNICIAN FUNCTION

## 2 - Aperture voltage

- Once entering the «APERTURE VOLTAGE» menu, connect a voltmeter between the ground (TP30) and TP19.
- Check that the value is **60V +2.8V -1.5V** (Diagram 1). The aperture voltage is not adjustable.
- Press any key to escape.

## 3 - Vacuum check



- Disconnect the tube from the vacuum/waste syringe coming from the valve  (see diagram 2).

- Follow the instructions given on the LCD screen :

«PLEASE PLUG BARFLEX ON AIR SYRINGE» (On the free nipple).

«CHECK PRESSURE : -200mB ± 10mB» (The piston has raised in order to create a vacuum in the syringe body).

- Check the stability of the vacuum during 30 secondes : The vacuum drop down must be  $\leq 2$  mbar.

### Diag.2

- If the results are not correct check the O ring and the tubing watertightness.

«PLEASE RECONNECT PREVIOUS TUBE» (disconnect the Barflex and replug the tube instead).

## 4 - WBC adjust

- Put the WBC latex to mix on a Vortex during 1min or shake thoroughly

## IMPORTANT

As the WBC gain is a factory adjustment it is mandatory not to readjust it without taking the following precautions :

Carry out previously an autoconcentrated cleaning to make sure of the cleanliness of the WBC counting circuit.

If necessary clean the WBC chamber aperture as described in RAS 169 A.

Make sure the Latex has been thoroughly mixed before.

- Run a blank cycle to check the cleanliness of the instrument.
- Enter the «WBC ADJUST» sub menu.

## TECHNICIAN FUNCTION



- Enter the «WBC ADJUST» sub menu and close the door of the piercing mechanism (See Diag 3) : A CBC's cycle starts.

During the cycle measuring phasis (around 1 minute) the **Lymphocyte** and the **Granulocyte** volumes are displayed on the screen every 3 seconds as shown below :

Diag.3

LYM <57 +/- 1> 57	GRA <180 +/- 2> 180
----------------------	------------------------

- Wait for several results to be displayed and check the stability of both values.


## IMPORTANT

It is mandatory not to operate the gain adjustment as long as the lympho and granulo values are not stable.

After one minute the last volume values displayed on the screen are printed out as well as the WBC, PLT and RBC histograms and the CBC's results. Check that these printed values correspond to the following target values :

- LYM =  $57 \pm 1$
- GRA =  $180 \pm 2$

- If not rerun a «WBC ADJUST» cycle after having previously mixed the latex vial again. By means of R74 (see diagram 4) adjust the volumes to the target values during the measuring phasis.



Diag.4

# TECHNICIAN FUNCTION

## NOTE

Both sub menus 2 - # LYM < > and 3 - # GRA < > allow the technician to change the Latex target values if the latex run on the instrument is different from the one recommended above.

### 5 - RBC PLT adjust

- Put the RBC and PLT latex to mix on a Vortex during 1min or shake thoroughly

## IMPORTANT

As the RBC/PLT gain is a factory adjustment it is mandatory not to readjust it without taking the following precautions :

Carry out previously an autoconcentrated cleaning to make sure of the cleanliness of the RBC/PLT counting circuit.

If necessary clean the RBC/PLT chamber aperture as described in the procedure RAS 169 A.

Make sure the Latex has been thoroughly mixed before.

- Run a blank cycle to check the cleanliness of the instrument.
- Enter the «RBC PLT ADJUST» sub menu.
- Present the vial of Latex to the open probe as shown on diagram 3 and press the sampling bar located behind the sampling needle : the needle directly delivers the latex sample in the RBC chamber dilution (for a usual analysis cycle the sample is first delivered to the mixing chamber) and a measuring phase begins.
- During the cycle measuring phase (around 1 minute) the **Platelet** and the **Red Blood cell** volumes are displayed on the screen every 3 seconds as shown below :

RBC <74 +/- 1> 74	PLT <59 +/-1> 59
----------------------	---------------------

- Wait for several results to be displayed and check the stability of both values.

## IMPORTANT

It is mandatory not to operate the gain adjustment as long as the platelet and RBC values are not stable.

## NOTE

From the latex lot # 980311 included, balls having a different size, a drift of MPV peak has been noticed, i.e. a modification of the PLT gain target value : it becomes 64 instead of 59.

The program default value will be modified in the next MICROS version.

Both sub menus 2 - MCV < > and 3 - MPV < > allow the technician to change the Latex target values. If the lot #  $\geq$  980311, modify the target values and proceed the same way to adjust the PLT gain.

MICROS 60 CT

- If not, rerun a «RBC PLT ADJUST» cycle after having previously mixed the latex vial again.
- Adjust the volumes to the target values during the measuring phasis by means of the following potentiometers (see Diagram 4) :

- RBC gain with R81
- PLT gain with R87

### 6 - Sensor adjust

#### • Needle sensor

If the needle sensor (diag 6) has been replaced by a new one or dismantled for any reason, it is mandatory to re-position it at the right height. Proceed as following :




Diag.5

- Install the piece of plastic (diagram 5) underneath the needle rinsing block.

- Once entering the menu «SENSOR ADJUST», enter the sub menu «NEEDLE SENSOR».
- Push the sampling needle downward until it stops against the piece of plastic and press any key in order to raise the needle back in its upper position.

The current number of steps, the mini and maxi values are displayed as well as the way to remove the sensor (shown by an arrow) if the current value is out of ranges (see below).

CURRENT : XXX      MIN : 422      MAX : 432      ▲



Diag.6

- For a current number of steps out of ranges, loosen the 2 cell fixation screws (diagram 6) and gently move the sensor
  - upward if the current value is too low
  - downward if the current value is too high.

- Tighten the screws and rerun a «NEEDLE SENSOR» cycle. Check that the current value is correct.

### NOTE


10 steps correspond to around 1 mm. The target number of steps is  $70 \pm 5$

# TECHNICIAN FUNCTION

## • Carriage sensor

This function allows the adjustment of the carriage sensor (diagram 8) position.

- Proceed as described below :




Diag.7

- Install the piece of plastic against the left side of the tube holder compartment (see Diagram 7).

- Once entering the menu «SENSOR ADJUST» (Selection 6), enter the sub menu «CARRIAGE SENSOR».

Move the carriage on the left until the piercing needle stops against the piece of plastic. Press any key : the carriage comes back in its initial position.

The current number of steps, the mini and maxi values are displayed as well as the way to remove the sensor if the current value is out of ranges (see below).



Diag.8

- For a current number of steps out of ranges, unloosen the 2 screws (see diagram 8) and gently move the sensor

- towards the right if the current value is too low

- towards the left if the current value is too high.

## NOTE


10 steps correspond to around 1 mm. The target number of steps is  $350 \pm 5$

- Tighten the screws (see diagram 8) and rerun a «CARRIAGE» cycle. Check that the current value is correct.

**7 - Position adjustment****• Needle height adjustment**

When replacing a needle, it is mandatory to adjust the height of the needle in the chambers. Proceed as follows :

- Remove the WBC/HGB chamber cover.
- Enter the menu 7 - POSITION ADJUST / 1 - NEEDLE / 2 - ADJUSTMENT.
- Manually pull down the needle until it comes into contact with the edge of the RBC chamber (Diag. 9).
- Press a key to continue : the needle moves up to the initial position ; the value is stored.

*Diag.9**Diag.10***• Needle motion adjustment**

The needle position in the WBC chamber can be adjusted as follows :

- Enter the menu 7 - POSITION ADJUST / 2 - CARRIAGE / 2 - ADJUSTMENT.
- Position the tool (P/N : GBC 218 A) over the RBC and WBC chambers (Diag 10).
- Manually lower the needle into the WBC chamber.
- Press a key to continue : the needle moves up, and the carriage comes back to the initial position ; the value is stored.
- Carry out a 3 - CHECK cycle : the needle comes down to the WBC chamber.
- Check that the needle is centered in the hole.

If not, enter the 2 - CARRIAGE / 1 - CARRIAGE < > menu ; the display shows the current number of steps carriage motion.

If the needle goes too far on the right, add 1 step to the current value for 0.1mm.

If the needle is too much on the left, decrease the current value of 1 step for 0.1mm.

**CARRIAGE ? :**  
**CURRENT : 893**

**EXIT : ESC**  
**SAVE : ENTER**

Carry out a 2 - CARRIAGE / 3 - CHECK cycle again to control the needle position.



# TECHNICIAN FUNCTION

## 8 - Bubbling adjustment

An overflow protection tank is installed on the drain circuit of each chamber. This one prevents from polluted liquid overflow during bubbling phasis (Diag .11).



Diag.11

Two bubbling phasis are adjustable :

- "BUBBLING 1" is the first dilution (WBC/HGB chamber) bubbling value.
- "BUBBLING 2" is the second dilution (WBC/HGB chamber + LYSE) value and RBC chamber bubbling value.

Both values correspond to a number of steps carried out by the waste/vacuum syringe.

Default values are        BUBBLING 1 : 175  
                                     BUBBLING 2 : 120

## CAUTION

These values are factory adjusted (and may be different from the default values shown above) and should be modified only when hematologic erroneous results are given by the instrument : If values are too important, liquid overflows can occur or if bubbling is too low homogeneity of the dilution can be decreased.

Ranges :    150 < BUBBLING 1 < 200  
               80 < BUBBLING 2 < 140


To modify the bubbling values, enter the menu :

- 2 - ADJUSTMENTS / 8 - BUBBLING ADJUST. / BUBBLING 1 < >
- 2 - ADJUSTMENTS / 8 - BUBBLING ADJUST. / BUBBLING 2 < >

and type in new step value.

Carry out a 3 - CHECK to control the adjustment.

8 - Needle sampling position





Diag.13

- Move the tube holder to one of the four sampling positions.

- Move the cursor to the function 9 - SAMPLING POSITION / 1 - NEEDLE ADJUSTMENT.

- Close the sample holder door. The sampler assy moves to the upper position (except in the Needle 6 position).

- Manually push the sampling needle to the bottom of the required sampling position and press  : The MICROS 60 CT adjust automatically the needle depth whatever the tube holder position can be. The sampler assy comes back to the initial position and the sampler door is opened.

**WARNING !** 


OEM instruments equipped with specific sample tube holders (see RAS 198 A procedure) must be used with the specific tubes or blood controls they were intended for.

- Turn the tube holder to another needle position and carries out the same procedure to adjust the needle depth.

It is possible to enter directly the required number of steps for each sampling position. Proceed as following :

Enter the sub menu that corresponds to the number of the needle : the current number of steps is displayed. Enter the new value.

Increase the number of steps to move the needle deeper or decrease the value to raise the needle.

Confirm the new value with  . The minimum and maximum step values are as follow :

NEEDLE	NUMBER OF STEPS		
	MINI.	DEFAULT	MAXI.
1	1	788	1100
2	1	661	1100
3	1	612	1100
4	1	948	1100
5	1	1003	1100
6	1	845	1100

# TECHNICIAN FUNCTION

## \* Check

The CHECK function allows to check the piercing operation on each sampling position. From the SAMPLING POSITION menu, move the cursor to the function <9> CHECK and press ENTER.


Select the required position on the tube holder and close the sample door. The piercing cycle is carried out and the number of steps for this position is displayed :

NEEDLE 2 CURRENT : 635	STANDARD : 630
---------------------------	----------------

## III - Temperature adjustment

When entering the «temperature adjustment» menu the following sub menus are displayed :

### 1 - Uncal. temperature

When pressing the  key the sensor temperature value **uncalibrated** is displayed. This value should be close to the diluent temperature.

### 2 - Calibrate :


The temperature must be calibrated according to the diluent temperature :

- Plunge a thermometer directly into the diluent container and leave it for a while until stabilization.


**NOTE** For a pack equipped instrument, the thermometer must be plunged in the WBC/HGB chamber and the temperature must be note as soon as possible.


- Run 2 diluent primes («SERVICE» menu, «PRIME» sub menu, «DILUENT» selection).

- Enter the «CALIBRATE» menu. Note the temperature of the diluent and type in the value (if it is different from the previous on the instrument).

- Press  to save the new value.

### 3 - Temperature :

**1 - Visible :** Press  to validate the temperature display.  
The value displayed is the calibrated temperature.

**2 - Invisible :** Press  to cancel the temperature display.  
(Temperature invisible by default)

## IV - Run mode

The instrument must be configured in the «USER» mode (configured by default in the «USER « mode), the «CONTROL» mode intends for a factory use.

## V - Reagent pack

This function is used to update the instrument from a bottle mode to a pack mode.

**VI - Serial #**

Displays the instrument serial number.

**VII - Cycle #**

Displays the :      Startup number since the first use of the instrument.  
                         Stand by number since the first use of the instrument.  
                         CBC number since the first use of the instrument.

**IMPORTANT**

The startup, stand by and CBC numbers are adjustable in this menu but it is mandatory to keep the initial values (useful for maintenance schedules).

**VIII - Burn-in**

**CAUTION**

This function which allows the burn-in of the instrument is intended for a factory use.

# LX300 PRINTER SETUP



Hématologie

01/04/98

## • CONCERNS

- 1 - Configuration
- 2 - Control panel
- 3 - Control LEDS and keys
- 4 - Printer description

## • REQUIRED TOOLS

- None

## • REQUIRED PRODUCTS

- None

## • INTERVENTION TIME

- 15 minutes

## • FREQUENCY

- On request

## • SPECIFIC KIT OR CONSUMABLES


- None

RAS 177 A Ind.A

M I C R O S 6 0

- PROCEDURE

**NOTE** On the MICROS 60, it is necessary to select the printer format RESERVED 1 of the "PRINTER" menu (function 4) accessible through the "OPTIONS" menu (function 5 of the main menu) then "RESULTS" (function 1).



*Diag.1*

# LX300 PRINTER SETUP

## 1 - Printer configuration :

The printer configuration is printed out when pressing the key <FONT> when the printer is switched ON. The configuration used for the MICROS 60 is the factory configuration :

```

*****
      << Current settings >>
Character spacing 10 cpi
Shape of zero    0
Skip-over-perforation Off
Character table  PC 437
Auto line feed   Off
Page length      12 inches
Auto tear off    Off
Tractor          Single
Interface        Auto selection (10 sec.)
Bit rate         9600 bps
Parity           None
Data length      8 bit
ETX/ACK         Off
Software         ESC/P
Auto CR         Off
    
```

Diag.2


Each parameter can be modified by the corresponding parameter chart. Each chart is accessible using the keys <PAUSE>, <FONT> and <LF/FF> according to the control LED combinations :

Table B				
Lights			Setting	Go to submenu
1	2	PAUSE		
BLINKS	OFF	OFF	Character spacing	Table C
BLINKS	ON	OFF	Shape of zero	Table D
OFF	BLINKS	OFF	Skip-over-perforation	Table E
ON	BLINKS	OFF	Character table	Table F
BLINKS	BLINKS	OFF	Auto line feed	Table G
BLINKS	OFF	ON	Page length	Table H
BLINKS	ON	ON	Auto tear off	Table I
OFF	BLINKS	ON	Tractor	Table J
ON	BLINKS	ON	Interface	Table K
BLINKS	BLINKS	ON	Bit rate	Table L
OFF	OFF	BLINKS	Parity	Table M
BLINKS	OFF	BLINKS	Data length	Table N
ON	OFF	BLINKS	ETX/ACK	Table O
ON	ON	BLINKS	Software	Table P
BLINKS	ON	BLINKS	Auto CR	Table Q

Diag.3

MICROS 60

2 - Control pannel : (See Diag.4)



Diag.4

The control pannel keys allow the user to set up the main functions of the printer : paper advance, paper ejection, and font selection. Control LEDS indicate the printer status.

### 3 - Control LEDS and keys :

**LED PAUSE** : The orange LED PAUSE lights when the printer stops printing. During each power ON, this LED blinks for few seconds and 4 audible beeps occur. When the printer runs out of paper, the LED blinks and 3 audible beeps occur. This LED lights also when the paper is in its tear off position. When a problem occurs, this LED lights ON and 5 audible beeps occur.

**LEDS FONT 1 and FONT 2** : These 2 green LEDS indicate the selected font. Refer to the printer user's manual to select the font.

**Key FONT** : During normal operation, the FONT key allows the font selection. For each pressure on this key, the selection is modified. Refer to the printer user's manual to select the font. When this key is pressed during the printer power ON, the printer setup menu is entered.

**Key LF/FF** : During normal operation, a quick pressure on this key allows a ligne feed of the paper. Keep the pressure on this key to feed a whole page. This key can be used to load or eject the paper. When this key is pressed during the printer power ON, the printing test starts.

**Key PAUSE** : When this key is pressed during the printing, the printout stops. Press again on this key to restart the printout.


**PAPER PARK** : If Z folded paper is used, the paper can be driven to its parking position when pressing simultaneously on the keys LF/FF and FONT.

**MICRO ADJUST** : This function allows to adjust the loading paper position. See the user's manual for details.




# LX300 PRINTER SETUP

## 4 - Printer description :



Diag.5



Diag.6

M  
I  
C  
R  
O  
S  
6  
0

# LIQUID SYRINGES



Hématologie

16/06/98

## • CONCERNS

Lyse dispenser O ring replacement  
Diluent dispenser O ring replacement  
Sampling needle dispenser O ring replacement  
Lubrication of the liquid syringes

## • REQUIRED TOOLS

- Hexagonal keys
- Dynamometric screw driver : A302 : MAG 019 A  
A301 : MAG 020 A

## • REQUIRED PRODUCTS

- Silicone grease : LAM 004 A
- Soft paper
- Grease for mechanical assemblies : XEA 381 AS

## • INTERVENTION TIME

- 30 min

## • FREQUENCY

- Once a year or on request

## • SPECIFIC KIT OR CONSUMABLES

- O ring kit : XEA 328 AS

RAS 178 A Ind.A


M I C R O S 6 0

### • PROCEDURE

**NOTE** In case of a leak on one of the three syringes it is recommended to replace the entire dilution block. If not proceed as following :


- Pull the piston assy in the upper position and clamp the diluent and lyse tubing as described on the diagram 1.

- Disconnect the diluent/lyse/sampling tubings from the 3 syringes and the tube on the sampling syringe side.



Diag.1

- Unscrew the 4 fixation screws (2) , and remove the dilution block (1) (Diagram 1).




Diag.2

- Unscrew the 6 screws (1) and (2) in order to remove the body cover (3) (Diagram 2).

- Pull out the pistons 1, 2 and 3 from the body with their respective o ring still around (Diagram 3).

- Replace the lyse and diluent O rings by new ones. Check the cleanliness of the piston and of the syringe bodies. If necessary clean with a soft paper.

# LIQUID SYRINGES



Diag.3

## • Lubrication



- Spread a little amount of silicone grease between two fingers and apply a very thin film of grease on the 2 new sampling syringe O rings. Replace the old ones (4) .

- Reinstall the dilution block assy in the reverse order.

Use the A302 dynamometric screw driver to tighten the screws (1) to **700mN.m** (Diagram 2).

Use the A301 dynamometric screw driver to tighten the screws (2) to **400mN.m** (Diagram 2).

- Disconnect the diluent and waste inputs located at the rear of the instrument (or remove the reagent pack)

- Move the liquid syringe by hand in order to have an access to the motor gearings.

- Spread a little amount of grease on the gearings and on the piston axis (Diag.4).

- Move by hand the syringe assembly to spread the grease on all parts of the gearings and piston axis.

- Re-install the instrument cover, reconnect the waste and diluent tubes, reconnect the power cable.

- Switch the instrument on and run several priming cycles.

## CAUTION

Place some absorbant paper at the instrument rear connections (diluent and waste) as some liquids may come out when the syringe is pushed.

M  
I  
C  
R  
O  
S  
  
6  
0

# VACUUM/WASTE SYRINGE MAINTENANCE



Hématologie

16/06/98

## • CONCERNS

- O ring replacement

## • REQUIRED TOOLS

- Hexagonal keys
- Dynamometric screw driver A302 : MAG 019 A

## • REQUIRED PRODUCTS

- Silicone grease : LAM 004 A
- Grease for mechanical assemblies : XEA 381 AS.

## • INTERVENTION TIME

- 20 min

## • FREQUENCY

- See maintenance chart table

## • SPECIFIC KIT OR CONSUMABLES

- O ring kit : XEA 328 AS

RAS 179 A Ind.A

M I C R O S 6 0


# VACUUM/WASTE SYRINGE MAINTENANCE

## • PROCEDURE

- Manually pull the syringe piston in order to freed it from the syringe body.

## NOTE

The syringe has to be linked to the atmosphere, pressing the valve **2**, to pull the piston out from the syringe body.  
Turn the cylindrical gearing **1** by hand to help the raising of the piston (See Diagram 1).



Diag.1

# VACUUM/WASTE SYRINGE MAINTENANCE

- Unscrew the fixation screws (5) in order to remove the syringe body.
- Unscrew the O ring tightening screws (3) and remove the O ring (4).
- Spread a little amount of silicone grease between two fingers and apply a very thin film of grease on a new O ring.
- Reinstall in the reverse order. Apply the following torque to the screws (5) : 700 mN.m.

## CAUTION

Push the piston back inside the syringe body before tightening the o ring fixation screws (3). Use the dynamometric screw driver to tighten the screws (3) to 400 mN.m

- Check the watertightness of the syringe running a «CHECK PRESSURE» cycle (see procedure : RAS 172 A for OT or RAS 176 A for CT).
- Run cycles and check for correct operations.

### • Lubrication



Diag.2

- Switch off the instrument and remove the power cable.
- Disconnect the diluent and waste inputs located at the rear of the instrument.
- Using a small and flat screwdriver, spread a little amount of grease on the gearings of the air syringe reductor plate (Diag.2).



Diag.3

- Spread a little amount of grease on the cogs of the piston axis (Diag.3).

M  
I  
C  
R  
O  
S  
  
6  
0

## VACUUM/WASTE SYRINGE MAINTENANCE



- Move by hand the piston axis up and down in order to spread the grease all around the gearings and along the axis (Diag.4).

*Diag.4*



# INSTRUMENT LANGUAGE



Hématologie

01/04/98

## • CONCERNS

Changing the instrument language

## • REQUIRED TOOLS

Pair of pliers

## • REQUIRED PRODUCTS

None

## • INTERVENTION TIME

5 minutes

## • FREQUENCY

On request only


## • SPECIFIC KIT OR CONSUMABLES

None

**RAS 180 A Ind.A**

**MICROS 60 S**


• PROCEDURE



# INSTRUMENT LANGUAGE

- Replug the jumpers on E1 according to the wished language :

## JUMPER CONFIGURATION ON E1



- Switch on the instrument.

M  
I  
C  
R  
O  
S  
  
6  
0

# REAGENT PACK



Hématologie

10/04/98

M I C R O S 6 0

## • CONCERNS

- Replacement of the waste connector O ring
- Replacement of the diluent/clean/lyse O rings

## • REQUIRED TOOLS

- Pair of pliers
- Torx keys

## • REQUIRED PRODUCTS

- None

## • INTERVENTION TIME

- 15 min

## • FREQUENCY

- O ring replacements : 1/year


## • SPECIFIC KIT OR CONSUMABLES

- O rings : FAA 036 A

RAS 181 A Ind.A

- PROCEDURE


- Disconnect the pack if this one is still connected.
- Dismantle the reagent pack front panel unscrewing the four torx screws as shown on the diagram 1.



*Diag.1*



*Diag.2*



*Diag.3*

- Use a pair of pliers to cut the O rings of the reagent connectors (diag 2) .
- Replace the O rings by new ones.
- Replace as well the waste connector O ring (diag 3).

# BARCODE READER



Hématologie

20/04/98

## • CONCERNS

Barcode reader installation & configuration

## • REQUIRED TOOLS

None.

## • REQUIRED PRODUCTS

None.

## • INTERVENTION TIME

10 minutes.

## • FREQUENCY

On request.

## • SPECIFIC KIT OR CONSUMABLES

Installation kit : XBA 379 AS

RAS 182 A Ind.A

MICROS 60

## • PROCEDURE

### Installation :


- Switch off the instrument.
- Open the instrument cover.
- Install the jumper (included in the XBA 379 AS kit) on **E10** as shown on the diagram 3.

## WARNING !



**Once the jumper installed on E10, the instrument data output receives 5 volts to supply the barcode reader. This voltage should cause damages on computer connection if this one is directly connected on the MICROS 60 data output. It is then mandatory to connect the computer connection only on the cable (DAC 023 AS shown on the diagram 2) intended for it.**

- Connect the BARCODE reader in the DIN plug of the RS adaptor wiring (diag 1 and 2).




Diag.1

### Setup :

- Switch on the instrument and within the 15 first seconds, read from the top to the bottom the 3 barcode labels located on the top left of the page 4 of this procedure. The audible beep occurs after each reading.

- After the 3 labels (the audible signal beeps 5 times) read all the labels from top to bottom and from left to right.

- When the last label is read, the signal beeps 5 times in order to indicate the end of the setup. Check on the test labels located on the last page that the reading is correct.




Diag.2

## NOTE

- In case of difficulties to read 2 consecutive 0, move away the barcode reader from the page after the first reading and then read the second 0.
- To carry out the *read test after setup* the instrument must be configured in a US mode in order to obtain the barcode identification in the ID field (menu 5 - SETUP / 3 - SPECIAL / 7- ID MODE).
- The barcode reader can be configured according to the type of barcode label in use (from the menu 5 - SETUP / 6 - BARCODE) and allows to enable the checksum or not.
- The barcode setup can be printed out by the function 5 - SETUP / 3 - SPECIAL / 5 - PRINT CONFIG.

# BARCODE READER



M I C R O S 6 0

Diag.3



Within 15 seconds after switch on :



(5 beeps)

Setup : default setting



Parity : even



1 bit stop



Code identifier : no



Codabar : yes



No start/stop code



No STF



No C93



C128



∅



∅



Enter



Exit



5 Beeps

# BARCODE READER

## READ TEST AFTER SETUP

### WITH CHECKSUM

EAN 8



1234567

EAN 8



1234567

EAN 13



123456789012

EAN 13



123456789012

### WITHOUT CHECKSUM

C39



TEST

2/5



9076543210

CODABAR



543210

C 128



98765

M  
I  
C  
R  
O  
S  
  
6  
0

## HYDRAULIC CYCLE CHECKUP



Hématologie

16/04/98

### • CONCERNS

Step by step control of the hydraulic cycle.

### • REQUIRED TOOLS

None

### • REQUIRED PRODUCTS

Blood samples

### • INTERVENTION TIME

15 minutes


### • FREQUENCY

On request

### • SPECIFIC KIT OR CONSUMABLES

None

• PROCEDURE




**Diag.1**

**1 - Cycle start condition**

- Needle (1) in the sampling position (diag 1).
- RBC chamber (2) filled with 2.5ml of diluent.
- WBC/HGB chamber (3) filled with 2.5ml of diluent.
- Liquid syringes (4) in standby position.
- Vacuum/waste syringe (5) in the lower position.

**2 - Sampling**

- Aspiration of 10µl of blood sample




**Diag.2**

Causes

The liquid syringes assembly moves down (diag 2) and pulls down the micro sampling syringe (1).

# HYDRAULIC CYCLE CHECKUP

## 3 - Outer sampling needle rinse



Diag.3

### Causes

- The sampling needle moves up.
- The liquid syringes ② (diag 2) send diluent for rinse through the rinsing block
- Polluted diluent is aspirated (from the lower to the upper tube) by means of waste/Vacuum syringe raise.




Diag.4

## 4 - WBC/HGB chamber rinse & HGB blank measure

- Sampling carriage transfer over the WBC/HGB chamber.
- Counting head rinse
- WBC/HGB chamber drain
- Needle motion downward of a few steps (Diag 5)
- Diluent is delivered from the rinsing block

### Causes

- Diluent is delivered by means of the liquid syringes raise
- Aspiration by means of the Vacuum/waste syringe raise (diag 4)
- The liquid syringes move up and a flow of diluent is delivered to the chamber via the outer needle.



Diag.5

MICROS 600T

## HYDRAULIC CYCLE CHECKUP

- WBC/HGB chamber drain (second time)
- Diluent is delivered from the outer needle (second time)
- HGB blank measure (beep triggered)
- RBC and WBC/HGB chamber drains.

### 5 - Dilutions

- Sampling needle moves down to the WBC/HGB chamber
- Injection of 1.7ml of diluent into the WBC/HGB chamber
- + Injection of 10µl of blood sample



#### Causes


Raise of the liquide syringes (diag 6) :

- delivers 0.5ml of diluent from the outer sampling needle
- injects 1.2ml of diluent + blood sample from the inner sampling needle

- Bubbling



Bubbling by means of vacuum/ waste syringe downward motion (Diag 7).




Diag.6




Diag.7

- Sampling needle moves up
- Outer needle short aspiration (dries the needle)
- Sampling needle moves back in the chamber
- Aspiration of 28.3µl of diluted blood (dilution 1/170)



#### Causes

The liquid syringes move down and pull down the micro sampling syringe (diag. 8)



Diag.8

- Sampling needle moves up
- Injection of 0.4 ml of diluent into the WBC/HGB chamber



#### Causes

Raise of the liquide syringes (diag.6) :

- delivers 0.4ml of diluent from the outer sampling needle

# HYDRAULIC CYCLE CHECKUP


- Outer needle aspiration
- Carriage motion over the RBC chamber
- Sampling needle moves down into the RBC chamber
- Lyse injection into the WBC/HGB chamber

Causes

The liquid syringes raise (diag 8) and the syringe ① delivers 0.52ml of lyse via the WBC chamber bottom

**FINAL DILUTION in the WBC CHAMBER :**

1.7ml + 0.4ml diluent + 0.52ml of lyse = 1/260



Diag.8



Diag.9

- Injection of 2.5 ml of diluent into the RBC chamber
- + Injection of 28.3 µl of diluted blood

Causes

The liquid syringes raise (diag 8) :  
 - delivers 0.5ml diluent via the outer sampling needle  
 - injects diluted blood from the inner needle + 2ml of diluent

**FINAL DILUTION in the RBC CHAMBER :**

28.3µl of diluted blood at 1/170 + 2.5ml diluent

$$\frac{1}{170} * \frac{28.33}{2500} = \frac{1}{15000}$$

- Bubbling (diag 10)
- Sampling needle moves up
- Carriage motion over the WBC/HGB chamber
- Counting head rinse



Diag.10

**6 - Counts**

- First counts (beep triggered)
- Counting head rinse
- Second counts (beep triggered)

**NOTE** A third count (C3) is carried out if the difference between first (C1) and second count (C2) is not within acceptable limits :

- WBC :

if C1 or C2 > 3000

C3 is carried out if difference between C1 and C2 > 7%

if Max C1 or C2  $\leq$  3000

C3 is carried out if difference between C1 and C2 > 9%

- RBC :

if C1 or C2 > 16000

C3 is carried out if difference between C1 and C2 > 5%

if Max C1 or C2  $\leq$  16000

C3 is carried out if difference between C1 and C2 > 8%

- PLT :


if C1 or C2 > 400

C3 is carried out if difference between C1 and C2 > 15%

if Max C1 or C2  $\leq$  400

C3 is carried out if difference between C1 and C2 > 20%

- Counting head rinse
- WBC chamber drain
- Diluent injection into the WBC chamber from the outer needle
- RBC chamber drain
- Carriage motion over the RBC chamber
- Diluent injection into the RBC chamber from the outer needle (diag 11)



**Diag.11**

- Carriage & needle motions back to the initial positions
- Results display and printed out



## HYDRAULIC CYCLE CHECKUP



Hématologie

02/06/98

### • CONCERNS

Step by step control of the hydraulic cycle.

### • REQUIRED TOOLS

None

### • REQUIRED PRODUCTS

Blood samples

### • INTERVENTION TIME

15 minutes

### • FREQUENCY


On request

### • SPECIFIC KIT OR CONSUMABLES

None

## Hydraulic cycle checkup

### • PROCEDURE



**Diag.1**

#### 1 - Cycle start condition

- Needle (1) in the sampling position (diag 1 & 2).
- RBC chamber (2) filled with 2.5ml of diluent.
- WBC/HGB chamber (3) filled with 2.5ml of diluent.
- Liquid syringes (4) in standby position.
- Vacuum/waste syringe (5) in the lower position.



**Diag.2**

#### 2 - Cap piercing

- Carriage motion over the WBC chamber.
- WBC/HGB chamber drain

#### Causes

Aspiration by means of the Vacuum/waste syringe raise.

- Diluent injection into the WBC/HGB chamber through the channels (A) and (B) (see diag.4)

Diluent is delivered by means of the liquid syringes raise

- HGB blank measure (beep triggered)
- Carriage return over the piercing device
- Atmosphere is provided inside the tube

#### Causes

Liquid valve #3 is activated


- Cap piercing

Sampling holder rises in the upper position.  
The needle pierces the tube cap

# Hydraulic cycle checkup

## 3 - Sampling

- Sampling needle (diag. 2) moves down to the lower position (inside the tube)
- Aspiration of 10µl of blood sample




### Causes

The liquid syringes assembly moves down (diag 3) and pulls down the micro sampling syringe ①.

Diag.3


- The sampling needle comes back in the upper position
- RBC chamber drain
- Counting head rinse
- Carriage motion over the WBC/HGB chamber
- Sampling holder door opens
- WBC/HGB chamber drain

## 4 - needle rinses



### Causes

- The liquid syringe ② (diag 3) sends diluent for rinse through the outer piercing needle (A) (diag 4) and inner piercing needle (B). The polluted diluent is sent to the WBC/HGB chamber.



Diag.5

Diag.4

## NOTE

The piercing needle inner rinse is equivalent to the sampling needle outer rinse.

- WBC/HGB chamber drain
- Second needle rinses
- WBC/HGB chamber drain
- The sampling needle moves down into the WBC/HGB chamber

## Hydraulic cycle checkup

### 5 - Dilutions

- Injection of 1.7ml of diluent into the WBC/HGB chamber  
+ Injection of 10µl of blood sample



#### Causes

Raise of the liquide syringes (diag 7) :



- delivers 0.5ml of diluent from the outer sampling needle (channel B diag.6)
- injects 1.2ml of diluent + blood sample from the inner sampling needle (channel C)

- Bubbling



Bubbling by means of vacuum/waste syringe downward motion (Diag 8).

- Sampling needle moves up
- Outer needle short aspiration (dries the needle)




Diag.6

Diag.7

Diag.8

- Sampling needle moves back in the chamber
- Aspiration of 28.3µl of diluted blood (dilution at 1/170)

- Sampling needle moves up
- Injection of 0.4 ml of diluent into the WBC/HGB chamber



#### Causes

Raise of the liquide syringes (diag.7) :

- delivers 0.4ml of diluent from the outer sampling needle (channel B)

- Carriage motion over the RBC chamber
- Sampling needle moves down to the RBC chamber

- Lyse injection into the WBC chamber + Bubbling



#### Causes

The liquide syringes raise (diag.7) and the syringe ① delivers 0.52ml of lyse via the WBC chamber bottom

**FINAL DILUTION in the WBC CHAMBER :**  
1.7ml + 0.4ml diluent + 0.52ml of lyse = 1/260

# Hydraulic cycle checkup


- Injection of 2.5 ml of diluent into the RBC chamber
- + Injection 28.3 µl of diluted blood



## Causes

The liquid syringes raise (diag.9) :

- delivers 0.5ml diluent via the outer sampling needle (channel B diag. 6)
- injects diluted blood from the inner needle + 2ml of diluent (channel C)



Diag.9



Diag.10

**FINAL DILUTION in the RBC CHAMBER :**  
 28.3µl of diluted blood at 1/170 + 2.5ml diluent

$$\frac{1}{170} * \frac{28.33}{2500} = \frac{1}{15000}$$

- Bubbling (diag 11)
- Sampling needle moves up
- Carriage motion over the WBC/HGB chamber
- Counting head rinse



Diag.11

MICROS 60 CT

### 6 - Counts

- First counts (beep triggered)
- Counting head rinse
- Second counts (beep triggered)

**NOTE** A third count (C3) is carried out if the difference between first (C1) and second count (C2) is not within acceptable limits :

- WBC :

if  $C1$  or  $C2 > 3000$

C3 is carried out if difference between C1 and C2  $> 7\%$

if  $\text{Max } C1$  or  $C2 \leq 3000$

C3 is carried out if difference between C1 and C2  $> 9\%$

- RBC :

if  $C1$  or  $C2 > 16000$

C3 is carried out if difference between C1 and C2  $> 5\%$

if  $\text{Max } C1$  or  $C2 \leq 16000$

C3 is carried out if difference between C1 and C2  $> 8\%$

- PLT :

if  $C1$  or  $C2 > 400$

C3 is carried out if difference between C1 and C2  $> 15\%$

if  $\text{Max } C1$  or  $C2 \leq 400$

C3 is carried out if difference between C1 and C2  $> 20\%$

- Counting head rinse
- WBC chamber drain
- Diluent injection into the WBC chamber from the outer sampling needle (channel B)
- RBC chamber drain
- Carriage motion over the RBC chamber
- Diluent injection into the RBC chamber from the outer sampling needle (channel B)
- Carriage & needle motions back to the initial positions
- Results display and printed out

# OVERALL MAINTENANCE



Hématologie

05/10/98

## • CONCERNS

Instrument maintenance step by step

## • REQUIRED TOOLS

Hexagonal keys	Voltmeter
Dynamometric screw driver A302, A301, A300	Flat screw driver
Clamps	Barflex
Scalpel	Thermometer
Cutting pliers	Torx keys
Pair of scissors	

## • REQUIRED PRODUCTS

Empty sample tubes	Flat piece of stiff plastic
Silicone grease	Latex WBC, RBC
Soft tissue	Felt pen
Liquid soap	syringe 5ml
Distilled water	Fresh blood samples
Micropipette tip	Calibrator
	Commercial control

## • INTERVENTION TIME

2 h 30

## • FREQUENCY

The yearly maintenance frequencies vary with the instrument output.  
According to the cycle numbers, 3 categories are created :

- < 6000 cycles/year -> 1 overall maintenances/year
- 6000 to 15000 cycles/year -> 2 overall maintenances/year
- > 15000 cycles/year -> 3 overall maintenances/year

## • SPECIFIC KIT OR CONSUMABLES

Spare parts kit : XEA 458 AS  
Needle position tool : GBC 218 A

RAS 191 A Ind.B

M I C R O S 6 0

### A - INSTRUMENT CHECKUP

#### 1 - Reagent check

- ✓ Type of reagent used
- ✓ Expiration dates
- ✓ Levels
- ✓ Pack : Number of cycles left

#### 2 - Operation check

- ✓ Blank cycle control
- ✓ QC control
- ✓ Fresh blood sample run
- ✓ Calibration coefficient checkup
- ✓ Leak control and general cleanliness of the instrument

### B - CLEANING AND MAINTENANCE

#### 1 - Chamber maintenance

- ✓ RBC & WBC/HGB chamber cleaning
- ✓ Aperture check
- ✓ Coaxial O ring replacement
- ✓ Aperture O ring replacement
- Aperture tightening torque : **100mN.m**

RAS 169 A

#### 2 - Liquid syringes

- ✓ Cleaning
- ✓ O ring replacement (Lyse, diluent, Micro syringe)
- Torque values :

RAS 178 A

① torx screws : **400mN.m**

② Hexagonal screws : **700mN.m**



Diag.1

- ✓ Lubrication of the gearings and piston axis

#### 3 - Sampling needle

- ✓ Cleaning
- ✓ Sampling needle O ring replacement
  - CT : O ring holder tightening torque : **700mN.m**
  - OT : O ring holder tightening torque : **100mN.m**
- ✓ Free pulley cleaning and lubrication
  - Pulley tightening torque : **400 mN.m**
- ✓ Piercing needle check/cleaning (CT)

RAS 167 A (CT)  
RAS 168 A (OT)

#### 4 - Air syringe

- ✓ Cleaning
- ✓ O ring replacement
  - O ring tightening torque : **400mN.m**
  - Syringe holding torque : **700mN.m**
- ✓ Lubrication of the gearings and piston axis

RAS 179 A

#### 5 - Piercing block maintenance

- ✓ Check and clean

RAS 198 A (CT)

#### 6 - Reagent pack (Option)

- ✓ Replacement of the connector O rings

RAS 181 A



# OVERALL MAINTENANCE

## C - MECHANICAL OPERATION CHECK

### 1 - Liquid syringes

- ✓ Operation check  
Menu 4 - SERVICE / 5 - MECHANIC / 4 - LIQ. SYRINGES

### 2 - Air syringe

- ✓ Operation check  
Menu 4 - SERVICE / 5 - MECHANIC / 5 - PRESSURE SYR.

### 3 - Liquid valves

- ✓ Operation check  
Menu 4 - SERVICE / 5 - MECHANIC / 6 - VALVES

### 4 - Sampling needle and carriage

- ✓ Operation check  
Menu 4 - SERVICE / 5 - MECHANIC / 2 - NEEDLE U/D  
Menu 4 - SERVICE / 5 - MECHANIC / 3 - CARRIAGE L/R

## D - MECHANICAL ADJUSTMENTS

### 1 - Needle height

- ✓ Menu 7 - TECHNICIAN / 2 - ADJUSTMENTS / 7 - POSITION ADJUST. / 1 - NEEDLE / 2 - ADJUSTMENT
- ✓ Needle on the edge of the RBC chamber, press a key.

### 2 - Needle motion

- ✓ Tool GBC 218 A over the RBC and WBC chamber.
- ✓ Menu 7 - TECHNICIAN / 2 - ADJUSTMENTS / 7 - POSITION ADJUST. / 2 - CARRIAGE / 2 - ADJUSTMENT
- ✓ Lower the needle in the WBC chamber , press a key.

## E - HYDRAULIC ADJUSTMENTS

- ✓ PRIME / ALL REAGENTS
- ✓ Check for leaks

### 1 - Check vacuum

- ✓ Barflex connected instead of the tube coming from the valve 8.
- ✓ Menu SERVICE / 7 - TECHNICIAN / 2 - ADJUSTMENTS / 3 - CHECK PRESSURE  
Value : -200mbar  $\pm$  10mbar
- ✓ Check vacuum drop down during 30 secondes  $\leq$  2 mbar

### 2 - Bubbling

- ✓ Menu 7 - TECHNICIAN / 2 - ADJUSTMENTS / 8 - BUBBLING ADJUST. / BUBBLING 1 < > (default value : 175)
- ✓ Menu 7 - TECHNICIAN / 2 - ADJUSTMENTS / 8 - BUBBLING ADJUST. / BUBBLING 2 < > (default value : 120)

### 3 - Drain detection

- ✓ Voltmeter between TP30 and TP6
- ✓ Drained sensor : adjust the voltage to **4,5V  $\pm$  0,3V** by means of **R6**
- ✓ Fill the sensor with diluent : voltage falls down below **1Volt**

Procedures

RAS 173 A

RAS 172 A (OT)  
RAS 176 A (CT)

RAS 172 A (OT)  
RAS 176 A (CT)

RAS 172 A (OT)  
RAS 176 A (CT)

RAS 172 A (OT)  
RAS 176 A (CT)

RAS 174 A

M  
I  
C  
R  
O  
S  
  
6  
0

### Procedures

#### 4 - Temperature

- ✓ Thermometer plunged into the diluent container (WBC chamber if pack)
- ✓ Menu 7 - TECHNICIAN / 3 - TEMPERATURE / 2 - CALIBRATE
- ✓ Temperature diluent = temperature displayed  $\pm 2^\circ$

RAS 172 A (OT)  
RAS 176 A (CT)

## IMPORTANT

Temperature adjustment must be done 20 minutes at least after the instrument has been switched on.

#### 5 - Photometer calibration

- ✓ Menu 7 - TECHNICIAN / 2 - ADJUSTMENTS / 1 - CAL PHOTOMETER
- ✓ Hgb channel displayed on the screen, adjusted by means of **R97**, according to the room temperature :

RAS 172 A (OT)  
RAS 176 A (CT)

ROOM TPT (°C)	CHANNEL			ROOM TPT (°C)	CHANNEL		
	Mini.	Nominal	Maxi.		Mini.	Nominal	Maxi.
18	238	243	248	24	233	238	243
19	237	242	247	25	232	237	242
20	236	241	246	26	231	236	241
21	235	240	245	27	230	235	240
22	234	239	244	28	229	234	239
23	234	239	244	29	228	233	238

#### 6 - Gains

##### • WBC

- ✓ Latex thoroughly mixed
- ✓ Blank cycle to check the cleanliness
- ✓ Menu 7 - TECHNICIAN / 2 - ADJUST. / 4 - WBC ADJUST / 1 - WBC ADJUST
- ✓ Run an analysis on Latex. Adjust with **R74** to obtain
  - LYM =  $57 \pm 1$
  - GRA =  $180 \pm 2$

RAS 172 A (OT)  
RAS 176 A (CT)

##### • RBC

- ✓ Latex thoroughly mixed
- ✓ Blank cycle to check the cleanliness
- ✓ Menu 7 - TECHNICIAN / 2 - ADJUST. / 5 - RBC PLT ADJUST / 1 - RBC PLT ADJUST
- ✓ Run an analysis on Latex.
- ✓ Adjust :
  - the RBC gain with **R81** to obtain RBC =  $74 \pm 1$
  - The PLT gain with **R87** to obtain PLT =  $59 \pm 2$

RAS 172 A (OT)  
RAS 176 A (CT)

# OVERALL MAINTENANCE

## F - VOLTAGE CHECKS

### 1 - Power supply

✓ Ground on TP31, TP30 or TP29

TEST POINTS	VOLTAGE
TP 20	-12V ± 0,5V
TP 22	24V + 1.5V - 0V
TP 23	5V + 0,3V - 0V
TP21	12V ± 0,5V

### 2 - RBC, PLT, WBC threshold

✓ Ground on TP31

THRESHOLDS	TEST POINTS	VOLTAGE	POTENTIOMETERS
WBC	TP 10	350 mV ± 7	R68
RBC	TP13	400 mV ± 7	R75
PLT	TP16	180 mV ± 3	R82

### 3 - Stepper motor voltage

✓ Ground on TP31

MOTORS	TEST POINTS	VOLTAGE	POTENTIO.
Air syringe	TP36	2.50 V ± 0.05 V	R195
Liquid syringe	TP35	2.00 V ± 0.05 V	R186
Horizontal carriage	TP37	1.50 V ± 0.05 V	R177
Vertical carriage	TP38	1.00 V ± 0.05 V	R168

### 4 - Aperture voltage

✓ Ground on TP30

✓ Menu 4 - SERVICE / 7 - TECHNICIAN / 2 - ADJUSTMENTS / 2 - APERTURE VOLTAGE

✓ Test point on TP 19 : check to have **60V +2.8V -1.5V**

Procedures

RAS 171 A

RAS 175 A

RAS 175 A

RAS 172 A OT  
RAS 176 A CT

M  
I  
C  
R  
O  
S  
  
6  
0

### G - BLOOD SAMPLE RUN

#### 1 - Preliminary

- ✓ Run a STARTUP cycle

#### 2 - Repeatability

- ✓ Fresh and normal blood sample
- ✓ Run 20 consecutives analyses

## CAUTION

For MICROS 60 CT, it is mandatory to remove the cap from the sample tube to prevent from piercing several times the same cap.

- ✓ Control to have variation coefficients within the following acceptable limits :

PARAMETERS	% CV	TEST LEVEL
- WBC :	< 2,5%	at 10.10 <sup>9</sup> /l
- RBC :	< 2%	at 5.10 <sup>12</sup> /l
- HGB :	< 1,5%	at 15 g/dl
- HCT :	< 2%	at 45 %
- MCV :	< 1%	at 90 fl
- PLT :	< 5%	at 300.10 <sup>9</sup> /l
- LYM :	< 5%	at 40%(16/18 param.)
- MON :	< 10%	at 10% (16/18 param.)
- GRA :	< 5%	at 50% (16/18 param.)

CV calculated by means of the below formula :

$$\bar{X} = \frac{\sum Xi}{n} \quad SD = \sqrt{\frac{\sum (\bar{X} - Xi)^2}{n-1}}$$

$\bar{X}$  : Mean

$Xi$  : measure value

$n$  : Measure number

$SD$  : Standard deviation

$$CV(\%) = \frac{SD}{\bar{X}} \times 100$$

### H - CALIBRATION

- ✓ Calibration passed :

Check that the calibration coefficients remain within the following ranges :


LIMITS	WBC	RBC	HGB	HCT	PLT	MPV	RDW	PDW
MINIMUM	0.80	0.65	1.10	0.83	0.86	0.75	0.75	0.75
MEAN	1.00	0.81	1.38	1.04	1.07	0.94	1.00	1.00
MAXIMUM	1.20	0.97	1.66	1.25	1.28	1.13	1.25	1.25

### I - QUALITY CONTROL

- ✓ Run a blood control
  - ✓ Option : QC smart card use

# OVERALL MAINTENANCE

# MICROS 60



# DECONTAMINATION



Hématologie

01/05/98

## • CONCERNS

- Instrument decontamination before maintenance operation in the following cases
- Instrument move out of the biologic risks area
- Maintenance intervention on contaminate suspected assemblies

## • REQUIRED TOOLS

Hexagonal keys  
Clamps  
Flat screw driver  
Torx keys

## • REQUIRED PRODUCTS

Fungicidal, bactericidal, virus killing detergent spray, non corrosive for metals,  
Non plastic altering.  
Bleach solution 12°CI  
Deionize water  
Protection gloves  
Absorbant paper  
Distilled water

## • INTERVENTION TIME

1h35min

## • FREQUENCY

On request

## • SPECIFIC KIT OR CONSUMABLES

Drain and rinse kit **XEA 349 AS** for Pack equipped instrument

RAS 192 A Ind.A

M I C R O S 6 0

## • PROCEDURE

### **WARNING !**



- Disposal gloves and white coat must be worn by the operator.
- Local or national regulations must be applied in all the operations .

#### **1 - Preliminary (20min)**

- Switch on the instrument
- Run a STARTUP cycle, then a SERVICE / AUTOCLEAN
- Switch off the instrument and remove the supplying cable
- Open the instrument cover
- Spray the bactericidal cleaner on all assemblies that may provide biologic risks and wait for 10 minutes (assemblies in contact with the operator such as instrument cover, tube holder, keyboard, start key, sampling needle neighbored assemblies...

#### **2 - Manual decontamination (20 min)**

- Remove the WBC/HGB chamber cover
- Dilute the 12°cl bleach to 1 part of bleach for 4 of deionize water (1/5).
- Instrument environment must be cleaned and decontaminated.
- No sponge, nor cloth must be used. Only absorbant paper, thrown after use, in contamination bins, can be employed. For small or weak assemblies use accurate drier papers.
- All assemblies that is suspected to have contact with biologic product must be disinfected with the diluted bleach (the stainless steel must be bleached below 30°Celsius).
- Blood stains or salt marks must be cleaned with spray detergent first.
- *Concerned assemblies*
  - Outer surfaces of the instrument (perpex, covers, LCD, reagent locations...)
  - Keyboards
  - Waste connector plug
  - Liquid valve push
  - Needle neighbored assemblies
  - Tube holder assy.
  - overflow trays

Reinstall all the assemblies and setup the instrument in its initial configuration.

#### **3 - Analysis circuit decontamination (30 min)**

##### • BOTTLE VERSION

- Prepare 1 bottle containing 1/2 litre of bleach diluted to 1 part of bleach for 9 parts of deionize water (1/10).
- Prepare 1 bottle containing 1/2 litre of distilled water.
- Switch on the instrument
- Replace the reagent bottles by the diluted bleach bottle.
- Run a SERVICE / PRIME / ALL REAGENTS cycle.
- Fill a sample tube with diluted bleach to 1 part of bleach for 4 of deionize water (1/5).
- Enter the TECHNICIAN / BURN-IN function, Type in 15 cycles and leave the instrument operating until it stops (On MICROS 60 OT run 15 manual cycles).

# DECONTAMINATION

## • PACK VERSION

- Prepare 1 bottle containing 1/2 litre of bleach diluted to 1 part of bleach for 9 parts of deionize water (1/10).
- Prepare 1 bottle containing 1/2 litre of distilled water.
- Prepare one empty bottle of 1 litre for waste.
- Switch on the instrument
- Replace the reagent pack by the Drain & Rinse kit (XEA 349 AS).
- Plunge the straws into the diluted bleach bottle and the waste tube into the empty waste bottle.
- Run a SERVICE / PACK / PRIME cycle.
- Fill a sample tube with diluted bleach to 1 part of bleach for 4 of deionize water (1/5).
- Enter the TECHNICIAN / BURN-IN function, Type in 15 cycles and leave the instrument operating until it stops (On MICROS 60 OT run 15 manual cycles).

## 4 - Drain and rinse (30 min)

- Remove the 3 reagent straws from the bottle containing the diluted bleach
- Wrap the straws in absorbant paper.
- Run two prime cycles : the bleach is drained.
- Replace the diluted bleach by the distilled water bottle and re-plunge the straws in distilled water.
- Run six PRIME / ALL REAGENTS cycles (Rinse).
- Remove the 3 reagent straws from the distilled water (Wrap the straws in absorbant paper).
- Run two PRIME / ALL REAGENTS cycles : the distilled water is drained.
- Run a STAND BY cycle.
- Re-install the reagent bottles and the straws (or re-install the Pack instead of the Drain & Rinse kit).
- Switch off the instrument.
- Close the instrument cover.

MICROS 60



# DRAIN & RINSE



Hématologie

05/10/98

## • CONCERNS

Instrument Rinse and drain before  
- an extended shutdown  
- an instrument removing

## • REQUIRED TOOLS

None

## • REQUIRED PRODUCTS

Distilled water

## • INTERVENTION TIME

35min

## • FREQUENCY

On request

## • SPECIFIC KIT OR CONSUMABLES

Drain & rinse kit : **XEA 349 AS** for Pack equipped instrument  
**ADVIA** Pack optional : Reagent output protections FFZ 015 A

RAS 197 A Ind.B

M I C R O S 6 0

## • PROCEDURE

**WARNING !**

- Disposal gloves and white coat must be worn by the operator.
- Local or national regulations must be applied in all the operations .

**1 - Preliminary (5min)**

- Switch on the instrument
- Run a STARTUP cycle, then a SERVICE / AUTOCLEAN

**2 - Drain and rinse (30 min)****• BOTTLE VERSION**

- Prepare one bottle containing 1/2 litre of distilled water
- Remove the 3 reagent straws MINIDIL, MINILYSE, MINICLEAN from the bottles.
- Wrap the straws in absorbant paper.
- Run two PRIME / ALL REAGENTS cycles : the reagents are drained.
- Replace the reagent bottles by the distilled water bottle and plunge the straws into distilled water.
- Run 6 PRIME / ALL REAGENTS cycles (Rinse).
- Remove the 3 reagent straws from the distilled water (Wrap the straws in absorbant paper).
- Run two PRIME / ALL REAGENTS cycles : the distilled water is drained.
- Run a STAND BY cycle.
- Check that the diluent syringe piston is in park position (upper position).
- Remove the distilled water and install the installation kit box instead.
- Install the black plastic carriage locking clip in order to block the needle carriage (see INSTALLATION procedures :
  - MICROS 60 CT : RAS 166 A
  - MICROS 60 OT : RAS 165 A)
- Clean the reagent stains from the instrument.
- Put an adhesive tape on the tube holder door (MICROS 60 CT) to prevent from opening it.
- Switch the instrument off.

# DRAIN & RINSE

## • PACK VERSION

- Prepare one bottle containing 1/2 litre of distilled water.
- Prepare one empty bottle of 1/2 litre for waste.
- Install the Drain & Rinse kit (XEA 349 AS) instead of the Reagent pack.
- Plunge the waste tube into the empty waste bottle.
- Wrap the 3 reagent straws MINIDIL, MINILYSE, MINICLEAN in absorbant paper.
- Run two SERVICE / PACK / PRIME cycles : the reagents are drained.
- Plunge the straws into distilled water.
- Run 6 SERVICE / PACK / PRIME cycles (Rinse).
- Remove the 3 reagent straws from the distilled water (Wrap the straws in absorbant paper).
- Run two SERVICE / PACK / PRIME cycles : the distilled water is drained.
- Run a STAND BY cycle.
- Check that the diluent syringe piston is in park position (upper position).
- Remove the Drain & rinse kit.
- Switch off the instrument.
- Install the black plastic carriage locking clip in order to block the needle carriage (see INSTALLATION procedures :
  - MICROS 60 CT : RAS 166 A
  - MICROS 60 OT : RAS 165 A)
- Clean the reagent stains from the instrument.
- Put an adhesive tape on the tube holder door (MICROS 60 CT) to prevent from opening it.
- Switch the instrument off.

## NOTE

ADVIA : Install the reagent output protections (see the warning sheet form : RAL 035 A and join it to the installation kit box)

- Install the installation kit box on the pack location.

M  
I  
C  
R  
O  
S  
  
6  
0

# Piercing block maintenance



Hématologie

15/06/98

## • CONCERNS

- Sampling position and available tubes
- Piercing block description
- Maintenance

## • REQUIRED TOOLS

- Torx keys

## • REQUIRED PRODUCTS

- See decontamination procedure : RAS 192 A

## • INTERVENTION TIME

- 30 min

## • FREQUENCY

- On request

## • SPECIFIC KIT OR CONSUMABLES

- None

RAS 198 A Ind.A

MICROS 60 CT


## Piercing block maintenance

### • PROCEDURE

#### 1 - Sampling position



Diag.1



Diag.2

According to the sample tube models (Vacutainers, Microtainers ...) used and to their position into the sample tube holder, the height of the needle in its lower position can be modified.

The 3 switches (diag 1) associated to the tube holder are able to detect the sampling position according to the following principle :

Two states 0/1 are possible for the 3 switches (see diag 2) :

- Switched OFF : 1
- Switched ON : 0

The binary codes obtained from the states of the switches gives the positions of the tube holder.

For each position of the sample tube holder (from 1 to 6) corresponds a position of the needle (from 1 to 6).

SWITCH 1	SWITCH 2	SWITCH 3	Sampling position	Needle
0	0	0	Bad position of the tube holder	
1	0	0	position 1	Needle 1
0	1	0	position 2	Needle 2
1	1	0	position 3	Needle 3
0	0	1	position 4	Needle 4
1	0	1	position 5	Needle 5
0	1	1	position 6	Needle 6
1	1	1	No tube holder	

**NOTE** The code «0 0 0» means that the tube holder has been turned in between two sampling positions. The code «1 1 1» means that the tube holder has been removed.

# Piercing block maintenance

The following chart tables give the tube positions available according to the models of sample holder used on the instrument

## MICROS 60 STANDARD TUBE HOLDER



Diag.3



Diag.4

Barrillet standard (GBC061AS) Standard Tube holder	Marque Trademark	Additif Additive	Volume Volume	Vide Vacuum	Volume mort (5) Dead volume	Ref. P/N	Modèle Model	Photo N° Photo #	Réglage (1) Adjustment	Perçage tube Tube piercing	Bouchon (2) Stopper
POSITION 6 CONTROLE	R&D system	na	2ml	na	250µl	na	na	8	standard	sans bouchon without stopper	
POSITION 5 STANDARD	Becton D.	EDTA-K3	5ml	nc	390µl	368452	Vacutainer	9	standard	avec bouchon with stopper	gomme rubber
	Becton D.	EDTA-K3	5ml	2ml	400µl	367651	Vacutainer	11	standard	avec bouchon with stopper	hemogard
	Becton D.	EDTA-K3	5ml	3ml	410µl	367856	Vacutainer	10	standard	avec bouchon with stopper	hemogard
	Becton D.	EDTA-K3	5ml	3ml	400µl	367652	Vacutainer	12	standard	avec bouchon with stopper	hemogard
	Becton D.	EDTA-K3	5ml	4.5ml	410µl	367654	Vacutainer	7	standard	avec bouchon with stopper	hemogard
	Terumo	EDTA-K2	5ml	3ml	430µl	VP-053SDK	Venoject II	15	standard	avec bouchon with stopper	Ultraseal
	Terumo	EDTA-K3	5ml	5ml	460µl	VT-050STK	Venoject	17	standard	avec bouchon with stopper	gomme rubber
	Terumo	EDTA-K3	5ml	3ml	460µl (4)	VT-053STK	Venoject	na	standard	avec bouchon with stopper	gomme rubber
	ABX	EDTA-K3	5ml	4ml	480µl	ABX-3004002	na	14	standard	avec bouchon with stopper	gomme (non recommandé) rubber (not recommended)
	Greiner	EDTA-K3	5ml	2ml	370µl (4)	454087	Vacuette	na	standard	avec bouchon with stopper	hemogard
	Greiner	EDTA-K3	5ml	3ml	370µl (4)	454086	Vacuette	na	standard	avec bouchon with stopper	hemogard
	Greiner	EDTA-K3	5ml	4ml	370µl	454036	Vacuette	13	standard	avec bouchon with stopper	hemogard
POSITION 4	LDM	EDTA-KE	5ml	4.5ml	480µl	nc	nc	2	standard	avec bouchon with stopper	hemogard
	Becton D.	EDTA-K3	3ml	na	30µl	6385	Vacutainer	6	standard	sans bouchon (3) without stopper	gomme (non recommandé) rubber (not recommended)
	Terumo	EDTA-K3	3ml	3ml	30µl	VT-030STK	Venoject	16	standard	avec bouchon with stopper	gomme (non recommandé) rubber (not recommended)
POSITION 2	Sarstedt	nc	0.5ml	na	30µl	901091	nc	4	standard	sans bouchon without stopper	bouchon imperdable unlostable stopper
	ABX	nc	0.5ml	na	30µl	ABX-3001001	nc	5	standard	sans bouchon without stopper	bouchon imperdable unlostable stopper

Table 1


# CAUTION

The needle sampling position must not be modified without referring to the procedure RAS 176 A.

MICROS 60 CT

# Piercing block maintenance

## ADVIA 60 STANDARD TUBE HOLDER



Diag.5




Diag.6

Barrillet standard (GBC216AS) Standard Tube holder	Marque Trademark	Additif Additive	Volume Volume	Vide Vacuum	Volume mort (5) Dead volume	Ref. P/N	Modèle Model	Photo N° Photo #	Réglage (1) Adjustment	Perçage tube Tube piercing	Bouchon (2) Stopper
POSITION 6 CONTROLE	STRECK	nc	nc	nc	nc	nc	nc	18	standard	avec bouchon with stopper	nc
	R&D system	na	2ml	na	250µl	na	na	8	standard	sans bouchon without stopper	
POSITION 5 STANDARD	Becton D.	EDTA-K3	5ml	nc	390µl	368452	Vacutainer	9	standard	avec bouchon with stopper	gomme rubber
	Becton D.	EDTA-K3	5ml	2ml	400µl	367651	Vacutainer	11	standard	avec bouchon with stopper	hemogard
	Becton D.	EDTA-K3	5ml	3ml	410µl	367856	Vacutainer	10	standard	avec bouchon with stopper	hemogard
	Becton D.	EDTA-K3	5ml	3ml	400µl	367652	Vacutainer	12	standard	avec bouchon with stopper	hemogard
	Becton D.	EDTA-K3	5ml	4.5ml	410µl	367654	Vacutainer	7	standard	avec bouchon with stopper	hemogard
	Terumo	EDTA-K2	5ml	3ml	430µl	VP-053SDK	Venoject II	15	standard	avec bouchon with stopper	Ultraseal
	Terumo	EDTA-K3	5ml	5ml	460µl	VT-050STK	Venoject	17	standard	avec bouchon with stopper	gomme rubber
	Terumo	EDTA-K3	5ml	3ml	460µl (4)	VT-053STK	Venoject	na	standard	avec bouchon with stopper	gomme rubber
	ABX	EDTA-K3	5ml	4ml	480µl	ABX-3004002	na	14	standard	avec bouchon with stopper	gomme (non recommandé) rubber (not recommended)
	Greiner	EDTA-K3	5ml	2ml	370µl (4)	454087	Vacurette	na	standard	avec bouchon with stopper	hemogard
	Greiner	EDTA-K3	5ml	3ml	370µl (4)	454086	Vacurette	na	standard	avec bouchon with stopper	hemogard
	Greiner	EDTA-K3	5ml	4ml	370µl	454036	Vacurette	13	standard	avec bouchon with stopper	hemogard
POSITION 4	Becton D.	EDTA-K3	3ml	na	30µl	6385	Vacutainer	6	standard	sans bouchon (3) without stopper	gomme (non recommandé) rubber (not recommended)
	Terumo	EDTA-K3	3ml	3ml	30µl	VT-030STK	Venoject	16	standard	avec bouchon with stopper	gomme (non recommandé) rubber (not recommended)
POSITION 2	Sarstedt	nc	0.5ml	na	30µl	901091	nc	4	standard	sans bouchon without stopper	bouchon imperdable unostable stopper
	ABX	nc	0.5ml	na	30µl	ABX-3001001	nc	5	standard	sans bouchon without stopper	bouchon imperdable unostable stopper

Table 2

# Piercing block maintenance

## OPTIONAL TUBE HOLDER



Diag.8



Diag.9

Groove

Barrillet optionnel (GBC217AS) Optional Tube holder	Marque Trademark	Additif Additive	Volume Volume	Vide Vacuum	Volume mort (5) Dead volume	Ref. P/N	Modèle Model	Photo N° Photo #	Réglage (1) Adjustment	Perçage tube Tube piercing	Bouchon (2) Stopper
POSITION 6 CONTROLE	R&D system	na	2ml	na	250µl	na	na	8	standard	sans bouchon without stopper	na
POSITION 5 STANDARD	Becton D.	EDTA-K3	5ml	nc	390µl	368452	Vacutainer	9	standard	avec bouchon with stopper	gomme rubber
	Becton D.	EDTA-K3	5ml	2ml	400µl	367651	Vacutainer	11	standard	avec bouchon with stopper	hemogard
	Becton D.	EDTA-K3	5ml	3ml	410µl	367856	Vacutainer	10	standard	avec bouchon with stopper	hemogard
	Becton D.	EDTA-K3	5ml	3ml	400µl	367652	Vacutainer	12	standard	avec bouchon with stopper	hemogard
	Becton D.	EDTA-K3	5ml	4.5ml	410µl	367654	Vacutainer	7	standard	avec bouchon with stopper	hemogard
	Terumo	EDTA-K2	5ml	3ml	430µl	VP-053SDK	Venoject II	15	standard	avec bouchon with stopper	Ultraseal
	Terumo	EDTA-K3	5ml	5ml	460µl	VT-050STK	Venoject	17	standard	avec bouchon with stopper	gomme rubber
	Terumo	EDTA-K3	5ml	3ml	460µl (4)	VT-053STK	Venoject	na	standard	avec bouchon with stopper	gomme rubber
	ABX	EDTA-K3	5ml	4ml	480µl	ABX-3004002	nc	14	standard	avec bouchon with stopper	gomme (non recommandé) rubber (not recommended)
	Greiner	EDTA-K3	5ml	2ml	370µl (4)	454087	Vacurette	na	standard	avec bouchon with stopper	hemogard
	Greiner	EDTA-K3	5ml	3ml	370µl (4)	454086	Vacurette	na	standard	avec bouchon with stopper	hemogard
	Greiner	EDTA-K3	5ml	4ml	370µl	454036	Vacurette	13	standard	avec bouchon with stopper	hemogard
LDM	EDTA-KE	5ml	4.5ml	480µl	nc	nc	2	standard	avec bouchon with stopper	hemogard	
POSITION 3	Becton D.	nc	0.5ml	na	30µl	365975	Microtainer	3	avec tube with tube	sans bouchon without stopper	Microgard (Equipé adaptateur: autre réglage aiguille) (Equipped with adaptor: other adjustment)
POSITION 1	Becton D.	nc	0.5ml	na	30µl	365973	Microtainer	1	standard	sans bouchon without stopper	na

Table 3

MICROS 60 CT



## Piercing block maintenance

na - not applicable.

nc - not communicated.


(1) - The standard adjustment positions the needle at the bottom of the tube holder, when mentioned "with tube", the presence of the tube is mandatory for the adjustment.

(2) - More information about sample tube is available in the user manual, section 3 "SPECIFICATIONS", point 3.4.2.

(3) - The thickness of the tube stopper blocks the tube into the piercing mechanism.

(4) - These volumes have been calculated; not measured.

(5) - The "dead volume" is determined after the manual adjustment of the sampling depth and increased by 20%, except in the pediatric tubes where a 30µl volume has been fixed arbitrarily for security.




Diag.10

**NOTE** The needle height adjustments is explained in the procedure RAS 176 A.

# Piercing block maintenance


## 3 - Description

Piercing block  
right hand side view



Diag.11


Piercing block  
left hand side view



Diag.12

- ① Door lock electro magnet
- ② Door detection switch
- ③ Sampling position switches
- ④ Gearings

## 2 - Maintenance



Diag.13

- Control the correct operation of the door lock electro magnet (diag 13).



Diag.14

- To clean the piercing block it is easier to dismantle the door front panel as shown on the diag 14.



Diag.15

- The sample tube holder can be pulled out from its location as shown on diagram 15.  
- Re-install the reverse order.

# CAUTION

Decontamination procedure (RAS 192 A) must be followed to clean the sample tube holder and its location.

## CONTENTS

1. ARGOS FORMAT PRINCIPLES .....	2
1.1. Introduction .....	2
1.2. Results characteristics .....	2
1.2.1. Key .....	2
1.2.2. Result format .....	3
1.2.3. End of communication key .....	4
2. STANDARD FORMAT .....	5
2.1 Message Structure .....	5
2.2. Details about the structure .....	5
2.3. Identifier list and their formats .....	5
2.3.1. Hematologyc numeric parameters .....	5
2.3.1.1. <i>Format description</i> .....	5
2.3.1.2. <i>Identifier list</i> .....	6
2.3.2. Pathology .....	7
2.3.3. Histograms and matrix .....	7
2.3.3.1. <i>Format description</i> .....	7
2.3.3.2. <i>Identifier list</i> .....	7
2.3.4. Patient result identification .....	8
2.3.4.1. <i>Format description</i> .....	8
2.3.4.2. <i>Identifier list</i> .....	8
2.3.4.3. <i>Analysis type (\$80)</i> .....	8
2.4. Heading title .....	8
2.5. Other identifiers .....	9
3. PIN ASSIGNMENTS .....	10

# 4 . OUTPUT FORMAT


3 RS output formats are available on MICROS 60 :

- Format 1 : ARGOS / HELIOS type
- Format 2 : ABX internal format (factory use)
- Standard format
- TR off : Transmission off

## 1. ARGOS FORMAT PRINCIPLES

### 1.1. Introduction

The ARGOS format is a fixed format (406 characters for one result) including a STX and a ETX. These characters are splitted into fields representing a transmitted item.



The fields have a fixed length separated by the **OD** character.

### 1.2. Results characteristics

#### 1.2.1. Key

Total ASCII characters emitted : 406

- (-) : Space \$20
- (I) : Carriage return \$0D

- CRC : exclusive "OR" of all the transmitted bytes except ETX and STX, then an inclusive "OR" with a \$40 value.

- zzzzz : numeric field completed by zeros on the left.  
ex : 04.55 (decimal separation with a period).  
When the analyser does not transmit parameters, the field (zzzzz) is put in place of (--.--).

- Y : Alphanumeric character from \$20 to \$7F.

- # : Space (\$20) if automatic sampling. Star (\$2A) if manual sampling.

# 4 . OUTPUT FORMAT

## 1.2.2. Result format

Line 1 :	STX (\$02)		Start of text	1
Line 2 :	R (\$52)		Character "R"	1
Line 3 :	zz]		Analyser No	2 + 1
Line 4 :	YYYYYYYYYYYYYYYY]		Identification No	16 + 1
Line 5 :	YYYYYYYYYYYYYYYYYYYYYYYYYYYYYYYY]		Identification	30 + 1
Line 6 :	zz/zz/zz-zzhzzmnzsz#]		Time & Date	20 + 1
Line 7 :	zzzzz-RN]	GB	rejection & limits	8 + 1
Line 8 :	zzzzz-RN]	LYC#	rejection & limits	8 + 1
Line 9 :	zzzzz-RN]	LYC%	rejection & limits	8 + 1
Line 10 :	zzzzz-RN]	MON#	rejection & limits	8 + 1
Line 11 :	zzzzz-RN]	MON%	rejection & limits	8 + 1
Line 12 :	zzzzz-RN]	GRA#	rejection & limits	8 + 1
Line 13 :	zzzzz-RN]	GRA%	rejection & limits	8 + 1
Line 14 :	zzzzz-RN]	NEU#	rejection & limits	8 + 1
Line 15 :	zzzzz-RN]	NEU%	rejection & limits	8 + 1
Line 16 :	zzzzz-RN]	EOS#	rejection & limits	8 + 1
Line 17 :	zzzzz-RN]	EOS%	rejection & limits	8 + 1
Line 18 :	zzzzz-RN]	BAS#	rejection & limits	8 + 1
Line 19 :	zzzzz-RN]	BAS%	rejection & limits	8 + 1
Line 20 :	zzzzz-RN]	ALY#	rejection & limits	8 + 1
Line 21 :	zzzzz-RN]	ALY%	rejection & limits	8 + 1
Line 22 :	zzzzz-RN]	LIC#	rejection & limits	8 + 1
Line 23 :	zzzzz-RN]	LIC%	rejection & limits	8 + 1
Line 24 :	zzzzz-RN]	RET*	rejection & limits	8 + 1
Line 25 :	zzzzz-RN]	RET*	rejection & limits	8 + 1
Line 26 :	zzzzz-RN]	RET*	rejection & limits	8 + 1
Line 27 :	zzzzz-RN]	RET*	rejection & limits	8 + 1
Line 28 :	zzzzz-RN]	RBC	rejection & limits	8 + 1
Line 29 :	zzzzz-RN]	HGB	rejection & limits	8 + 1
Line 30 :	zzzzz-RN]	HCT	rejection & limits	8 + 1
Line 31 :	zzzzz-RN]	MCV	rejection & limits	8 + 1
Line 32 :	zzzzz-RN]	MCH	rejection & limits	8 + 1
Line 33 :	zzzzz-RN]	MCHC	rejection & limits	8 + 1
Line 34 :	zzzzz-RN]	RDW	rejection & limits	8 + 1
Line 35 :	zzzzz-RN]	RET*	rejection & limits	8 + 1
Line 36 :	zzzzz-RN]	PLT	rejection & limits	8 + 1
Line 37 :	zzzzz-RN]	MPV	rejection & limits	8 + 1
Line 38 :	zzzzz-RN]	PCT	rejection & limits	8 + 1
Line 39 :	zzzzz-RN]	PDW	rejection & limits	8 + 1
Line 40 :	ABCDEFGHIJKLMNQRSTU]		WBC 5 DIFF Flags	21 + 1
Line 41 :	LMMGGG]		WBC LMG Flags	6 + 1
Line 42 :	PSM]		PLT Flags	3 + 1
Line 43 :	CRC			1
Line 44 :	ETX (\$03)		end of text	1

# 4 . OUTPUT FORMAT

## *1.2.3. End of communication key*

- (I) : Carriage return \$0D
- CRC : the exclusive «OR» of all the transmitted bytes, except ETX and STX, then the inclusive «OR» with a \$40 value.
- zz : Number of the analyser.

# 4 . OUTPUT FORMAT

## 2. STANDARD FORMAT

The **standard format** can have a different a different number of fields according to the transmitted items setup by the user (results, curves, flags, etc...).

### 2.1 Message Structure

STX  
Size + carriage return.

Identifier followed by heading title + carriage return.  
Identifier followed by the Information associated to the heading title + carriage return.  
Remainder of the other Identifiers and Informations associated to the heading title + carriage returns.

Other heading titles + associated Informations  
.....

Identifier followed by the Checksum + carriage return  
ETX

### 2.2. Details about the structure

Size : 5 bytes representing the total amount of the data except STX and ETX.  
Heading title : An 8 characters chain preceded by a space, indicating the associated data type.  
Identifier : 1 byte (moving about \$21 to \$FF, it describes the information type which follows this indicator).  
CheckSum : Sum modulo 65535 of all the characters except ETX, STX and all the informations linked to the checksum (identifier - space - checksum) in the hexadecimal format on 4 bytes, preceded by a space.

### 2.3. Identifier list and their formats

#### *2.3.1. Hematologyc numeric parameters*

##### *2.3.1.1. Format description*

###### • Numerical field

For all indicated parameters from \$21 to \$43, the format is a numerical field of 5 digits completed with zeros on the left side (ex. : 04.55) and preceded by a space.

The unit is the one chosen by the operator.

When the parameter cannot be calculated by the analyzer, the field is replaced with (—.—).

###### • Parameter status

Following the numerical field, a first digit gives the counting rejection status or the suspicion, a second one gives the parameter value status according to high and low normalities, and to the overloading capacities.

# 4 . OUTPUT FORMAT

First digit (letter)	correspondance
R	Parameter rejected for a counting default
S	Suspicious parameter value
'space'	No anomaly observed
Second digit (letter)	correspondance
l	Parameter < to the low normal value
'space'	Parameter normal value
h	Parameter > to the high value
O	Parameter exceeding the capacity

• *Example*

5.5 millions RBC with a counting error in the standard units :  
 \$32 \$20 \$30 \$35 \$2E \$35 \$30 \$52 \$68 \$0D or «**2 05.5Rh**» + carriage return.  
 The length is fixed and is worth 2+7+1, that is to say 10 bytes for one parameter.

2.3.1.2. Identifier list

Identifiers		Parameters	Units
\$21	!	WBC	Standard - SI g/dl - SI mmoles
\$22	"	Lymphocytes	(#)
\$23	#		(%)
\$24	\$	Monocytes	(#)
\$25	%		(%)
\$26	&	Granulocytes	(#)
\$27	'		(%)

Identifiers		Parameters	Units
\$32	2	RBC	Standard - SI g/dl - SI mmoles
\$33	3	HGB	
\$34	4	HCT	
\$35	5	MCV	
\$36	6	MCH	
\$37	7	MCHC	
\$38	8	RDW	

Identifiers		Parameters	Units
\$40	@	PLT	Standard - SI g/dl - SI mmoles
\$41	A	MPV	
\$42	B	THT	
\$43	C	PDW	



# 4 . OUTPUT FORMAT

## 2.3.2. Pathology

\* *Flags associated with parameters*

• Format description

Flags are transmitted in a comprehensive mode, preceded by a space (same presentation than on the screen, i.e. dependant from the language) 2 characters which are replaced with spaces when the flag has not been detected.

• Identifier list

Identifiers		Parameters	Formats	Length
\$50	P	WBC or LMG	L1M1M2G1G2G3	2 + 12 + 1
\$53	S	PLT	PcScMc	2 + 6 + 1

## 2.3.3. Histograms and matrix

### 2.3.3.1. Format description

• Histograms

Histograms are transmitted on 128 or 256 channels, preceded by a space. They are automatically rescaled to a 223 maximum amplitude value.

The zero amplitude value is \$20, the maximum amplitude value is \$FF.

### 2.3.3.2. Identifier list

Identifiers		Parameters	Formats	Length
\$57	W	WBC	Amplitude of each channel	2 + 128 + 1
\$58	X	RBC	"	"
\$59	Y	PLT	"	"
\$5D	I	WBC thresholds	5 thresholds	1 + 20 + 1
\$5E	^	RBC thresholds	2 thresholds	1 + 8 + 1
\$5F	_	PLT thresholds	1 threshold	1 + 4 + 1

# 4 . OUTPUT FORMAT

## 2.3.4. Patient result identification

### 2.3.4.1. Format description

All the described fields have a fixed size character chain type and are completed with spaces for the non significant informations.

### 2.3.4.2. Identifier list

Identifiers		Correspondance	Formats	Length
\$70	p	Analyzer number	01	2 + 2 + 1
\$71	q	Analysis date and time	94/06/06 13h15mn31s	2 + 19 + 1
\$73	s	Analyzer sequence number	0128	2 + 4 + 1
\$74	t	Sampling mode	'O' : open tube 'C' : close tube	2 + 1 + 1
\$75	u	Identification number	1450302154275-42	2 + 16 + 1
\$76	v	Identification	SMITH Ronald	2 + 30 + 1
\$80	ç	Analysis type	Defined on 1 character (see description) 'A' CBC analysis 'D' LMG analysis From 'G' to 'Z' can be configured by the user	2 + 1 + 1

### 2.3.4.3. Analysis type (\$80)

This identifier defines the analysis type CBC, LMG to carried out on the sample. It also provides the analysis of one or several specific parameters.

The CBC analysis includes the 12 parameters of the CBC's count.

The LMG analysis includes the CBC analysis and the % and # of the 6 WBC populations.

## 2.4. Heading title

The heading title provides the data of the whole message : current hematological results or results coming from statistics. This heading title is able to drive commands that can be interpreted by the analyser or by an external computer.

Analyzers being able to communicate in the bidirectionnal mode and supporting the remote control mode, can interprete the heading title and runs the corresponding actions.

This string is a 8 characters length, preceded by a space, containing data that follows, or the command type to be carried out. The identifier is \$FF.

Identifier	Correspondance	Format	Length
\$FF	type of data packet	String of characters	2 + 8 + 1

# 4 . OUTPUT FORMAT

## Data exported by the analyzer

Data packet string	Use
RESULT	Hematological result transmission on a routine mode

## 2.5. Other identifiers

Identifiers		Correspondance	Format	Length
\$FC		Number	On 8 bytes	2 + 8 + 1
\$FD		Checksum value	hexadecimal on 4 Bytes	2 + 4 + 1

Identifier \$FC : This identifier allows the number transmission that could be an error #, a position #, a burn-in sequence # or an hexadecimal status (see "Error list").



Identifier \$FD : Checksum value : see chapter on the message constitution.

## Error list

Error N°	event linked to the analyzer
1	Operating temperature out of limits

# 4 . OUTPUT FORMAT

## 3. PIN ASSIGNMENTS



# 5. TRAINING SLIDES